



Steven M. Ruben
Appl. No. 10/662,429

Department Protein Expression
Subject _____
Name Edward Wensdike-Johnson
Address _____



 43-648
Computation Notebook
Dennison Stationery Products Co., Framingham, MA 01701


75 Sheets
11 1/2" x 9 1/2"
4x4 Quad.

0 73333 43648 8

BEST AVAILABLE COPY

Ruben EXHIBIT #70

Department Protein Expression
Subject _____
Name Edward Mendelke, Jr.
Address _____
 43-648
Computation Notebook
Dennison Stationery Products Co., Framingham, MA 01701

75 Sheets
11 1/2" x 9 1/2"
4x4 Quad.
0 73333 43648 8

Ruben EXHIBIT 2070
Ruben v. Wiley et al.
Interference No. 105,077
RX 2070

REDACTED

Transfection

5/9

127/94

Cells

Prepare dishes 1x35 mm² with 1 to 1.5×10^6
Sf9 cells/ml in Grace's media + 5% H₂O. Inactivated
FBS. Incubate overnight on Sf900

obs. my cells were prepared with 10% H₂O FBS

Check the plate to see if the cells are ok
Remove the media carefully

Add 3 ml of media without FBS 1 or 2 hrs

Prepare the transfection mix

A: 10 μ l of lipofectin
90 μ l of Media without FBS

B: 50 μ l of Media without FBS
5 μ l of plasmid SC3 (A2) at 1 μ g/ μ l
SC4 (A2)

5 μ l of Unmodified Baculogold DNA [0.1 μ g/ μ l]

Mix the contents of tube A and B and
allow the mixture to stand at least 15
minutes at room temperature before
adding it to the cells

Remove the culturing dish from the incubator
and add the transfection mixture to it drop
wise around the plate

Shake the plate forward and sideways
and incubate again 20% for 4 to 6 hours

After this time add 1 ml of media
containing 10% H₂O FBS and antibiotic
and incubate for 3 to 4 days or until
some polyhedra are visible

1.3094 Add 4 ul Xgal (stock solution 20mg/ml)
each plate (dish)

- Harvest 1 ml of supernatant into a
sterile epp

Spin at top speed to pellet any
cells and debris

- Transfer the supernatant into labeled
sterile epp tube and store at
4°C
first stock

Clone Purification

Prepare a Flat bottom 96 well with
100 ul/well of S19 cells at 1.5×10^5 cells/ml

in S1900 media + Serum H.1 or
Grace's media 10% H.1

Add 40 ul/ml (200 ul would be 10 ml)
X-Gal

- Add 100 ul of the transfected viral
supernatant to well A1 in the multi
well plate

A1 OOO OOO OOO OOO OOO OOO OOO OOO 200 ul

- Add 100 μ l of media to well A1 to A12
- transfer 100 μ l to well A to B, C etc with multichannel pipet
- Return the plates to incubate at 21°C for about 4 to 5 days or until the blue color becomes visible across the plate

make a photo copy of the plate
 choose a well at the highest virus dilution which gives a strong blue staining reaction

21.94 2 prepared again 2 microtiter plates 96 well

SF900 media + 2% HI Serum + X gal (200 μ l/10ml)
 dilution according to two protocols

Virus \rightarrow SC3 A2
 SC4 n70

21.94 6 well plates \rightarrow prepared with HI Serum for labelling

Grace's 101

SC3	WT
SC4	A2

Sf900 2%

SC3	WT
SC4	A2

2 ml of media each well
 50 μ l of each virus/well

24.94 First labelling 35 S cysteine
 \rightarrow Scott's media 35 S met

Procedure:

Remove media (complete) from 6 well plates prepared 21

add 1.5 ml SF 900 II (-Eys met - serum)
 without

incubate 1 hour

add 10 ul ³⁵S cysteine each well
 10 ul ³⁵S met " "

after 2 hours remove

50 ul → add sample buffer → freeze
 200 ul → stock -20°C

after 5 hours do the same like above
 after 24-28 h remove the rest

- 50 ul → + sample buffer -20°C
 - 800 ul → stock
 + 100 ul → harvest cells

Protein Gyl → Scott
 Sample
 Preparation

2994 Prepared 6 wells for labelling → 2nd 12/94
 SF 900 II + 2. H.I. Serum → Heimer did it

REDACTED → SF 900 II
 → Phages 101 H.I.S

cells seeded in 2 ml media + 40 ul of each virus
 from stock

2/3/94 Picked samples:

	1	2	3	4	5	6	7	8
g. cells	wt sep.	" "	SC3	Phages	SC4	SC4	Cells	" "
10:	SF 900	" "	SF 900	Phages	SF 900	wt SF 900	" "	" "
11:	SC3	SC3	SC3	Phages	SC4	SC4	SC4	SC4
12:	SC4	SC4	SC4	SC4	SC4	SC4	SC4	SC4

Heimer
 did it

2.9.94 Prepared 6 well dish for infection

REDACTED

REDACTED

media: SF 900 II + 2% FBS H.S. - 2ml

2.18.94 Collect supernatant ^{virus} Stock

REDACTED

2.18.94 - Add new media to the plate above

2.18.94 Prepared 2 dishes 6 well for

labelling: infected with ^{virus} Stock collect 2.18.94

SF900II
media

REDACTED

Another infection: in Grace media with virus
Stock from 2.18.94

REDACTED

2.23.94 → Collected the third stock of virus

REDACTED

2.23.94 → Infection: big flask 0.5 ml of C5 stock virus from 2.23.94
5 gacs for H.S.

2.23.94 → third labelling → Meiner did it

2.24.94 → Samples labelled (2.23) were collected - spin 30'
- collect SPN
- freeze -20°C

REDACTED

2.25.94 Collected stock of VIRUS CS Infected 22394
in Grace's 10% H.F.S.

STEP

3.1.94 PLAQUE PURIFICATION OF BACULOVIRUS STOCK

CS stock from 2.25.94

E12 from 51900 plate Nov 2.1.94

6. Sequentially remove supernatant medium from each plate with a pipette and replace with 1 ml of a dilution of the virus. Incubate at room temperature for 1 h, initially tipping the plates back and forth every 5 min to insure even adsorption of the virus.
7. While plates are incubating, melt the 4% Agarose Gel completely in a 70°C water bath (~10 min) and maintain at 37°C.

Protocol -> from Gibco - BEVs methods

PLAQUE PRODUCTION

Materials:

In addition to the materials listed in *General Information*, page 5, you will need the following:

Exponentially growing Sf9 cells (viability >90%)

GIBCO BRL 4% Agarose Gel (or GIBCO BRL 4% Agarose Gel with Bluo-gal)

GIBCO BRL Sf-900 Medium, 1.3X or GIBCO BRL Grace's Insect Medium, 2X; HI-FBS; and sterile, deionized, distilled water

Sealable plastic container (4" x 8" x 8")

Automated pro-pipette or other suction/reservoir system

60-mm cell culture dishes

27°C incubator

70°C water bath

37°C water bath

Protocol:

1. Prepare >1 ml of each of the log serial dilutions of the virus to be titered or purified from 1×10^{-2} through 1×10^{-6} .
2. Prepare a 50-ml stock suspension at 5×10^5 cells/ml using Sf9 cells from an exponentially growing culture at >95% viability.
3. Dispense 4 ml of this suspension into 60-mm cell culture dishes, swirling gently while dispensing to insure even dispersal of cells.
4. Allow cells to settle to bottom of plate for 1 h. Note: In serum-containing media, transport the plates gently because cells do not adhere tightly to the plate surface.
5. Inspect the cultures on an inverted microscope to confirm a 50% confluency. Adjust density as necessary.

To prepare the Grace's plaquing overlay:

8. Place a bottle of sterile, deionized, distilled water and an empty, sterile, 100-ml glass container in a 37°C water bath. Add, under sterile conditions, 20 ml of HI-FBS to a bottle of Grace's Insect Medium, 2X. Mix by pipetting, and move the mixture to 37°C water bath.
9. Combine 25 ml of Grace's Insect Medium, 2X; 12.5 ml of the sterile water; and 12.5 ml of the melted 4% Agarose Gel in the 37°C container to make the plaquing overlay. Maintain at 37°C to 42°C. Note: This example provides for a final agarose gel concentration of 1%. The final concentration may be varied from 2% to 0.2% by adjusting the water complement of the overlay solution.

Continue with protocol:

10. Sequentially remove inoculum from each plate with a Pasteur pipette and add 4 ml of the desired plaquing overlay. Note: This must be done quickly to avoid drying of the monolayer and premature gelling of the overlay solution.
11. Allow overlays to solidify (10 to 30 min) at room temperature.
12. Carefully place the plates in a sealed container with a damp cloth (to provide humidity) and incubate at $28^\circ\text{C} \pm 0.5^\circ\text{C}$ for 4 to 6 days. Note: Wild-type virus produces highly refractile, near-white plaques to the naked eye. Recombinant virus produces milky-gray plaques of slight contrast without staining or

other detection methods. Bluo-gal produces a deep-blue precipitate in the immediate area of the β -galactosidase product.

**This Page is Inserted by IFW Indexing and Scanning
Operations and is not part of the Official Record**

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

- ☐ BLACK BORDERS
- ☐ IMAGE CUT OFF AT TOP, BOTTOM OR SIDES
- ☒ FADED TEXT OR DRAWING
- ☒ BLURRED OR ILLEGIBLE TEXT OR DRAWING
- ☐ SKEWED/SLANTED IMAGES
- ☐ COLOR OR BLACK AND WHITE PHOTOGRAPHS
- ☐ GRAY SCALE DOCUMENTS
- ☒ LINES OR MARKS ON ORIGINAL DOCUMENT
- ☒ REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY
- ☐ OTHER: _____

IMAGES ARE BEST AVAILABLE COPY.

As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.