

Steven M. Ruben
Appl. No. 10/662,429

Department MOL. Biol.
Subject 615/95 - 9/14/96
Name ANN KIM # 10
Address _____

Nationel™ Brand

Computation Notebook

11 3/4" x 9 1/4", 4 x 4 Quad., 75 Sheets

43-648



0 73333 43648 8



AVERTY
DENNISON

Office Products
Chicopee, MA 01022

BEST AVAILABLE COPY

Ruben EXHIBIT #93

Department MOL. Biol
Subject 615195 - 9/14/96
Name ANN KIM # 10
Address _____

Nationel® Brand

Computation Notebook

11 3/4" x 9 1/4", 4 x 4 Quad., 75 Sheets

43-648



0 73333 43648 8



**AVERTY
DENNISON**

Office Products
Chicopee, MA 01022

Ruben EXHIBIT 2093
Ruben v. Wiley et al.
Interference No. 105,077
RX 2093

Inoculate 100ml LB + Amp Kan

3ml A2 uninduced

Incubate 37°C until $OD_{600} = 0.4-0.6$

2 1/2 hours w/ aeration

Add 100mM IPTG to 2mM

Incubate 37°C w/ aeration 4 hours

Spin culture 15 min 3K

Resuspend pellet 10ml 1M Gm HCl pH 8

Store 4°C over weekend

6/26/95

Isolate HTPANX8 S04 51bp ATG + PD10

in 6M Gm HCl pH 5 in Acrylamide
preparation Gel

in 450ul H₂O

Add 50ul Protein

50ul 0.15% Na-Dox

75ul 50% TCA

Mix well

Spin 10 min

Resuspend pellet in 25ul 0.2M NaOH

Combine into 2 Tubes

Add equal Volume 2x Dissociation

Buffer

Heat 100°C 5 min


Run 100V 1 1/2 hours

Stain + Destain marker + part of gel

Cut desired band from gel

Cut up gel slice - ready to ship

12/26/95

pg 141
book 20
pnc 29Send for Ab
production

to Pocono farm

Cindy Haub

~ 0.75mg Acrylamide
~ 0.5mg eluted on
DunndigoleSpin Crude extract of HNSAF220A2
8K 20 min

Transfer Supernatant to fresh tube

Prepare NaSO_4 Column

2ml Resin to make 1ml bed

Wash Column 10ml H_2O Equilibrate with 20ml 0.1M Gm HCl pH 8

Add supernatant. Collect Flow

Wash column 30ml 0.1M Gm HCl pH 8

Collect pH 8

Wash Column 30ml 0.1M Gm HCl pH 6

Collect pH 6

Elute Protein 5ml 0.1M Gm HCl pH 5

Collect eluted

Strip Column 30ml 0.1M Gm HCl pH 2


Collect pH 2

Prepare to run on 12% Stacking
PAGE Gel

6/26/9

450 μ l H_2O
 50 μ l Protein
 50 μ l 0.15% NaDOC
 75 μ l 50% TCA

Mix well
 Spin 10 min
 Remove Supernatant
 Resuspend pellet 240 μ l 0.2N NaOH
 Add 20 μ l 2x Dissociation Buffer
 Heat 100°C 5 min
 Run 20 μ l on gel with MW Marker
 150V 1 1/2 hrs
 STAIN / DESTAIN



Does not look
 like anything
 induced

Remake insert
 fragment

PCR Fragment

11749 + 11848

10X #2

10X BSA

H_2O

Xba I

20 μ l

6

6

27

1

60

Incubate 37°C
 Overnight

6/27/95

Inoculate 5ml LB + Amp + Kan
 with frozen stock of
 HIPAN08 51bp + PD10 D5
 from frozen glycerol stock
 Incubate 37°C w/ aeration 2 1/2 hours
 Inoculate 5ml LB + Amp/Kan w/ culture
 Add 100ul 100mM IPTG to culture
 Incubate all 37°C w/ aeration
 Inoculate 50 ml LB + Kan/Kan
 to do large scale induction

PCR Fragment
 Kba Digested

PPT with Ethanol and NaAcetate
 Spin 10min
 Remove Supernatant
 Wash pellet 70% Ethanol
 Spin 5min
 Remove Supernatant
 Allow pellet to dry
 Resuspend pellet 44ul TE
 Add 5ul Buffer M
 1ul Sph I
 Let Digest 37°C 6/12

6/28/95

Inoculate 600 ml LB + Amp/Kan
 with 35ml overnight culture
 of HIPAN08 504 51bp ATG + PD10
 D5
 OD₆₀₀ = 0.15
 Incubate 37°C w/ aeration for 2 hrs
 until OD₆₀₀ = 0.4-0.6

6/28/95

①

$OD_{600} = 0.45$
 Add 100mM IPTG to 2mM - 12ml
 Incubate 37°C with aeration
 4 hours
 Spin cultures 4K 20min
 Pour off supernatant
 Resuspend pellet in a total of
 80 ml Cell Gen HCl pH 8
 Store 4°C

Run Gels for Westerns
 12% PAGE Gels - Stacking
 15well combs

Samples: 1 HTPAN08S04 51bp ATG+PD10 D5 - UNINDUCED
 2 " " " " " " - INDUCED
 3 HNSAF22 + PGE70 H10 UNINDUCED
 4 HTPAN08S04 51bp ATG+PD10 Purified
 pH 5

Run 15 μ l each on gel

(3) (2) (1) - Spin 3ml Culture
 Resuspend in 200 μ l H₂O
 Add 200 μ l 2x Dissociation
 Buffer
 Heat 100°C 5min

Spin 5min
 (2A) - 900 μ l H₂O
 70 μ l Protein in pH 5
 100 μ l 0.15% Na DCL
 150 μ l 50% TCA

Mix well
 Spin 10 min
 Resuspend pellet 150 μ l 0.2N NaOH
 Add 150 μ l 2x Dissociation Buffer
 Heat 100°C 5min

6/28/95

Run

- 3X
- 1 - Rainbow Marker
 - 2 - HTP408504 51bp ATG + PDI0 UNINDUCED
 - 3 - " " " " INDUCED
 - 4 - HASAE22 + PDE UNINDUCED
 - 5 - HTP408504 51bp ATG - pH5.

Run 150 V 1 1/4 hours.


Set up transfer.

Transfer Buffer:

25 mM Tris pH 8.3
 192 mM Glycine
 20% METHANOL.

4°C.

Set up Sandwich



Run 100 mV
 40 min
 Run 200 mV
 20 min

Store blots in Blocking Buffer

WESTERN BLOCKING BUFFER:

3% BSA Fraction V
 0.02% Na AZIDE
 in PBS

At 4°C

34.
 PM

6/29/95

Run 2ul of Fragment with 1kb ladder



fragment looks good

Set up ligations at 21°C

5' kb / 5' Sph I

	1	2	3	4	5	6
DNA Fragment	2	2	2			
PQE70 Nco / Sph I	1			1		
PQE70 Sph I / Nco		1			1	
10X Buffer	2	2	2	2	2	2
T4 ligase	1	1	1	1	1	1
H ₂ O	14	14	15	16	16	17

Incubate reactions 16°C overnight

HTPAW08504 81bp ATG + PD10 DS protein

Spin 8K 30 min

Transfer Supernatant to fresh tubes

Crude Extract

Prepare Ni-NTA Column

6/29/92

Prepare Buffer - 3ml
 Strip Column - 30ml 0.2N NaOH
 Wash 50ml H₂O
 Charge 50ml 0.1M Na₂SO₄
 Wash 50ml H₂O - Equilibrate 30ml 0.1M HCl pH 2
 Add Supernatant to column
 Collect - flow
 Wash 40ml 0.1M HCl pH 8
 Collect - pH 8
 Wash 40ml 0.1M HCl pH 6
 Collect - pH 6
 Elute 10ml 0.1M HCl pH 5
 Collect - pH 5
 Strip - 40ml 0.1M HCl pH 2
 Collect - pH 2

Run fractions on Gel

450ul H₂O
 50ul Fraction
 50ul 0.15% NaDOC
 75ul 50% TCA

Mix well

Spin 10 min

Remove supernatant

Resuspend pellet 10ul 0.2N NaOH

Add 10ul 2X dissociation Buffer

Heat 100°C 5min


Spin 5min

Run all on 12% Stacking

PAGE Gel 150V 1hr

Stain 30min / DESTAIN 30min

6/30/95



looks like most of
the protein did
not stick onto
columns -
Try Reapplying
crude extract -
Flaw was
over with freshly
changed resin

Transform ligations

Thaw Chemically Competent M15 Cells
on ice
To 100ul thawed cells add 10ul of
ligation
let sit on ice 1 hour
heat 42°C 45 sec
Put on ice 2 min
add 400ul LB to tube
Incubate 37°C 1 hour
plate 100 ul onto LB + Amp/Kan
plates
incubate at Room Temp over
weekend.

7/3/95

Received primers for HTPAN08504
 into pCDNA 3' HA TAG
 5' HA TAG
 +Noting

#12305 HTPAN08504
 5' Bam HI + 32 5' Untranslated Region
 for 3' HA TAG in pCDNA.
 Bam HI
 GCG GGC GGA TCC TGC CTG GCT GAG TTA CAG
 CAG TC

#12306 HTPAN08504
 3' Xho I for 3' HA TAG in pCDNA
 No Stop Codon
 Xho I
 GCG GGC CCG GAG GCC ACC AAC TAA AAA GGC
 CCC GAA AAA ACT G

#12307 HTPAN08504
 5' Bam HI for 5' HA Tag and regular pCDNA
 Bam HI
 GCG GGC GGA TCC GGT ATG ATG GAG GTC CAG GGG
 GGA C

#12308 HTPAN08504
 3' Xho I for 5' HA Tag and regular pCDNA
 has authentic stop
 Xho I
 GCG GGC CCG GAG TTA TTA GGC AAC TAA AAA
 CCG CC GAA AAA AE

PCR products

8/3/95

Set-up PCR

3' HA TAG

		10x
12305	6	60
12306	6	60
10x dNTP	10	100
10x PCR	10	100
TAG	0.4	4
H ₂ O	60.4	60.4
DNA	1	10
	100.2	

5' HA TAG

		10x
12307	37	37
12308	6	60
10x dNTP	10	100
10x PCR	10	100
TAG	0.4	4
H ₂ O	68.9	68.9
DNA	1	10
	100.2	

PCR

95°C 5 min
 95°C 30 sec
 65°C 30 sec } 25x
 72°C 1 min
 72°C 9 1/2 min
 4°C Hold

Run Tylon gel



looks good
 Combine Tubes

Precipitate with equal volume PEG/xal
 Spin 10 min
 Remove & Supernatant
 Wash pellet 100% 70% Ethanol
 Spin 5 min
 Remove Supernatant
 Allow pellet to dry at RT
 10 min

7/3/95

Resuspend pellet in a total of
100 μ l TE
Rem. 1 μ l on gel

HYDRODOL 501



Store -20°C Fragment #3
Box

10 μ l

Pick colonies from transformation
(6/30) into LB + Amp/Kan. - HMSA F22

- 46 μ l #1 - 2 μ l #3
- 46 μ l #2 - 2 μ l #5

Inoculate 37°C w/ampicillin 4 hours
Set up PCR

		100x
5' PDE Vector	0.2	2.0
3' PDE Vector	0.14	1.4
10x dNTP	3.2	32.0
10x PCR	3.2	32.0
H ₂ O	23.1	231.0
Taq	0.16	1.6
Culture	2	
	32.0	300.0 - 30 fold

		100x
5' SpHMSA F22	1.6	16.0
3' PDE Vector	0.14	1.4
10x dNTP	3.2	32.0
10x PCR	3.2	32.0
H ₂ O	21.7	217.0
Taq	0.16	1.6
Cult.	2	
	32.0	300.0 30 fold

7/3/95


PCR Program # 106

95°C	5 min	} 30x
95°C	20 sec	
75°C	20 sec	
72°C	1 min	
72°C	7 1/2 min	
4°C	Hold	


7/4/95 Holiday

7/5/95 Run 10 out of PCR reactions
with 1K ladder - HUSAF 22

5' PGE + 5' TGE Vector Specific





5' PGE + 5' TGE




5' PGE + 5' TGE

5' PGE + 5' TGE



7/5/95

5' + 3' PCR primers



Does not look like the ligation worked.

Try Redtransforming
10/7/95

Set up Digestion for pc-DNA HTPAN08504
(p339)

5' HA TAG	
Fragment	20
10x #2	15
Eam HI	1
Xho I	1
H ₂ O	23
	50

3' HA TAG	
Fragment	20
10x #2	15
Eam HI	1
Xho I	1
H ₂ O	23
	50

Incubate 37°C overnight

7/6/95



Pick more colonies from 4/9/95 - H1MSAF 22
into LB + Amp / Kan
Incubate 37°C 4 hours

Set up PCR

5' Sense	1.6	160
3' PDB	0.14	14
10x dNTP	3.2	320
10x PCR	3.2	320
Taq	0.16	16
H ₂ O	226.217	2170

5/6/95

PCR Prod 16 HMASAP22
 Run 10 million Gel w/ 1 kb ladder




looks good.
 incubate 200ul of
 LB + Amp/Kan
 with (+) and (-) AC(2)
 incubate 37°C w/
 aeration overnight

LMP Gel Purification
 HTPANOTSO 3' HA Tag
 5' HA tag

Load onto 0.8% LMP Gel.
 Cut out gel slice
 Ex. Take Picture
 Here Clean fragment

HCPANOSOM

4/6/95



Gene Clean

Add 100 μ l NBT

Heat 55°C 5min

Add 10 μ l Glass milkLet incubate at R.T.
5min

Spin 10sec

Remove Supernatant

Resuspend pellet in
500 μ l Wash Buffer

3X Spin 10sec

Remove Supernatant

Spin 10sec

Remove Supernatant

Resuspend pellet 40 μ l TE

Heat 55°C 5min

Spin 1min

Transfer Supernatant to fresh tube

Resuspend pellet 20 μ l TE

Heat 55°C 5min

Spin 1min

Transfer Supernatant

Run 1 μ l of Fragment on gel with

1 Kb ladder

Store -20°C in Fragment #5



Will need Vector to do ligations
None premade w/ correct
Restriction. So set up digests

pCDNA3 HA	10 μ l
10x2	5
H ₂ O	23
Prim / Xho I	1/1
	50


pCDNA3 HA	17
10x	5
H ₂ O	26.0
Prim / Xho	1/1
	50

7/6/95 pCDNA I 1.1
 LOX #2 5
 H₂O 46.9
 Bam HI 1.1
 50

incubate 37°C overnight

7/13/95

Remainder of Digest on Gel
 17.6 kb ladder



pCDNA Should be
 24.7 kb

pCDNA 5'43' HA
~~Tag~~
 Race DNase insert
 2000 bp

pCDNA 3' Nucleotides

looks like not completely
 digested.
 Add 1 µl more. Enzyme.
 incubate 37°C. Seal
 on day 2
 incubate 37°C over
 weekend.

Rec'd Clone from Sazie

Clone I.D. HE2PM 21 # 2(RB)
 Requested By Ann Kim
 Transfer Folder 07/05/95 transfer 2
 Processed By Laurie Riedel


Start Ah
 7/13/95

7/7/95

Small scale inductions
 - from overnight cultures from pg 42.
 Inoculate 200 μ l CB7 Amp Kan with
 15 μ l of Culture -
 incubate 37°C w/ aeration 2 hrs.
 $OD_{600} = 0.4 \pm 0.6$
 Add 100 ml IPTG to 2 mM - 4 μ l
 incubate 37°C w/ aeration 4 hrs.
 Spin culture 10 min
 Remove supernatant
 Resuspend pellet in 10 μ l H₂O
 Add 10 μ l of Dissociation Buffer
 Heat 100°C 5 min
 Spin 2 min
 Run 5 μ l on gel with marker and
 uninduced culture as 15 well
 of 12 to stacking gels.
 Run 150 V 1.3 hrs
 Stain over weekend at RT.

7/10/95

Run pCDNA Digests on 0.8% LMP Gel
 with 1 kb ladder




Cut out Gel Slices
 Gene Clean

Resuspend Gel Slice
 800 μ l NaI
 Heat 55°C 5 min
 Add 5 μ l glass beads
 let sit at RT 2-3 min
 w/ occasional mixing
 Spin 10 sec

7/10/95

Remove Supernatant
 Resuspend pellet 100ul Wash Buffer
 Spin 10 Sec
 2x { Remove Supernatant
 Add 100ul Wash Buffer
 Spin 10 Sec
 Remove Supernatant
 Spin 10 Sec
 Remove Supernatant
 Resuspend pellet 15ul TE
 Heat $+55^{\circ}\text{C}$ 1 min
 Spin 10 Sec
 Transfer Supernatant to fresh tube
 Resuspend pellet 15ul TE
 Heat 50°C 1 min
 Spin 10 Sec
 Transfer Supernatant to New tubes
 Run 2ul on gel with 1 kb ladder



Vectors look good
 Etc

store -20°C
 Vector #2 box


Set up ligation
 10ul HPRAN08504
 5' + 3' HA Tags
 Fragments
 (pg 15)

7/10/95

	1	2	3	4	5	6	7	8	9
HTPANO6504 5' HA Tag Bam/Kho	1	1		1					
HTPANO6504 3' HA Tag Bam/Kho			1		1				
PCDNA Bam/Kho	4					4			
PCDNA 5' HA Bam/Kho		3					3		
PCDNA Bam/Kho 3' HA			2	1				2	
10x Buffer	2	2	2	2	2	2	2	2	2
H ₂ O	12	12	12	16	16	13	14	15	17
T4 Ligase	1	1	1	1	1	1	1	1	1

Incubate 16°C overnight

DESTAIN Gels - 37°C - 2 hrs - HMSA P22 PDE70



7/11/95

	5' HA	3' HA	pCDNA
DAADNA	17	27.8	1.8
10X BSA	5	5	5
10X BSA	5	5	5
H ₂ O	22	11.2	37.2
Enzyme	11	1	1
	50	50	50.0

Digest with Xho I and Bam HI

Incubate 37°C O/N

HTPAND8504 pCDNA constructs.
 digestion from 7/10/95 pg 47

Transform ligations (10 µl) into
 100 µl of Chemically Competent
 * DH5α cells.

Thaw cells on ice
 Aliquot 100 µl into sterile tubes
 Add 10 µl of ligation
 Mix by pipetting
 Incubate on ice 1 hour
 Heat 45°C 45 sec
 Place on ice 2 min
 Add 400 µl LB
 Incubate 37°C 1 hour
 Plate 100 µl + 200 µl into
 1 B₂ Amp plates
 Incubate 37°C overnight

7/12/95

Pick transformed cells into LB+amp
HTPANO8504 + pCDM1 constructs

- 1 HTPANO8504 + pCDM1
- 2 HTPANO8504 + 5'HA Tag pCDM1
- 3 HTPANO8504 + 3'HA Tag pCDM1

Anneal 37°C w/ 1 Overhang 4 hrs
Setup PCR.

for #1 & 2

for #3

H2O	2
FPI3	0.01
10x dNTP	3.2
10x PCR	3.2
Taq	0.15
H ₂ O	21.4
Culture	2
	<u>32</u>

H2O	2.5
FPI3	0.01
10x dNTP	3.2
10x PCR	3.2
Taq	0.15
H ₂ O	20.9
Culture	2
	<u>32</u>

PCR Prog 6L

95°C	5 min
95°C	20 sec
55°C	20 sec
72°C	1 min
72°C	7 1/2 min
4°C	Hold

⊕ HTPANO8504
Plasmid
⊖ H₂O only

Run 10ul on Gel with 1kb Ladder

No Positives seen - Not even ⊕
plasmid PCR'd. Re-do Tomorrow
with other primers

Carry out 07/13-07/17

Add 13% PEG 8000 / 1.6M NaCl - 850, l
 Mix well.
 Store -20 °C O/N

7/13/95

Ratio of Galactokinase

T ₁ 4°C 30-180 sec		T ₂ 4°C 30-180 sec		T ₃ 4°C 30-180 sec	
1-1	$\frac{0.6582 - 0.6551}{150 \text{ sec}}$ $-6 \times 10^{-6} / \text{sec}$	1-2	$\frac{0.6919 - 0.6752}{150 \text{ sec}}$ $1.1 \times 10^{-4} / \text{sec}$	1-3	$\frac{0.7737 - 0.7784}{40 \text{ sec}}$ $-3 \times 10^{-5} / \text{sec}$
2-1	$\frac{0.6665 - 0.5709}{150 \text{ sec}}$ $6.4 \times 10^{-4} / \text{sec}$	2-2	$\frac{0.6664 - 0.4669}{150 \text{ sec}}$ $9.3 \times 10^{-4} / \text{sec}$	2-3	$\frac{0.8013 - 0.7769}{40 \text{ sec}}$ $2.8 \times 10^{-4} / \text{sec}$
3-1	$\frac{0.7006 - 0.7033}{150 \text{ sec}}$ $3.5 \times 10^{-5} / \text{sec}$	3-2	$\frac{0.6917 - 0.6701}{150 \text{ sec}}$ $1.3 \times 10^{-4} / \text{sec}$	3-3	$\frac{0.7730 - 0.7730}{40 \text{ sec}}$ $-5.5 \times 10^{-9} / \text{sec}$
4-1	$\frac{0.6325 - 0.5931}{150 \text{ sec}}$ $6 \times 10^{-4} / \text{sec}$	4-2	$\frac{0.6203 - 0.4792}{150 \text{ sec}}$ $9.4 \times 10^{-4} / \text{sec}$	4-3	$\frac{0.7334 - 0.7557}{40 \text{ sec}}$ $3 \times 10^{-4} / \text{sec}$
5-1	$\frac{0.6841 - 0.5678}{150 \text{ sec}}$ $7.8 \times 10^{-4} / \text{sec}$	5-2	$\frac{0.6273 - 0.4796}{150 \text{ sec}}$ $9.0 \times 10^{-4} / \text{sec}$	5-3	$\frac{0.9937 - 0.7835}{40 \text{ sec}}$ $5.8 \times 10^{-5} / \text{sec}$
Avdh GalK	$\frac{0.6105 - 0.1350}{150 \text{ sec}}$ $3.2 \times 10^{-3} / \text{sec}$	Blank	$\frac{0.7250 - 0.7212}{150 \text{ sec}}$ $4.5 \times 10^{-5} / \text{sec}$	Avdh GalK 10% G4	$\frac{0.533 - 0.1695}{40 \text{ sec}}$ $4 \times 10^{-3} / \text{sec}$

7/14/95

Ri-PCR HTA1308 SOL pCDNA constructs
 Run on gel with 1 Kb ladder

See pg 54

7/17/95

H5/HQ1	0.7765 - 0.6000	HQ DT	0.8019 - 0.6740
DT - 20°	150 sec	PS - 20°	150 sec
H1/H95	$1.8 \times 10^{-3}/\text{sec}$	H1/H95	$8.5 \times 10^{-4}/\text{sec}$
Galk + 10%	0.6771 - 0.1012	Blank	0.8696 - 0.0653
Galk - 20°	150 sec		150 sec
H1/H95	$3.8 \times 10^{-3}/\text{sec}$		$2.9 \times 10^{-3}/\text{sec}$

Inoculate 5ml TB + Amp
with HTPANOS pCDNA Cultures

- 1: A1, A2, A3, A4
- 2: D1, D2, D3, D4, E1, E2, E5
- 3: E9, E10, F2, G1, G2, H3, H4

Incubate 37°C w/ aeration o/n

7/18/95

Boiling Menu pups - HTPANOS pCDNA Construct

Spin 2ml Cultures

Remove Supernatant

Resuspend pellet 800 μ l STET + RNase

lysozyme

Heat Samples 70°C 1 min

Spin 5 min

Remove pellet

Add - 800 μ l 13% PEG 8000 / 1.0M NaCl - mix well

Spin 10 min


Remove Supernatant

Wash pellet 1 ml 70% EtOH

Spin 5 min

9/18/95

Remove Supernatant
 Allow pellet to dry at RT 10 min.
 Resuspend pellet in 200 μ l TE
 Run gel on gel with 1 kb
 ladder




Set up Digestions:

	μ l	Bar	
DNA	10		10x
NotI	3		50
H ₂ O	10.8		302.4
Bar	0.1		1.8
μ l	0.8		1.8
	30		

Incubate Digests at 37°C O/N

Run HMAF22 + PDE
 with 1 kb ladder


Digests on 1% Gel
 Buff II / Sph I



Cannot tell
 anything with
 this digest
 - Too much chromosomal
 - Try Digests with
 EcoRI and HII

7/19/95

Rem. 1/2 of fragment on gel with
1 Kb ladder



HGDSH43 T7+ATG
ready for T7
Store -20 Fragment
#2 Box


4/44

HTUSB02 Fragment 5
Set up Digest

3' HA 12408 + 12409					RASHA 12406 + 12407					PAZ 12404 + 12402					PGE 12403 + 12405				
DNA	10	DNA	10	DNA	10	DNA	10	DNA	10	DNA	10	DNA	10	DNA	10	DNA	10	DNA	10
10x	5	10x	5	10x	5	10x	5	10x	5	10x	5	10x	5	10x	5	10x	5	10x	5
H ₂ O	34	H ₂ O	34	H ₂ O	34	H ₂ O	34	H ₂ O	34	H ₂ O	34	H ₂ O	34	H ₂ O	34	H ₂ O	34	H ₂ O	34
Boon		Xho	1	Boon	1	Xho	1	Boon	1	Xho	1	Boon	1	Xho	1	Boon	1	Xho	1
	50		50		50		50		50		50		50		50		50		50

Incubate BPPC 4 hrs - Add Alkaraile
enzyme and incubate ON at 37°C

HTRANDS - pCDNA Constructs - Boon/Xho Digest



7/20/95

HTPANO'S Mini preps

Spin 50ml culture

8K 15 min

Resuspend pellet in Resuspension

Buffer + RNase - 10ml

Add 100ml Lipid Buffer

mix gently

let sit 15 min at RT

Add 10ml Neutralization Buffer

mix well

let sit on ice 20 min

Spin 8K 30 min

Transfer Supernatant to fresh tube

Add 10ml 30% DCC and 10ml

5M NaCl

mix well

Store 4°C 20 min

Spin 8K 30 min

Pour off Supernatant

Resuspend pellet 10ml 70% Ethanol

Spin 8K 15 min

Pour off Supernatant

allow pellet to dry slightly

Resuspend in 2 x 400ul TE

precipitate DNA

1/10 vol 3M Na acetate 2x vol 100% Ethanol

Store -20°C O/N

No 1/10/95

7/21/95

Pick Colonies into 1 B + Amp Medium

1-24

6-12

9-48

2-24

10-12

10-46

3-12

7-24

4-12

8-24

7/21/95

Clones look good -
Monday..... inoculate for miniprep

HTPANO6

Spin DNA 10 min

Pour off Supernatant

Wash pellet 2x 70% Ethanol

Spin 5 min

Pour off Supernatant

let pellet Dry at RT

Resuspend in 100 μ l of TE Total

Ready to clean up

7/24/95

Inoculate 5ml TB + Amp w/
Cultures from HT4SBO2 Legation

1- HT4SBO2	3' HATAG B/A	- A1, A2, A3, A4
2	3' HATAG B/A/B	- C1, C2, C3, C4
3	5' HATAG B/A	- E1, E2, E3, E4
4	5' HATAG X/B	- F5, F6, F7, F8
5	pcDNA B/A	- G1, G2, G3, G4
6	pcDNA X/B	- H5, H6, H7, H8
7	PA2 B/A	- A1, A2, A3, A4
8	PA2 X/B	- C8, C10, C12, D1

4 pcDNA	3' HATag	15
pcDNA	5' HATag	16
pcDNA		17
PA2		18

Inoculate 37°C with Aeration
over night

7/24/95

Transform ligations from 7/20/95
 ligations # 90 10 - HT45802 + PCE60

Thaw M15 sup 4 Cells on ice
 Add 100 μ l of Chemically Competent
 cells to 100 μ l of ligations
 incubate on ice 1 hour
 heat 42°C 45 sec.
 place on ice
 Add 400 μ l LB
 incubate 37°C 1 hour
 plate 100 + 250 μ l onto LB Amp/Km
 plates
 incubate at 37°C O/N.

HTPAND8 cDNA Constructs

2x Phenol / Sarg (1:1) Extract
 2x Sarg extract
 PPT DNA

1/10 vol 3M Na Acetate pH 5.3
 2 vol 100% ethanol

Spin 10 min

Pour off supernatant

Wash pellet 1000 μ l 70% ethanol

Spin 5 min


Pour off supernatant

Allow pellet to dry at R.T 10 min

Resuspend pellet 200 μ l TE

Run gel on gel with 1 kb ladder

7/24/95



Plasmid looks good
Set up Dye 15

DNA	4
10x #2	3
H ₂ O	22
Bam HI	0.5
Xho I	0.5
	<hr/> 30

Inoculate 37°C O/N.
Run on gel tomorrow

7/25/95

RADIATION PROBLEM HOUND

- VARIABLE to travel freely around 3rd floor 9670
- Take plate out of 37°C HT45B02 + P0B00
- Take culture tubes out of 37°C HT45B02 + P0A2 + P0B00
- Boiling minutes:
 - Spin 2ml culture 5 min
 - Remove supernatant
 - Resuspend pellet 1000 μ l STEA
 - Boil and lysozyme
 - Heat 100°C 2 min
 - Spin 10 min

7/25/95


Remove Pellet
Add equal volume 13% PEG/NaCl
mix well
Store 20°C

7/26/95

NO WORK - Cleaning of
of Radioactivity by RSO

7/27/95

Spin B. Mon. 5 min
Remove Supernatant
Wash pellet 1.5 ml 70% ethanol
Spin 5 min
Remove Supernatant
Allow Pellet to Dry at RT 15 min
Resuspend in 200 µl TE
Resuspend on gel with 1 kb ladder



1-8 HT4S602 + 3'HA CAN
A1, A2, A3, A4, C1, C2, C3, C4
9-16 HT4S602 + 5'HA pCDNA
E1, E2, E3, E4, F1, F2, F3, F4
17-24 HT4S602 + pCDNA
G1, G2, H1, H2, I1, I2, J1, J2
25-32 HT4S602 + PAZ
A1, A2, A3, A4 C8, C9, C10, C11
33-34 3'HA pCDNA
35-36 pCDNA
37-38 PAZ

9/27/95

Set up Digests


1, 24 + 33, 34, 35 Digest Bam/Xho.

25 - 32 + 36 Digest with Bam/Xba.

	10	29X		10	10X
DNA	3	971	DNA	3	30
IDV#2	16.8	2487.2	IDV#2	0.1	1
H ₂ O	0.1	2.9	Bam H7	0.1	1
Bam HI	0.1	2.9	Xba I	0.1	1
Xho I	0.1	2.9	H ₂ O	16.8	16.8
	30ul	20ul/tube		30ul	20ul/tube

Incubate at RT - 3 hrs

Run HPAAG Digests on gel w/ 1Kb ladder




Digests look good

Submit for Seq.
with internal
primers

8/10/95

Run ~~11~~ ¹ ~~11~~ of plasmid on gel w/ 1.0x
 plasmid looks good
 Set up Digests



Bam / Xba I

DNA	1
10x	3
H ₂ O	25.6
Bam HI	0.2
Xba I	0.2

30ul

Incubate 37°C 5/1h

Transfections

Remove Spent media
 2x Wash Cells in 3ml DMEM without Serum


make transfection Cocktail

300 ml	
3 ml	Nutrodome
0.75 ml	DEAE Dextran (20mg/ml stock)
0.8 ml	Chloroquine (5mg/ml stock)
25.8 ml	DMEM

To 4 ml of ~~same~~ Transfection Cocktail
 Add 10ug DNA

HT45602	pcDNA	Constructs
HT16N08	pcDNA	Constructs

8/16/95



Remove Transfection Cocktail

Add 1ml ~~clear~~ shadow media:
 DMEM + 10% FCS + 10% DMSO
 + 100 units - 100 units
 2.5 min - 10 min
 10 min - 10 min

Remme Shale med

Add 1 ml corn Em


Remove media

Add 3ml complete DMEM

Incubate 30°C for 72 hrs
in CO₂ incubator.

Ren Dixon on oil
w/ 100 Ladder

Boose good.
Submit HTR/BBZ
for Signing



8/17/95

8/30/

Digest with Bam / Xba - HEZPM21 + PAZ

DNA	2
H ₂ O	15.0
10x#2	2
BamHI	0.2
Xba I	0.2
	<u>20ul</u>

Incubate 37°C 2hrs

Run 10ul on gel w/ kb ladder



Submit for Sequencing with internal primers and with PAZ specific primers.

HT4SBO2 + PAZ

Construct sent to Protein Expression
Submit for Sequencing to Confirm

FAS ligand - HTPANXX

Antibody - Rec'd Rabbit Ab
from Pocoro farms

Ready for Western Detection

FAS LIG	RABBIT #11940	→ #1
	#11941	→ #2

Pre Bleed - 6/30/95

Test Bleed #1 - 8/5/95

8/30/95

Western blots from 8/28/95
 (See pg 31-32)

Remove 24 blots from Blocking Buffer
 To 10 ml from blocking Buffer
 add 100 μ l Anti Serum - 1:100
 Dilution
 incubate ON at RT Shaking

#1 = Test bleed #1 - #11941
 #2 = Test bleed #1 - 11940
 3 = Prebleed #11940
 4 = Prebleed #11941


8/31/95

FAS LIG Western

Cells -
 Pour off 1^o Ab
 Rinse filters in 1x PBS
 Wash 1x in 10 ml 1x PBS - 5 min
 at RT w/ Shaking
 Add 2^o Ab - Rabbit Anti Rabbit
 Peroxidase - Dilute 1:2000 in 1x PBS
 Incubate at RT w/ Shaking for
 1 1/2 hrs
 Rinse filters 1x PBS
 Wash 1x in 10 ml 1x PBS - 5 min
 at RT w/ Shaking
 Rinse filters in 20 ml
 50 mM Na₂HPO₄ pH 7 for
 5 min at RT w/ Shaking

8/30/95

Add Substrate:

12.5 ml 50 mM $\text{Na}_2\text{H}_2\text{PO}_4$ pH 725.0 mg β -NADH25.0 ~~ml~~ ^{ul} - Phenol equilibrated
in TRIS pH 7.0.8.5 ~~ul~~ ^{ul} 30% H_2O_2 solution375 ~~ul~~ ^{ul} NBT - Nitro Blue Tetrazolium
(10 mg/ml)Incubate at RT till color change
develops.Stop Rxn with dH₂O.Dry and label blots
FAS Lig - WESTERN.Prebleed
#11941Prebleed
#11940

1° Ab.

Test Bleed #1
#11940Test Bleed
#11941

lane #1 - Rainbow Marker

#2 - HTPAN08504

51bp ATG in PD10

UNINDUCED

#3 -

#4 - HUSAF22 + PQE70

INDUCED

#5 -

HTPAN08504 51bp ATG in PD10

purified

NTA Column - pH 5.0.

See Pg 141

8/31/95

PCR HIPANOS Constructs

POE60

12975	5bp	2
2865		0.4
10x		10
10x		10
H ₂ O		76.3
Taq		0.3
DNA	5µl	1
		100

12986	185bp	2
2865		0.4
10x		10
10x		10
H ₂ O		76.3
Taq		0.3
DNA	5µl	1
		100

PCR HELPM21 + PA2 constructs

T7 Bac		0.1
12805	3' Nco I	2
10x		10
10x		10
H ₂ O		76.6
Taq		0.3
DNA	5µl	1
		100


PCR PRTA # 42

95°C	5 min	
95°C	30 sec	
55°C	30 sec	25x
72°C	1 min	
72°C	7 min	
4°C	10 min	

Total Volume PCR Run
500µl

8/3/95

Run 5ul of PCR product on
Gel with 1kb ladder



Add equal Volume
to 13% PEG / A.M.H. (500ul)
per Merscoel
Precipitate at 4°C O/N.

Resuspend 5ul in LB + Amp/Kan
for ~~over~~ inductions
HYPER 516 + Kan + PD10
HECM21 + PCE60 #1
Incubate 37°C w/aeration 10/N

9/1/95

Spin DNA fragments 10min
in 70% Ethanol Wash.


Spin 5min

in 70% Ethanol

Allow pellet to dry 10min in RT

Run 5ul on gel with 1kb ladder

Resuspend pellet in 10ul TE



fragments are ready
to digest

Submit HF2 PM21 / HPA2

TA + 3.4kb fragment
cfr 70/7 to

9/1/95

HTPANOS 516pA14 + PD10
HE2PM21 + POEGG #1

Inductions

To 300 ml LB + Amp + 0.1%
DN culture till OD \approx 0.8
Incubate w/ aeration at 37°C
for 2 hrs till OD₆₀₀ = 0.4-0.6

Add 100 mg IPTG to 2 mM \approx 6.5 ml
incubate w/ aeration 4 hrs
at 37°C

Spin 5K, 15 min
Remove supernatant
Resuspend pellet in 60 ml H₂O pH 8
HE2PM21 \approx 40 ml
HTPANOS \approx 20 ml
Store 4°C over weekend.

9/4/95 Labor Day


off.

9/5/95

Spin Cultures 5K, 20 min
Transfer supernatant to fresh
tube
Prepare fresh Ni column
2 ml Resin Bed
20 ml 0.1 M Ni SO₄
20 ml d H₂O

9/5/95

Run 10 μ l on gel with 11 kb ladder



looks like one of the primers is not working
 - probably #12496 - T7 PGE60
 made the New PGE60 T7 primer
 & redo.

Digest New HTUS302 51bp + 18.5bp
 PGE60 constructs

PCR Fragment	10
10x	5
H ₂ O	34.0
Nco	0.5
BspFI	0.5
	50 μ l

Incubate
 37C O/N.

9/6/95

Precipitate Digests
 Add 150 μ l TE - 200 μ l 13% PEG/NaCl
 mix well
 Spin 10 min



9/16/95

Wash pellet 70% Ethanol 5000
Spin 45 min
Remove Supernatant
allow pellet to Dry at R.T. - 10 min
Resuspend pellet in 100 µl TE
Set up ligation

	1					
5bp NIB	4	—	4	—	—	—
185bp NIB	—	4	—	4	—	—
PGE NIB	2	2	2	2	2	2
10X	2	2	2	2	2	2
T4 ligase	1	1	1	1	1	1
H ₂ O	11	11	13	13	15	17

Incubate 16°C O/N.

DE STAIN Gels



Computer Work

9/7/95

Transform ligations

Thaw 20 ml sup 5 Cells on ice
 100 μ l of Chemically Competent
 cells add 10 μ l ligation
 Let sit on ice 1 hour
 Heat 42°C 45 sec.
 Place on ice
 Add 400 μ l LB
 Incubate 37°C 1 hour
 plate 100 + 200 μ l onto LB + Amp/Km
 plates
 Incubate 37°C O/N

- HTPACK 51bp - Ab from Peconofarms + Block O/N

1:200 Dilution in Blocking Buffer
 Computer Work -

ORF for
 HNFAR64
 HTUSB02
 HSEBN09
 HNSAF22
 HE2PM21

9/8/95

Pick clones into 200 μ l LB + Amp/Km
 in 96 well dish
 N 90 of 185bp ATG and 51bp ATG
 Incubate 37°C 4 hrs O/N
 Ready to do PCR to check for
 inserts.

9/8/95

PCR to check for H. TRANS. 504
51bp & 185bp ATG

51bp			185bp		
		90x			90
5' Primer H ₂ O	2	190	5' N ₂ O	2	1
3' Primer PDE	0.1	9	3' Primer PDE	0.1	1
10x	3.2	288	10x	3.2	28
10x	3.2	288	10x	3.2	28
H ₂ O	21.3	1917	H ₂ O	21.3	19
Taq	0.2	18	Taq	0.2	18
cutt	2	—	cutt	2	—
	3.2	—		3.2	—

PCR program lab


95°C 5min
 95°C 20sec
 55°C 20sec } 30x
 72°C 1min
 72°C 7.5min
 4°C Hold

Run 10ul on gel with 1/16 ladder

H. TRANS. 504 51bp ATG + PDE

H. TRANS. 504 51bp ATG + PDE

9/8/95



Looks like
transformations
worked very
well

A-H = HTPAN088x
51bp ATG + PGE

A'-H' = HTPAN08804
185bp ATG + PGE

Inoculate 3ml LB + Amp/Kan with
10ul of 4hr culture
incubate at RT over weekend
1-2 loops both 51bp + 185bp ATG

FAB Legend AB. See pg 8/31/95 (p. 131)
#11940 & #11941 large plasmid

Pour 100 Buffer + 1° Ab
Wash 1x PBS 10 min
Add Anti Rabbit Peroxidase (2° Ab)
at 1:2000 in 1x PBS
Incubate at RT w/ shaking 1 1/2 hrs
Wash filter 1x PBS 10 min
Rinse 50mM Na2H2PO4 pH 7 for
5 min
Add Substrate 12.5ml 50mM Na2H2PO4 pH 7
25 mg B-NADPH
25 mg Phenol
8.5ml 30% H2O2
0.375 ml NGT

9/8/95

Shape cont'd color appears to desired
 Distances.
 Stop Reaction by adding dH₂O
 Rinse filters 2x dH₂O
 Dry on Whatman paper



Lane

- 1 HTPAN08 51bp ATG - 1u
- 2 HTPAN08 51bp ATG - 1u
- 3 #MSAF22 HPG70 - 1u
- 4 HTPAN08 51bp ATG -
purified protein

9/11/95

To the cultures left at RT merge
 weekend add 100mM IPTG to
 10mM - 4 hours
 incubate 37°C
 Spin 10000g 5 min
 Resuspend pellet 75ul H₂O
 Add 75ul 2x dissolving
 Buffer
 Run 12 ul on 10% Gel with
 Rainbow Color

9/11/95

Run 100V 1/2 hours

STAIN Overnight

#P001- HE2PM21 5' Bam HI-Pac

GGG CGC GGA TCC GCC ATC ATG GCG GCA GCA GTG GTC C

PCR New HE2PM21 PA2 Construct
with old New 5' Bam primer
#1301

4-Undo

4-Undo

2-Undo

4

in

13.101

2

3 x dNTP

2

10 x dNTP

10

10 x PCR

10

H₂O

74.8

Taq

0.2

DNA (original)

1

100

total of 100 µl Run

PCR Program FS

95°C 5min

95°C 30sec

55°C 30sec

72°C 1.5min

72°C 7.5min

4°C Hold



25x

9/12/95


DESTAIN
HTPACR

GELS

9/12/95




HTPAN08 51bp ATG + PEG60



HTPAN08 185bp ATG + PEG60

Looks like HTPAN08 185 bp ATG induced
 put w/ 51bp.
 Do Boiling preps & Digest
 in 500ul 500ul TB + Amp / 100ul
 with 1-18 of 51bp + 185bp
 incubate 37°C 0.1N w/ 100ul

HEZPM21 PA2 Construct
 Run 10 min on gel with
 1kb ladder




looks good
 PPT with equal volume
 13% PEG / 100ul
 Spin 10 min
 Pour off Supernatant
 wash pellet

146

9/12/95

TNT RESULTS 9/12/95			
INVESTIGATOR	SAMPLE NAME	EXPECTED PRODUCT SIZE KDa	OBSERVED PRODUCT SIZE KDa
MARIO CEPEDA	HTSFH75	30 OR LESS	NO PRODUCT OBSERVED
LAURIE INSCORE	HHFFK18	44	48 (major product), 44, 31, 23
LAURIE INSCORE	HHFFK18 PA2 2/25	33	NO PRODUCT OBSERVED
LAURIE INSCORE	HHFFK18 PA2 3/18	33	NO PRODUCT OBSERVED
ANN KIM	HE2PM21 PA2	44	36 (faint band)
T7 POSITIVE CONTROL	DNASE 02-105	33	33
T7 NEGATIVE CONTROL	WATER	NONE	NONE
T3 POSITIVE CONTROL	HCAC193	33	33
T3 NEGATIVE CONTROL	WATER	NONE	NONE




Set up PCR of DNA TNT constructs
using the 0.125 μl T7 PCR 100/70
primers (see pg 136 9/5/95).

HE2PM21		HHFFK18		HT45B02		HTPND
4042	0.1	4042	1.5	4042	1.5	4042
3' RNA	1.5	3' RNA	1.5	3' RNA	2.5	3'
10x dNTP	10					
10x PCR	10					
H ₂ O	76.2					
Taq	0.3					
DATA	1					

9/12/95

PCR Prep lab.

Run 10 μ l on gel w/ 1 kb ladder

looks like
only HmSAF22
worked well
- Try different
primers?

T₂ Promoter is very long
4042 - 3' end


PCR PPT 1 fmSAF22

9/13/95

To Rxns from 9/12/95

HE2PM21	HTMSB02	516 ₂ HTPAB08	1856 ₂ HTPAB08
4032 0.1	4032 0.1	4032 0.1	4032 0.1
3' 1.5	3' 1.5	4.5	1.5
10x 1	10x 1	1	1
10x 1	10x 1	1	1
H ₂ O 6.1	H ₂ O 6.1	6.1	6.1
Taq 0.3	Taq 0.3	0.3	0.3
10	10	10	10
10 μ l / Tube			

9/13/95

Run 10 μ l on Gel w/ 1 kb ladder

all PCR Expts
seem to have
worked

PEG PPT
- Add Equal Vol
PEG / NaCl
- 70% ethanol wash
- Rept 450 μ l TE

GET Cos Cells from Peter Hudson (PLH)
Frozen Stock - Plate into T25 flask

Boiling mini Prep.

Spin 2ml culture 5min
Resuspend pellet 750 μ l
STE + RNase / Lysozyme

Boil 5min

Spin 10min

Remove Pellet

Add 750 μ l PEG/NaCl

Mix well

Spin 15min

Remove Supernatant

Wash pellet 70% ethanol
1ml

Spin 5min

Remove Supernatant


Allow Pellet to Air Dry 15min

Resuspend pellet 200 μ l TE

Run 2 μ l on gel w/ 1 kb ladder

150

9/13/95



Set up Digests


DNA	5	36x
10x	3	10x
H ₂ O	21.6	77.6
Nco I	0.2	7.2
Bgl II	0.2	7.2
	30 μ l	

Incubate 37°C O/N

9/14/95

Run 2ul of PPT'd PCR fragments in Gel w/ 1 kb ladder

Run Boiling Municipal Digests on
and water 1 lb. each -
10:15



Nco / Boy II
Diplo

All Look
Like They
Disputed
Contract
Should Induce

Rem Western of 5/6p Induction
(see pg 144 9/12/95)

Run 150 V
Transfer blot 100 V 1 hr

Block 2 hours with Blocking
Buffer

Pour 100 lbs
Add 10 lbs at 1:20 - Fac
Test Blood # ~~1199~~ 1199/4240

Uncut ^{Beading} 0/1000 w/ Shalimar

See pg 3 Book # 11
Lab Note book 405

**This Page is Inserted by IFW Indexing and Scanning
Operations and is not part of the Official Record**

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

- ☐ BLACK BORDERS
- ☐ IMAGE CUT OFF AT TOP, BOTTOM OR SIDES
- ☐ FADED TEXT OR DRAWING
- ☐ BLURRED OR ILLEGIBLE TEXT OR DRAWING
- ☐ SKEWED/SLANTED IMAGES
- ☐ COLOR OR BLACK AND WHITE PHOTOGRAPHS
- ☐ GRAY SCALE DOCUMENTS
- ☐ LINES OR MARKS ON ORIGINAL DOCUMENT
- ☒ REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY
- ☐ OTHER: _____

IMAGES ARE BEST AVAILABLE COPY.

As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.