



Steven M. Ruben
Appl. No. 10/662,429

Department MOL. BIOL. - PROT. EXPR.
Subject 12/95 - 6/5/95
Name ANN KIM # 9
Address _____
 43-648
Computation Notebook
Dennison Stationery Products Co., Framingham, MA 01701

75 Sheets
11 1/4" x 9 1/2"
4x4 Quad.
0 73333 43648 8

BEST AVAILABLE COPY

Ruben EXHIBIT #92

Department MOL. BIOL. - PROT. EXPR.
Subject 12/95 - 6/5/95
Name ANN KIM # 9
Address _____



43-648

Computation Notebook

Dennison Stationery Products Co., Framingham, MA 01701



0 73333 43648 8

75 Sheets
11 1/4" x 9 1/4"
4x4 Quad.

Ruben EXHIBIT 2092
Ruben v. Wiley et al.
Interference No. 105,077
RX 2092

2

HIPANOS 185 bp + PDE60

pg 151 Book 8 #236

1/24/95

Spin HIPANOS 185 bp + PDE60 in
6M GdnHCl pH 8 - 8K 20min.Equilibrate N. 504 column with pH 8
6M GdnHCl

Apply Supernatant to Column - Collect Flow
 Wash 30 ml pH 8 6M GdnHCl - Collect pH 8
 Wash 30 ml pH 6 6M GdnHCl - Collect pH 6
 Elute 5 ml pH 5 6M GdnHCl - Collect pH 5
 Strip 30 ml pH 2 6M GdnHCl - Collect pH 2.

Add 50 μ l of ~~eluate~~ collected material
 to 450 μ l H₂O
 50 μ l 80% 0.5% NaDOC
 75 μ l 50% TCA.

Mix well

Spin 5 min

Remove supernatant


Resuspend pellet in 15 μ l 0.2N NaOH / 15 μ l 2X Buff

Heat 100°C 5 min

Run 20 μ l on gel with Rainbow Markers
on 12.5% gel with at 150V 1 1/2 hrs.

Stain 1 hr

De Stain 1 1/2 hr.



1	Uninduced
2	Crude Extract
3	flow
4	pH 8
5	pH 6
6	pH 5
7	pH 2

After 2 Send some
 Dialyze to protein exp
 20 min column

185 bp

1/24/95

11-6 PQECO

HTPAW08185 bpt PQECO

pg 748 Book 8 #236

2/2/95

11-6 - Reapply to column + Strip again
do try & Purify more
Start Dialyzing ~~the~~ ^{some} in Dialysis
Tubing

3M Gn HCl / Hepes 5 hrs
1.5 M Gn HCl / Hepes Over night
Reapply 43ml to Ni Sep Column - to Renature
over column of HTPAW08185 bpt PQECO.

2/3/95

Change Buffer.

1 M Gn HCl / Hepes 10 hrs.
0.5 M Gn HCl / Hepes Over Weekend.

Carrie will finish

2/6 - 2/10 Vacation

2/3/95

HTPAW08185 bpt + PQECO

Strip Column that has ~~leak~~
unatured Protein in 2
Strips in Imidazole Elution
Buffer: 50 mM NaPO₄ pH 6
250 mM Imidazole
300 mM NaCl
10% Glycerol.

2 Strips at 2.5 ml each.

Run on Stacking gel with LNK
Marker - 12.5% gel 150 V


Stain 1/2 hr.
DeStain 1 hr.

4

1L6 PQE60 / HIPANOS 185bp + PQE60


2/3/95

HIPANOS 185bp + PQE60



Store 4°C.

1L6 PQE60 + PDD 2nd Column Extraction



43

29

18.4

H1PAN08/H1PB411 - PD10 - PAZ

5

pg 152 Book 8 # 236

1/26/95

Received Primers for
H1PAN08 & H1PB411
Reprotect 55°C O/N.

1/27/95

PCR Fragments

H1PB411 PAZ

9130 5' Bam.	1
3' Xba New.	1
10x dNTP	50
10x PCR	50
Taq	2
H ₂ O	397
DNA (10ng/ul)	1
	<hr/> 500ul

100ul / Reaction

H1PAN08

	51 PAZ	185 PAZ	51 PD10	185 PD10
3146	2	2	2	2
9111:	20	9112: 20	9113: 24	9114: 20
10x dNTP	50	50	50	50
10x PCR	50	50	50	50
Taq	20 2	20 2	20 2	2
H ₂ O	375	375	371	375
DNA 10ng/ul	1	1	1	1
	<hr/> 500	<hr/> 500	<hr/> 500	<hr/> 500
	100ul / Rxn			


6 HTPA008 | HTPB411 PDB / PAZ

1/27/95

Run PCR:

95°C	5min	} 25X
95°C	30 sec	
55°C	30 sec	
72°C	1 min	
72°C	7 1/2 min	
4°C	Hold	

Run 5ul of rxn on gel with 1kb ladder



1	HTPB411	5' Bam / 3' Xba PAZ
2	HTPA008	51bp PAZ
3		185bp PAZ
4	↓	51bp PAZ
5	↓	185bp PAZ

Precipitate Reactions
Add equal Vol
13% PEG / NaCl

Spin 10min
Remove Supernatant

Wash pellet 1000ul 70% Ethanol
Drop Remove Supernatant
Dry pellet 5min at RT.
Resuspend pellet in 100ul TE.
Set up Digests

DNA (PCR rxns)	10 ul
10X #2 Buffer	5
H ₂ O	34
Bam	0.5
Xba	0.5
	<hr/> 50ul

Incubate 37°C 4hrs.
Run on 0.8% LMP gel with
1kb ladder.


HTPANO8 / HTPBY11

PD10 / PAZ

7

1/27/95

Cut out bands
Take picture:



1	HTPANO8	51 bp	PAZ
2	↓	185 bp	PAZ
3		51 bp	PD10
4	↓	185 bp	PD10
5		HTPBY11	PAZ

Gene Clean fragments

- Resuspended in 40ul TE

Set up ligations.

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17
9111 + 3146	6	6							6								
9112 + 3146			6	6						6							
9113 + 3146					6						6						
9114 + 3146						6						6					
HTPBY11 (PAZ)							6	6					6				
10X Buffer	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2
H ₂ O	9	9	9	9	9	9	9	9	11	11	11	11	11	15	15	15	17
T ₄ Lig (14ul)	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
PD10 B/x 10ul					2	2								2			
PAZ B/x 1/6	2		2				2								2		
PAZ B/x 1/27		2		2				2								2	

10X Buffer
T₄ ligase
H₂O

2
1
9
12ul

18X
1836
18
162

12ul / Tube

Add Appropriate
of Vector Fragment
or H₂O.

8

Hunters 11/18/94

PAZ / PDR

1/27/95

Incubate ligations 16°C overnight

1/27/95

1/30/95

Transform ligations

100ul of Chem Competent Cells
10ul of ligation rxn.

PAZ Constructs DH5α

⊕ Control PAZ

PDR Constructs M15 rep 4

⊕ Control PDR

Incubate on ice 1 hr

Heat 42°C 45 sec

Place on ice

Add 400ul LB

Incubate 37°C off 1 hr.

Plate 200ul + 300ul onto

LB+Amp plates for all ligations
1-8 onto 150mm plates
for ligations 9-17 plate 100ul
onto LB+Amp 100mm plates.

Incubate 37°C O/N

1/31/95

Plates look good

No colonies on ⊖ Control plates

Colonies on ⊕ Control plates

- Inoculate plates colonies onto LB+Amp
for PAZ Constructs -
200ul of LB+Amp in 96 well dish- Inoculate colonies onto LB+Amp Kan
for PDR Constructs
200ul of LB + Kan Amp in 96 well Dish

HTPB411 / 1/27/95 PA2 / ~~PD10~~ PD10

9

1/31/95

B+Amp:

① 9111/3146 + PA2 1/6 48
 ② 9111/3146 + PA2 1/27 13
 ③ 9112/3146 + PA2 1/6 35 } Plate # 1

③ 9112/3146 + PA2 9/6 12
 ④ 9112/3146 + PA2 1/27 13
 PA2 1/6 2
 PA2 1/27 1
 ⑤ HTPB411 + PA2 9/6 48
 ⑥ HTPB411 + PA2 1/27 12 } Plate # 2

LB+Amp/Ran:

⑤ 9113/3146 + PD10 48
 ⑥ 9114/3146 + PD10 48 } Plate # 3

Incubate plate 4 hrs 37°C. with aeration
 Setup PCR's.

9111	1	70x	9112	1	55x
3146	0.1	70	3146	0.1	55
10xPCR	3.2	7	10xPCR	3.2	5.5
10xdNTP	3.2	240	10xdNTP	3.2	176
H ₂ O	22.35	240	H ₂ O	22.35	176
Taq	0.15	150.5	Taq	0.15	1229.25
Cult.	2	10.5	H ₂ O	0.15	8.25
	<u>32</u>	30ul/tube		<u>2</u>	
				32	30ul/tube
9113	1.0	(50x)	9114	1	(50x)
3146	0.1	50	3146	0.1	50
10xPCR	3.2	5	10xPCR	3.2	5
10xdNTP	3.2	160	10xdNTP	3.2	160
H ₂ O	22.35	160	H ₂ O	3.2	160
Taq	0.15	117.5	Taq	22.35	117.5
Cult.	2	7.5	Cult.	0.15	7.5
	<u>32</u>	30ul/tube		<u>2</u>	
				32	30ul/tube

10

HTPB411 / HTPB411

PD10 / ~~PD10~~ PA2

1/31/95

HTPB411 5'Bam	0.1	65X
3'Xba	0.1	6.5
10X dNTP	3.2	6.5
10X PCR	3.2	208
H ₂ O	23.25	208
Taq	0.15	1.511.25
Cult.	2	9.75
	<u>32ul</u>	<u>30ul / tube.</u>

PCR.

95°C	5min	
95°C	20sec	
55°C	20sec	} 30X
72°C	1min	
72°C	7 1/2 min	
4°C	Hold.	



PA2 controls
for 9111, 9112
9 HTPB411

2/1/95

Run Reactions on 1% TAE Agarose
gel with 1kb ladder.

9113 A1-D12

9113/9114 E1-H12



HTPAND8 / HTPB411 + PD10/PAR



11

9114
EI-H12

9113 9111
A1-F1

9111
A1-F1

2/195





9111
A1-F1

9112
F2-H12

9112
F2-H12 A1-C1

PAR Controls


HTPB411
DI-H12



HTPB411
DI-H12

HTPB411
DI-H12

H



12

HTPA008 / HTPB411

PA2 / PD10

1/1/95

Inoculate 200 ml LB + Amp/Kan
with 20 μ l of ⑤ clones to do
small scale inductions (micro)
incubate at 37°C w/aeration
2 hrs.

Add 5 μ l of 100 mM IPTG To 2 mM IPTG
final conc.

incubate 4 hrs 37°C w/aeration

Spin 10 min

Resuspend pellet in 5 μ l H₂O

Add 15 μ l 2X Buffer

Store -20°C till tomorrow.

Inoculate 5 ml TB + Amp
with PA2 ⑤ clones.

①	-	11	9N1 / 31218 + PA2 9/6
②	-	8	9N1 / 31416 + PA2 1/27
③	-	10	9N2 / 3146 + PA2 9/6
④	-	10	9N2 / 3146 + PA2 1/27
⑤	-	1	HTPB411 + PA2 9/6

Incubate 37°C w/aeration O/N.

2/2/95

Run Protein gels 15% Stacking
15 μ l of Samples + 1 MW marker.

150V 1 1/2 hrs

Stain 30 min 37°C

DeStain 1 hr 37°C


HTPA008/HTPB411

PA2/PD10

18

HTPA00851bp + PD10

2/2/95




HTPA00851bp

Should produce
a protein.

32.5 kD

HTPA008 51bp + PD10




HTPA008185bp

Should produce
a protein: 27.7 kD

HTPA008
185bp + PD10

HTPA008
185bp + PD10




looks like HTPA008
185bp + PD10 -
done
has induction

14

HTRAND8 | HTRBY11 PDI0 / PAZ

2/2/95



looks like
all induced
except \ominus control.

On large scale
prep...

HTRAND8S04 (SSD) + PDI0

Chk3

Do Binding Mini Preps. of PAZ Constructs

2ml culture Spin 2 min

Remove Supernatant

Resuspend pellet in 750ul STET +

RNase + Hypoxanthine

Heat 100°C 2 min

Spin 10 min

Remove pellet

Add 750ul 13% PEG-8000 / 1.6M NaCl

Mix well

Spin 10 min

Remove Supernatant

Add 1ml 70% Etanol to wash pellet

Spin 5 min

Remove Supernatant.

Allow Pellet to dry at RT 10 min

Resuspend pellets in 150ul TG

Run 2ul on 1% TBE gel with 1 kb


Add ladder

HTPAN081 HTPB411

PD10/PAZ

15

2/2/95



Mini preps all
look good.

Check for inserts.
By digesting w/
Bam / Xba.

DNA	5 μ l
IOx	30 μ l
H ₂ O	21.8 μ l
Bam	0.1 μ l
Xba	0.1 μ l
	30 μ l

Incubate 37°C
over night

2/3/95

Run 10 μ l on gel w/ 1 Kb ladder.

39422340 HTPB411
3PAZ



All look correct
- 1 kb
Select 2 to sequence
with internal
primers.

(1)	1, 2	} RPO6 FPI6
(2)	1, 2	
(3)	1, 2	
(4)	1, 2	

HTPB411 + PAZ - FPR3

Carrie checked Seq. Some looked good so did Mark

16

HTPANO8/HTPB411 + PAZ PD10

1/13/95

- Inoculate 3ml TB+Amp w/ cultures 1-1, 1-2, 2-1, 2-2, 3-1, 3-2, 4-1, 4-2 & HTPB411 + PAZ. Make Glycerols.
- Inoculate 200ml LB+Amp w/ HTPB411 + PAZ to do Maxi Prep
- Inoculate 10ml LB+Amp Kan w/ induced culture of HTPANO8 185bp + PD10 (#12) - Do large scale induction

1/14/95

- Made Glycerol Stocks - 80°C Protein Expression Box #1
- HTPANO8 185bp + PD10 #12 -
 Inoculate 300ml LB+Amp/Kan w/ 300ml of O/N culture
 Incubate 37°C 3hrs - until $OD_{600} = 0.4-0.6$
 Add 100mM IPTG to 2mM - 6mM
 Incubate 37°C 4 1/2 hrs w/ aeration
 Spin 5K 15 min
 Resuspend Pellet in 30ml of 6M GuHCl pH8
 Store at 4°C O/N.
- HTPB411 + PAZ Maxi Prep

Run Maxi
 Spin culture 6K 20min
 Pour off supernatant

(pg 27)

HPRANUS / HTPBY11

27

(page)

2/15/95

HTPB11 FPAZ

Resuspend pellet in 10ml of P1 + RNase

Let sit RT 5min

Add 10ml P2 while mixing

Add gently 10ml P3 while mixing

Let sit on ice 20 min

Spin 20 min 8K

Equilibrate Tip-500 with 10ml QBT

Apply Supernatant to Equilibrated

Column

Wash Column 2x with 30ml

QF

Elute DNA ca 15ml QF

Add 0.7X (10.5ml) of isopropanol

Mix Well

Spin 9K for 25min

Pour off Supernatant

Wash pellet with ice cold 70% Ethanol

(10ml)

Spin 9K 10min

Pour off 70% Ethanol

Allow Pellet to Dry at RT

Resuspend pellet in 100ul TE

Run 1ul on gel

Read OD 260/280 at 1:200 Dilution

abs	abs	bkg abs	260.0 nm	280.0 nm
260.0 nm	280.0 nm	320.0 nm	280.0 nm	260.0 nm
-0.0063	-0.0028	-0.0022	6.9003	0.1449
0.1519	0.0964	0.0232	1.7591	0.5685

1.52 ug/ul + 608 ug
Total


Store 4°C Plasmid Box #2

Sequence w/ internal primers to confirm sequence

28

HTPB4/11 & HTPAND8

2/14/95



looks good.
- See if sequence is good

2/15/95

Digest @ 1mg DNA w/ 1u Bam/Kla
to see if insert "pops" out.

DNA (250ng)	4
10X #2	3
H ₂ O	22.6
Bam	0.2
Kla	0.2
	<hr/> 30

Incubate 37°C

Submit for sequencing w/ internal
primers. HTPB4/11 PA2 RA/FP

RP50A
2nd sequence
primer

RP01A
RP03A
RP04A
RP05

RP06A
RP07
FP08
RP09

RP10
RP11
PA2
RP13

FP14
FP15
FP16
FP17

FP18
FP19
FP20
FP21A

FP22A
FP23C
RP50

HTPAN08 / HTPBYN

29

Submit for Sequencing w/ internal

2/15/95

HTPAN08 51bp + PAZ

HTPAN08 185bp + PAZ.

Clone Plasmids to Steve to submit
to Protein Expression for baculovirus.

51bp.

RP12	FP14	RP01
FP13	RP05	FP18
RP10	FP08	RP06
FP16	FP17	RP50.

185bp

FP16	FP17	RP50.
FP14	RP01	
RP05	FP18	
FP08	RP06	

HTPAN08PA51RP/FP

HTPAN08P18SRP/FP.

HTPAN08 185bp + PD10 #12 large scale

~~inductions~~ inductions.

Spin Culture 20 min 8K.

Prepare NiSO₄ Column.

Wash Apply 2ml Resin to Column

Wash 20ml H₂O.

Add 30ml 0.1M NiSO₄ to Change

Wash 30ml H₂O.

Equilibrate with 30ml 0.1M Gn HCl pH8

Apply Supernatant - Collect Flow

Wash 45ml pH8 - Collect pH8

Wash 45ml pH6 - Collect pH6

Elute 6ml pH5 - Collect pH5

Strip 45ml pH2 - Collect pH2.

Run on 15% Acrylamide Stacking gel.

490ul H₂O

20ul 0.1M Gn HCl prep-

75ul 50% TCA

50ul 0.15M NaDOC.

mix well

Spin 10 min

30

HTPB411 & HTPA0508


2/15/95

Remove Supernatant
 Resuspend pellet 10ul 0.2N NaOH
 Add 10ul 2X Dissociation Buffer
 Boil 5min - lost pH 8 + pH 2 < samples
 mixed with water.
 Spin 5min
 Run 20ul on gel with 1MW
 Protein markers
 150V 1 hour.

Blain 30 min. 37°C
 DESTAIN overnight.

2/16/95

HTPA008 185 PDD #12



Protein
 looks good

Reapply 5ml of
 pH 5 to Column
 add 30ml pH 8
 and to Protein
 Expression &
 have renatured
 over column

Inoculate 500ul of LB + Amp/Kan
 with 20ul of HTPB411 + PDD

H1PBIII & H1PANO8

31

2/16/95

Incubate 37°C w/ aeration
2 hrs.

Add 100 mM IPTG to 2 mM
10 µl.

Incubate 37°C w/ aeration
overnight.

2/17/95

Spin Cultures 2 min

Remove Supernatant.

Resuspend culture 20 µl H₂O.



Add 20 µl 2X Dissociation Buffer.

Heat 100°C 5 min

Spin 5 min.

Run 10 µl on 10% Stacking
gel.

Accidentally ran 1 kb Marker instead
of Rainbow Marker.



H1PBIII + PQE60

Run 150V

1/2 hr

Stain
DESTAIN


Ab + Rb

2/17/95

32

HTPAND8 / HTPB411

2/2/95



HTPB411 + PDE60

Try growing up 1 for induction -
large scale

2/21/95

Inoculate 5ml LB + Amp/Kan
with HTPB411 + PDE60.
2-2 & 3-4
Incubate 37°C

Transform - HTPAND8 5/6p + PD10 9/13/95
HTPAND8 5/6p + PDE60 1/5/95
into M15 Chemically Competent Cells

Thaw M15 on ice.
To 100 µl of Cells add ligations
Incubate on ice 1 hour.

43

11.1.1.1 + SV
Pg 26


33

2/17/95

Mix well
Spin 10min
Remove Supernatant
Resuspend in 10ul 0.2N NaOH
Add 10ul 2x dissociation Buffer
Heat 100°C 5min
Run all on gel.
150V 1 hour.

STAIN 30min 37°C
DESTAIN 30min 37°C

Take Picture.



looks like I have protein.

Looks slightly contaminated

Try all elutions
to try and clean
up prep?
Need RT in PDE60

2/21/95

pick 48 more clones from the
2/13/95 into 200ul of LB Amp/kan

H1PAN08 | H1PB411

43

pg 32

2/21/95

Heat 48°C 45 sec.
Place on ice
Add 400 μl of LB
Incubate 37°C 1 hour.
Plate 300 μl onto LB+Amp - 750mm plate
100 μl onto LB+Amp/Kan 100mm plate
Incubate 37°C O/N.

2/22/95

PCR H1PAN08 ST + PQE60 48 AEI-H12
H1PB411 ST + PD10 48 AI-D12
into 200 μl LB+Amp/Kan
Incubate 37°C w/ aeration O/N.

H1PB411 + PQE60.
10 300 μl LB+Amp/Kan. Add
3 μl of O/N culture. 2-2/3-4
Incubate w/ aeration 3 hrs until
 $\text{OD}_{600} \sim 0.4 - 0.6$
Add 100 mM IPTG to 2 mM (10 μl)
Incubate 37°C 4 1/2 hours
Spin cultures 7K 20 min
Remove supernatant
Resuspend pellet in 30 μl CoM6mHCP
pH 8
Store 4°C O/N.

2/23/95

Do PCR of H1PAN08 clones.

44

HTRAN08504 | HTR0611

2/23/95

HTRAN08504 51bp + PD10
(50X)



9113	1.5	75
3146	0.1	5
10x dNTP	3.2	160
10x PCR	3.2	160
H ₂ O	21.9	1095
Taq	0.1	5
Culture	2	
		30ul/tube

HTRAN08504 51bp + PD10

9113	1.5	50X
3146	0.1	75
10x dNTP	3.2	5
10x PCR	3.2	160
H ₂ O	21.9	160
Taq	0.1	1095
Culture	2	5
		30ul/tube

PCR program # 66

Run 10ul on gel w/ 1 Kb ladder



HTRAN08 PD10 A1-D12
A2, A7, B12, C5, C7, D5, D9, D10
Inoculate 2ml LB+amp/kan
with 20ul O/N culture

HTRAN08 PQE60 E1-H12
E1, E2, F1, F2, G1, G3, H1, H2
Inoculate 2ml LB+amp/kan with
20ul O/N culture
Incubate at RT Overnight

HTPBX/HTPBII

45

2/23/95

HTPBII + PGE₂ 2-2 + ~~3~~ 3-4
Spin Culture 8K 10 min.
Transfer Supernatant to fresh tube
(Crude Extract).

Prepare Column.
Prepare Column with 2ml NTA Resin.
Wash 30ml H₂O
Change 30ml 0.1M NaSO₄
Wash 30ml H₂O
Equilibrate 30ml 6M GnHCl pH 8.

Apply Supernatant to Column.
Collect as flow.
Wash 30ml 0.1M GnHCl pH 8
Collect as pH 8.
Wash 30ml 6M GnHCl pH 6
Collect as pH 6.
Elute protein 6 5ml 6M GnHCl pH 5
Collect as pH 5.
Strip Column 30ml 6M GnHCl pH 2
Collect as pH 2.

Store 4°C O/N.

2/24/95

HTPBII + PGE₂.


400ul H₂O
20ul of eluted Protein in 6M GnHCl
50ul 0.8% NaDOC
75ul 50% TCA
Mix well
Spin 5 min
Remove Supernatant
Resuspend pellet in 0.2N NaOH - 10ul

46

HTPAN08 | HTPB411

2/24/95

Add 10 μ l 2X Dissociation Buffer
 Heat 100°C 5 min
 Spin 1 min
 Run 15 μ l on gel with Uninduced
 and Rainbow marker
 12.5% Acrylamide Stacking gel
 150V 1 hour



STAIN 37°C
 30 min
 DESTAIN 30 min
 37°C

Take Picture

looks like isolated
 protein

Store at 4°C
 Ask Steve about
 how to get
~~SDS~~ signature

HTPAN08 516p + PDE/PD10

Place tubes at 37°C w/aeration

for 2 hrs
 Add 100mM IPTG to 2mM 4 μ l
 incubate 37°C 4 hours

Spin 1ml 2 min
 Remove supernatant

Resuspend pellet in 30 μ l H_2O

Add 30 μ l 2X Dissociation Buffer
 Run 15 μ l on gel with Rainbow Marker
 150V 1 hour

HTPA N08 | HTPBY11



47

2/24/95

Stain 30 min 37°C
DESTAIN over weekend at RT

2/27/95

Take picture



A2
C7
D5
D9

Induced

Show up to
make Glycer.
Pick 1 to 10
Large Scale
Prep.

C7 - large Scale
Inoculate 5ml
2B + Amp/Kan w/
C7. incubate O/N
w/ aeration. Make
Glycerols of A2, C7, D5, D9
Does not look
like anything
Induced -)

Try running
again - need only
10ul this time

48

HTPAN08 HTPB4

2/27/95

Re Run HTPAN08 51bp + PGE6e
10ul.

Run 10ul each of the 1st + 2nd
imidazole elution of HTPAN08 185 + PD10
#2

2/28/95



HTPAN08 185bp + PD10 #12


looks like both fractions
have protein.

~~HTPAN08 185bp + PD~~

inoculate 300ml LB +
Amp Kan with
HTPAN08 51bp + PD10 #7
to do Oagor Map
incubate 37°C 24h/overnight
O/N

Re Run HTPAN08 51bp + PGE6e

15ul 1 lane
Stain / DESTAIN



Does not look like
there was any
induction
- try other clones
- shut it up

E 4, 5, 6, 8
F 3, 4, 5, 6
G 4, 6, 7, 8
H - B, C, D, E, 4, 5, 6, 7

inoculate 2ml LB +

HTPAN08 / HTPB11

49

2/28/95

Wash 20ml of culture from 423
log cell. D
let sit at RT O/N

HTPAN08 51bp + PD10 CF Oxygen
maxi - See Ref 42 - Ded
along side O/L to SV POE(0).

HTPAN08 51bp + PD10 CF
Large Scale Induction
Inoculate 300ml LB + Amp Kan
5ml of O/N culture.
Incubate 37°C w/ aeration
until $OD_{600} = 0.4 - 0.6 - 2\frac{1}{2}$ to
3 hours
Add 100mM IPTG to 2mM (10ml)
Incubate 37°C w/ aeration
4 1/2 hours.
Spin culture 5K 15 min.
Pour off supernatant.
Resuspend pellet in 30ml
6M Gdn HCl pH 8
Store 4°C O/N.

3/1/95

Incubate HTPAN08 51bp + POE(0)
in 2ml LB + Amp Kan 37°C
w/ aeration
Incubate until ~ 2 hours
Add 100mM IPTG to 2mM (4ml)
Incubate 37°C 5 hours
Spin 1ml culture
Remove supernatant
Resuspend pellet in 40ul H₂O
Add 40ul 2X Dissociation Buffer
Aber 1/1/95

50

H118W 1H1864

3/1/95

Heat 100°C 5 min

Spin 2 min

Run 100V on gel w/ 2 MW / Marker

180V

40 min

12.5% gel

Stain O/N -



~~3/2/95~~

3/2/95

~~DESTAIN~~

DESTAIN

Take Picture



Nothing induced -

Page

60

New Clono - HTA Screen

3/7/95

Spin through G25 Column.
Count

HTAI33	41	1	899759.00	0.21	1.00
HTHCM40	42	1	965819.00	0.20	1.00
HT40B94	43	1	1254146.25	0.20	0.80
HTABG94	44	1	1446598.62	0.20	0.70

Add 100ul of Salmon Sperm DNA
Heat 100°C 5min
Quick Chill.

Pour off prehyb.
Add 1500ul Hyb Buffer
(2x PIPES, 10% Dextran Sulfate, 50% Form.)
Add probe to hyb buffer
Incubate 42°C O/N

3/8/95

Wash filters.
Pour off hyb
Rinse filters 0.2xSSC / 0.1% SDS.
Wash filters 65°C
0.2xSSC / 0.1% SDS.
Wash 3x
Put on film
- 1-20A - HTAI33
- 21-40B - HTABG94
The other 3 leave washing at 65°C.
Not enough cassettes.

3/9/95

Develop film
Place remaining filters on film
develop 10 min at 18°C
No Clono 10 min at 18°C HTABG94

KC8V / HPRN08 51bp + PD10

61

(pg 42)

(pg 50)

3/1/95

Add 0.7 volumes (sopropanol) 10.5ml.
 Mix well.
 Spin 8K 30min
 Wash pellet 15ml 80% Ethanol
 -20°C.
 Spin 8K 10min
 Pour off allow pellet to Dry at RT
 20min.

Resuspend pellet in total of 400ul
 TE and transfer to eppendorf
 tube.

Read OD_{260/280} - 1.200 Dilution.

Sample ID	abs 260.0 nm	abs 280.0 nm	bkg abs 320.0 nm	260.0 nm 280.0 nm	280.0 nm 260.0 nm	
1 KC8V PD10	0.1502	0.0995	0.0219	1.6537	0.6047	1.5ug/ml
2 HPRN08 51bp PD10	0.1091	0.0715	0.0144	1.6582	0.6031	1.1ug/ml

Run 2ul on gel w/ 1kb ladder

Plasmid looks good
 Strong plasmid #2.



abs 260.0 nm	abs 280.0 nm	bkg abs 320.0 nm	260.0 nm 280.0 nm
0.0942	0.0638	0.0201	1.6947

0.94ug/ml

3/8/95

Start Culture to do any induction of
 HPRN08 51bp + PD10

Inoculate 30ml LB Amp Kan with
 C-7
 incubate 37°C ON

62

HIPANOS 51bp PD1D -

3/9/95

inoculate 300ml LB+Amp/Km
 10ml of O/N culture of
 HIPANOS 51bp + PD1D @ C50
 incubate 37°C w/aeration
 until OD₆₀₀ ~ 0.4-0.6
 Add 100mM IPTG to 2mM (10ml)
 incubate 37°C 4 hours.
 Spin 8K 20min
 Pour off supernatant
 Resuspend pellet in 30ml 1M Gm HCl
 pH 8
 Store 4°C O/N.

3/10/95

Spin Culture 8K 20min
 Transfer supernatant to fresh tube.
 Run on gel.

400ul H₂O
 20ul Protein in Gm HCl
 50ul 0.15% Na DOK
 75ul 50% TCA.

Mix well


Spin 10 min

Remove supernatant

Resuspend pellet in 10ul 0.2N NaOH

Run on 8.5% gel.

Stain DeStain



looks good -
 - Load Pure
 Column

(PD) 21

pg 60

HT4 Screen / HTA Screen

69

3/10/94

Develop film for

21A-40A HT4CM40
1-20B HT4CB44

3/13/94

Pick positive clones of

HT4A133

HT4CM40

HT4CB44

into 20ul SM in 96 well dish

Pick 48 from each clone
Incubate samples at RT O/N

Plate HTA for screening of
HTA BF 94

Dilute 1/1000 up 20ul into
1ml 1E392 cells $OD_{600} = 1.0$
Incubate $37^{\circ}C$ 15 min
Plate into 150mm NZM Plates con
Final LB + 0.75% Amp. let set
Incubate $37^{\circ}C$ 5 hours
Store $4^{\circ}C$ O/N

3/15/95

for the 48 clones in SM of

HT4A133

HT4CM40

HT4CB44

To 50ul 1E392 $OD_{600} = 1.0$ - add 20ul
of phage
Incubate $37^{\circ}C$ 15 min
Add 15ul NZM Broth

pg 6Z

3/13/95

Prepare NiSO₄ Column.

2 ml Resin

Wash 30 ml H₂O

Strip 20 ml 6M Gn HCl pH 2.

Wash 30 ml H₂O

Change 30 ml 0.1M NiSO₄

Wash 30 ml H₂O

Equilibrate 30 ml 6M Gn HCl pH 8.

Apply supernatant - Collect flow

Wash 30 ml pH 8 6M Gn HCl - Collect


Wash 30 ml pH 6 6M Gn HCl - Collect pH 6

Elute 6 ml 6M Gn HCl pH 5 - Collect

Strip 30 ml 6M Gn HCl pH 2 - Collect

3/14/95

Run 20 µl on gel with uninduced induced cultures.



Protein looks good.

HTPAN08 51bp + PD10

3/15/95

prepare 12.5% Papanicolaou gel
15 mm

72


HTRANDS 88 5167 +PD10

3/15/95

Precipitate 80 μ l of pH 5 in H_2O / 50% KA
 Add 6.15% NaDOC
 Add 400 μ l 0.5N NaOH
 Add 400 μ l 2X Dissociation Buffer
 Heat 100°C 5 min
 Spin 2 min
 Run on gel, 100V.

Cut off Marker a Port of gel.

Spin DESTAIN



Align w / corresponding gel.

Cut out from unstained gel.

Place in 15ml conical
Ready for AT Production


133

HT4/HTA screen

75

3/16/95



HT4CB44



HT4CM40

HT4CM40

HT4CB44



HT4CB44

HT4CB44

Wash HTA & HT4 filters
0.2xSSC / 0.1% SDS. - 3x 65°C
Put on film
-80°C O/N.

3/17/95

Develop film.

He M. Puh
3/17/95

Random Prime Probe - HSEN37, HTUSB02, HSKB09

4/6/95

Mix by Flicking
Quick Spin
Incubate 37°C 10 min

For HSEN37 use Primer if

Primers	10ul
DNA	4ul
H ₂ O	30ul
	34ul

Heat 100°C 5min

Quick Chill

Quick Spin

Add 10ul 5x dCTP Buffer
5ul d₃₂P dCTP
1ul Klenow

Incubate 37°C 10min

Put through G 25 Column

Count probes

* Did Not put ~~HTUSB02~~ ~~HSKB09~~ HSKB09
HTUSB02 through Column.

CPM	2516%	
2391106.50	0.19	HSKB09
3005574.25	0.19	HTUSB02
1012483.00	0.20	HSEN37

Add 100ul Salmon Sperm DNA

Heat 70°C 5min

Quick Chill

N. M. R. h
4/6/95

Screen HOAAH/HSEEN (H5KBN)/HT4A4/HT4CB

4/27/95

Inoculate 30 ml TB+Amp
with HTAB 94, S01, S02, S03, S04, S05 & S06
Incubate 37°C O/N.

From plated resus - pick 6 white
Colonies into each 200 µl LB+Amp
Incubate 37°C 4 hours
Do PCR

HOAAH02		HSEEN37		H5KBN09	
FP50	1	FP50	1	FP50	2
M13R	0.05	M13R	0.05	M13R	0.05
10x	3.2		3.2		3.2
10x	3.2		3.2		3.2
H ₂ O	22.4		22.4		21.4
Taq	0.15		0.15		0.15
Cult	2		2		2
	<u>32</u>		<u>32</u>		<u>32</u>

PCR Profile

HT4A480	
FP01	1.2
M13R	0.05
10x	3.2
10x	3.2
H ₂ O	22.2
Taq	0.15
Cult	2
	<u>32</u>

HT4CB44	
FP02	1
	0.05
	3.2
	3.2
	22.4
	0.15
	2
	<u>32</u>



95°C 5min	
30X	95°C 20sec
	55°C 20sec
	72°C 1min
	72°C 7 1/2 min
4°C Hold	

Alan R. R.
4/27/95

TNT - HTPAN08 51bp Protein Prep

183

5/5/95



1972

5/8/95

Inoculate LB + Amp Kan with
HTPAN08 51bp A74 in PD10
into 100ml.
Incubate 37°C O/N.

Inoculate 250ml LB + Amp with
~~HTPAN08~~ H.TABG 94506.
for May Prep.

Inoculate 30ml LB + Amp with
HOAATH04 5002

Incubate 37°C O/N w/ aeration

5/9/95

Inoculate 1 lot LB + Amp Kan
with 50ml O/N culture
of HTPAN08 51bp A74 in PD10

Incubate 37°C until OD₆₀₀ = 0.4-0.6
Add 100mM IPTG to 2mM → 20ml

134

Maxi HTABG94506 | Medi HOAAH62502

5/9/95

Incubate 37°C 4 1/2 hours

Spin 5K 20 min

Pour off Supernatant
Resuspend pellet in a total of 100 ml
0.1M Gm HCl pH 8
let sit O/N at 4°COxygen Maxi Prep
by HTABG94506
1:200 Dilution

abs	abs	bkg abs	abs	abs
260.0 nm	280.0 nm	320.0 nm	260.0 nm	280.0 nm
0.1548	0.1071	0.0400	1.7092	0.5851

1.55 ug/ml

Run 0.5 ug on gel

HOAAH62502 Alkylin Lysis #1

Spin Culture

Pour off Supernatant

Resuspend pellet in P1 + RNase
(from Qiagen) - 2 ml

let sit RT 5 min

Add 2 ml P2

mix gently

let sit RT 10 min

Add 2 ml P3

mix well

let sit on ice 10 min

Spin 20 min 8K

Transfer Supernatant to fresh tube

Add isopropanol 0.7 Volumes - 4.2 ml

mix well

Spin 30 min 8K

HTRAN08504 51bp ATG + PD10

135

pour off Supernatant
wash pellet 70% Ethanol
Spin 10 min 8K
Allow pellet to Dry at RT 2W

5/9/95

5/10/95

H0AANK02502

Resuspend pellet in 1ml TE
Transfer to 2 microfuge tubes
Add equal volume 13% PEG/1.6M KCl
Mix well
Spin 15 min
1X 70% Ethanol Wash
Resuspend pellet in a total of 200ul of TE
Run 1ul on gel -
Read OD_{260/280} 1:200 Dilution

OD₂₆₀ OD₂₈₀ OD₃₂₀ 2nd 1st 2nd 1st

0.0279

0.0168

0.0034

1.8322

0.5458

0.26ug/ul.

HTRAN08 51bp ATG + PD10.

Spin 8K 20 min
Transfer Supernatant to fresh tube.
Prepare NiSO₄ Column.

Pour 3ml Resin into Column.
Wash 20ml H₂O.

Equilibrate 20ml 6M Gn HCl pH8
Pour on Supernatant. (crude extract)
and Collect Flow through.

Wash Column 60ml 6M Gn HCl pH8

Collect - pH8

Wash Column - 60ml 6M Gn HCl pH6

Collect pH6

Elute BSA Protein - 5ml 6M Gn HCl pH5

Collect pH5

136

HTPAN08504 51bpATG + PD10

5/10/95

Strip Column 30ml 6M Gm HCl pH 2.
collect pH 2.

5/11/95

Run Protein Samples on 12.5% Gel
To Samples in 10 M Gm HCl450 μ l H₂O30 μ l of Sample50 μ l 0.15% Na DOK35 μ l 50% TCA

mix well


Spin 10 min

Remove Supernatant

Resuspend in 10 μ l 0.2 N NaOHAdd 10 μ l 2x Dissociation Buffer

Heat @ 5 min 100°C

Spin

Run Samples on @ SDS PAGE Stacking
gel. 150 V 1 1/2 hours

1 - Unduced
 3 - Crude Extract
 4 - flow through
 5 - pH 8
 6 - pH 6
 7 - pH 5
 2 - Rainbow Marker

Does not look good try reprecipitating
 culture and do another
 induction

~~HTPAN08504~~ 51bpATG + PD10

Fragment prep. - Northern Blots

137

5/17/95


Set up Digests of Unknowns to give
to Brent Kieder - He'll do
Northern.

Clone ID	Conc	5ug DNA	H ₂ O	10X ¹ 2	Xho I	Eco RI
HSKBNO9.	2.06 µg/µl	2.5	42.1	5 µl	0.2	0.2
HNBAAZ6	0.64 µg/µl	7.5	37.1		0.2	0.2
HILBY30	0.5 µg/µl	10	34.6			
HT4A133	PCR Product.	20	24.6			
HT4CM40	0.49 µg/µl	11	33.6			
HT4CB44	1.2 µg/µl	4.2	40.4			
HNFAAG4	0.73 µg/µl	6.8	37.8			
HT4A180	0.34 µg/µl	15.6	29			
HTABG894	0.97 µg/µl	5.2	38.4			

Digest 37°C O/N.

5/12/95

Run 5 µl on gel with 1 kb marker




- Xho I/Eco RI Digests
- 1 - HSKBNO9.
 - 2 - HNBAAZ6
 - 3 - HILBY30
 - 4 - HT4A133
 - 5 - HT4CM40
 - 6 - HT4CB44
 - 7 - HNFAAG4
 - 8 - HT4A180
 - 9 - HTABG894

North
5/12/95

Fragment preps - Northern Blots

5/12/95

Run on 0.8% LMP Gel with 1 kb ladder.
 Run 80V 2 1/2 hours.



Cut out fragments
 - place into 1.5 ml
 microfuge tube
 Store 4°C over
 weekend.

Need to digest
 HILBY30 Again.


201

- H TRANS 5.1 kb ATG in PDI
- INDUCTION
 - Resuspend in 40 µl ml COM Gm HLP/H8
 - Store 4°C over weekend.

5/15/95

Gene Clean fragments
 Resuspend in 3 µl TE.

Run back on gel with 1 kb ladder



- 1 H5KBA009
- 2 HNBAA26
- 3 HT4A133
- 4 HTUCM00
- 5 HTUCB44
- 6 HNFPA66
- 7 HT4A180
- 8 HT4B694
- 9

Send to Brent Knieder - Kind Kovacs

5/15/95

Clone ID	Libraries Expressed	Size in Kb RI/Xho I	Eco	Approx. Conc(ng/ul)	
HSKBN09	HT4, HSK	~1.4	AMK	~100	15ul
HNBA26	HBJ, HBM, HCA, HNB, HRG, HTO	~0.800		~100	↓
HILBY30	HIL, HLM	~0.60		Will Send 5ul	
HT4AI33	HT4	~0.90		~150	15ul
HT4CM40	HT4	~1.7		~150	
HT4CB44	HT4	~2.5		~150	
HNFAA64	HNF, HSI, HSK, HTA	~1.80		~150	
HT4AY80	HT4, HTX, HT3	~0.85		~50	
HTABG94	HTA	~1.7		~150	
HT4CI56	HT4	~1.7		~100	
HT4AI55	HT4	~1.7		~140	
HT4CL32	HT4, HT5, HT3	~1.1		~150	
HT4CA46	HT4, HT3, HCE, HGO, HTA, HL3, HMW	~1.7		~250	↓
HMSAF22	HMS, HOS, HHF, HSR HTN	~2.0		~100	25ul
HT3SB70	HNF, HT3, HT4, HT5, HTA	~1.5	CLF	~300	5ul
HT4SB02	HET, HGL, HHF, HSU, HT4, HTE, HTP	~1.3	AMK	~200	10ul

Set up Digest of HILBY30 PCR product

DNA	20ul
10x B2	5
H ₂ O	24.6
EcoRI	0.2
Xho I	0.2
	<hr/> 30ul

Digest at 37°C
Overnight

HIPAN08 51bp ATG + PD10
Put over column
Collect Flow through
pH 8
pH 6
pH 5

140


HTPANOS 504 516p A14 + PD18

(my 130)

Run

gel

HTPANOS Samples on 12% Acrylamide




- 1 Rainbow Marker
- 2 Crude Extract
- 3 flow through
- 4 pH8
- 5 pH6
- 6 pH5
- 7 pH2.

Store at 4°C

Carrie Did Gene clean of HILBY30
Xho RI

Run full on gel with 1Kb ladder

Give to Brent
~600 bp
~100 ng full.

HIPANOS8504 51bp ATG +PO10

141

5/18/95

Give

Reapply pH5 of HIPANOS851bpATG to
fresh Column. (from 3/31/95 - pg 711)
Add 3ml pH5 + 2ml pH8.

Wash 3ml pH8.

Wash 3ml pH5.

Give to protein expression for
renaturation over column.

5/19 & 5/22

Computer


5/23/95

Received Column Back

Elute 2ml Immidazole elution Buffer.

2 times - Run on 12% Gel

with Marker + HIPANOS851bpATG



Computer

off.

5/23, 5/24
5/25

pg 28
Prot # 215
DNA # 10

5/26/95

HIPAND8504 5bp ATG +AD10

141

5/18/95

Good

Reapply pH5 of HIPAND8516pATG to
fresh Column. (from 3/31/95 - pg 31)

Add 3ml pH5 + 2ml pH8.

Wash 3ml pH8.

Wash 3ml pH5.

Give to protein expression for
renaturation over column.


5/19-5/22

Computer

5/23/95

Received Column Back

Elute 2ml Immobilized elution Buffer.
2 times - Run in 12% Gel
with Marker + HIPAND8516pATG



Computer

off.

5/23, 5/24
5/25

pg 28
Book #115
DMK #10

5/26/95

**This Page is Inserted by IFW Indexing and Scanning
Operations and is not part of the Official Record**

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

- ☐ **BLACK BORDERS**
- ☐ **IMAGE CUT OFF AT TOP, BOTTOM OR SIDES**
- ☐ **FADED TEXT OR DRAWING**
- ☐ **BLURRED OR ILLEGIBLE TEXT OR DRAWING**
- ☐ **SKEWED/SLANTED IMAGES**
- ☐ **COLOR OR BLACK AND WHITE PHOTOGRAPHS**
- ☐ **GRAY SCALE DOCUMENTS**
- ☐ **LINES OR MARKS ON ORIGINAL DOCUMENT**
- ☒ **REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY**
- ☐ **OTHER:** _____

IMAGES ARE BEST AVAILABLE COPY.

As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.