





Department MOL. BIOL. - PROT. EXPR.  
Subject VE/95 - 6/5/95  
Name ANN KIM # 9  
Address \_\_\_\_\_  
43-648  
 National Brand  
**Computation Notebook**  
Dennison Stationery Products Co., Framingham, MA 01701  
75 Sheets  
11" x 9"  
4x4 Quad.  
  
0 73333 43648

BEST AVAILABLE COPY

Ruben EXHIBIT #92

Department MOL. BIOL - PROT. EXPR.  
Subject 12/95 - 6/5/95  
Name ANN KIM # 9  
Address \_\_\_\_\_  
 43-648  
**Computation Notebook**  
Dennison Stationery Products Co., Framingham, MA 01701  
 75 Sheets  
11 1/4" x 9"  
4x4 Quad.  
0 73333 43648 4

Ruben EXHIBIT 2092  
Ruben v. Wiley et al.  
Interference No. 105,077  
RX:2092

1410502 19500 100000

pg 15 | Book 8 #236

1/24/95

Spin HIPANOS 8.85 bp + PDEGO 1u.  
with Hcl pH 8 - 8K 20min.

Equilibrium N. 504 column with pH8  
Cation HCl

Apply Supernatant to Column - Collect Flow  
 Wash 30ml pH 8.0 6M HCl - Collect pH 8  
 Wash 30ml pH 1.0 6M HCl - Collect pH 6  
 Elute 5ml pH 5 6M HCl - Collect pH 5  
 Strip 30ml pH 2 6M HCl - Collect pH 2

add 50ml of ~~distilled~~ collected material  
to 450ml H<sub>2</sub>O  
50ml 80% C. 15% NaDOC  
75ml 50% TCA.

mil wheel

Spinn 5 m

Remove supernatant

Resuspend pellet in 15ml 0.2M NaOH / 15ml 2x Buff

Net 100°C 5 min

Pen 20 und magel w. in Rahmen Mark

100% 10.5% gel fraction at 150V 1 1/2 hrs.

[illegible]

De Staat 1 1/2 hr.

1		Unbuffered
2	1	Crude Extract
3	2	flow
4	B <sub>1</sub>	pH 8
5	B <sub>2</sub>	pH 6
6	B <sub>3</sub>	pH 5
7	B <sub>4</sub>	pH 2

$pH \approx 2$        $TDS$   
 Filter      & Send some  
 Dialyze      to protein exp  
 do protein exp      to measure  $OD_{600}$   
                                  column

10260

Handwritten signature: *Handwritten signature*

11-6 PQEGO

HTAN08185 bot PQEGO

(pg H8 Box 8 #236)

2/2/95

11-6 - Reapply to column + Strip again  
do try to Purify more.  
Start Dialyzing ~~at 20 ml~~ in Dialysis  
Tubing

3M 6N HCl / Hepes 5 hrs  
1.5 M 6N HCl / Hepes Over night  
Reapply 4.3 ml to Ni Sep Column to Reattach  
over column of HTAN08185 bot PQEGO.

2/3/95

Change Buffer:

1 M 6N HCl / Hepes 10 hrs  
0.5 M 6N HCl / Hepes Over Weekend

Carrie will finish

(2/6 - 2/10 Vacation)

HTAN08185 bot + PQEGO

2/3/95

Strip Column that has cracked  
material - Protein in 2  
strips in imidazole Elution  
Buffer: 50 mM NaPO<sub>4</sub> pH 6  
250 mM Imidazole  
300 mM NaCl  
10% Glycerol

2 strips at 2.5 ml each

Run on Stacking gel with 1 M HCl  
Marker - 0.125% gel 1.50 V

Stain 1/2 hr  
de stain 1 hr


4

1L6 PQE60 / HTPANOS 185bp + PQE60


2/3/95

HTPANOS 185bp + PQE60

Store 4°C.



1L6 PQE60 + PQE60 2nd Column Extraction



43

29

184

H1PAN08/H1PB411 - PD10 - PAZ

5

pg 152 Book 8 # 236

1/26/95

Received Primers for  
H1PAN08 & H1PB411  
Reprotect 55°C O/N.

1/27/95

PCR Fragments

H1PB411 PAZ

9130 5' Bam.

3' Xba New.

10x dNTP

10x PCR

Tag

H<sub>2</sub>O

DNA (10 µg/µl)

1  
1  
50  
50  
2  
397  
1

500 µl - 100 µl / reaction

H1PAN08

3146

9111:

10x dNTP

10x PCR

Tag

H<sub>2</sub>O

DNA (10 µg/µl)

51 PAZ

2  
20  
50  
50  
2  
375  
1  
500

185 PAZ

2  
9112: 20  
50  
50  
2  
375  
1  
500

51 PD10

2  
9113: 24  
50  
50  
2  
371  
1  
500

185 PD10

2  
9114: 20  
50  
50  
2  
375  
1  
500

100 µl / rxn

HTPB411 / HTPB411 VDB / PAZ

1/27/95

Run PCR:

95°C	5 min	] 25X
95°C	30 sec	
55°C	30 sec	
72°C	1 min	
72°C	7 1/2 min	
72°C	Hold	

Run 5ul of Run on gel with 1 kb ladder



1	HTPB411	5.0um/3.0um	PAZ
2	HTPB411	51bp	PAZ
3		185bp	PAZ
4	↓	51bp	PAZ
5		185bp	PAZ

Precipitate Reactions

Add equal Vol

13% PEG-NaCl

Spin 10min

Remove Supernatant

1000ul 70% Ethanol

Wash pellet

Digo Remove Supernatant

Dry pellet 5min at RT

Resuspend pellet in 100ul TE

Set up Digests

DNA (PCR Run)

10X #2 Buffer

H<sub>2</sub>O

Bam

Xba

10 ul

5

34

0.5

0.5

50ul

Incubate 37°C 4hrs


Run on 0.8% LMP gel with 1 kb ladder

HTP808 / HTPB11

PD10 / PAZ

7

1/27/95

Cut out bands  
Take picture:

1	HTP808	51 bp	PAZ
2	↓	185 bp	PAZ
3		51 bp	PD10
4		185 bp	PD10
5	HTPB11	PAZ	

Gene Clean fragments

- Resuspended in 40ul TE

Set up ligations.

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17
9111 + 3146	6	6							6	6							
9112 + 3146			6	6						6							
9113 + 3146					6						6						
9114 + 3146						6						6					
HTPB11 (PAZ)							6	6					6				
10X Buffer	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2
H <sub>2</sub> O	9	9	9	9	9	9	9	9	11	11	11	11	11	15	15	15	17
T4 Lig (14ul)	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
PD10 8/x 9/16					2	2								2			
PAZ 8/x 9/16	2		2				2								2		
PAZ 2 1/2 1/27		2		2				2								2	

 10X Buffer  
 T4 ligase  
 H<sub>2</sub>O

 2  
 1  
 9  
 12ul

 18X  
 183g  
 18  
 162  
 12ul / Tube

 Add Appropriate  
 to Vector Fragment  
 or H<sub>2</sub>O



Hunt 11/18/64

VALUATION

$$1/27/95$$

Membrane ligations 16% over weekend

At A/Ry  
127/95

1/30/95

Transform ligotensis

100ul of Chem. Competent Cells

10. ul of degeneration Rhen

PA2 Constructs DH5 $\alpha$

P.D10 Constructs M.15 rep4.  $\oplus$  control P.D10

Incubate on ice 1 hr

Heat  $42^{\circ}\text{C}$  45 sec

Place on ice

Add 400 ml LB

Uncultivate 37°C off / hr

Plate 200ul + 300ul onto

Lb + Amp plates for all lysates

1-8 onto 150mm plate

for, heptanoic 9-17 plaster 100ml

onto 2.8t Amp 100min Plates

Incubate 37°C O/N

1/31/95

Plates look good

No colonies in control states

Colonies on (+) Control plates

- Inoculate plasmid colonies into LB + Amp for pA<sup>+</sup> constructs - 200  $\mu$ l of; 100  $\mu$ l in 96 well dish.

Inoculate colonies into LB + Amp Kan  
for PAID constructs

200ml of LB + Amp in 96 well plate

HTPB411 / 1/17/95 PAZ / ~~PAZ~~ PD10

9

1/31/95

B+amp:

① 9111/3146 + PAZ 1/6 48  
② 9111/3146 + PAZ 1/27 13  
③ 9112/3146 + PAZ 1/6 35 } Plate # 1

③ 9112/3146 + PAZ 9/6 12  
④ 9112/3146 + PAZ 9/27 13  
PAZ 1/6 2  
PAZ 1/27 1  
⑤ HTPB411 + PAZ 9/6 48  
⑥ HTPB411 + PAZ 1/27 12 } Plate # 2

LB+amp/Ran:

⑤ 9113/3146 + PD10 48  
⑥ 9114/3146 + PD10 48 } Plate # 3

Incubate plate 4 hrs 37°C. with aeration  
Setup PCR's.

9111	1	70x
3146	0.1	70
10xPCR	3.2	7
10x dNTP	3.2	240
H <sub>2</sub> O	22.35	240
Tap	0.15	1504.5
Cult.	2	10.5
	<u>32</u>	30ul/tube

9112	1	55x
3146	0.1	55
10xPCR	3.2	5.5
10x dNTP	3.2	176
H <sub>2</sub> O	22.35	176
Tap	0.15	1229.25
H <sub>2</sub> O	2	8.25
	<u>32</u>	30ul/tube

9113	1.0	(50x)
3146	0.1	50
10xPCR	3.2	5
10x dNTP	3.2	160
H <sub>2</sub> O	22.35	160
Tap	0.15	117.5
Cult.	2	7.5
	<u>32</u>	32ul/tube

9114	1	(50x)
3146	0.1	50
10xPCR	3.2	5
10x dNTP	3.2	160
H <sub>2</sub> O	22.35	160
Tap	0.15	117.5
Cult.	2	7.5
	<u>32</u>	30ul/tube

10

HTPB411 (HTPB411)

PD10 / ~~PD10~~ 0x2

1/31/95

HTPB411 5' Bam  
3' Xba

10x dNTP

10x PCR

H<sub>2</sub>O

Taq

Cult.

0.1

0.1

3.2

3.2

23.25

0.15

2

32ul

65X

6.5

6.5

208

208

151.25

9.75

32ul / tube.

PCR.

95°C

95°C

55°C

72°C

72°C

4°C

5min

20sec

20sec

1min

7 1/2 min

Hold.

30X


PA2 controls  
for 9111, 9112  
4 (HTPB411)

2/1/95

Run Reactions on 1% TAE Agarose  
gel with 1kb ladder.

9113 A1-D12

9113/9114 E1-H12




HTPB02 / HTPB11 + PD10/PAR

11


911A  
FI-H12

911B  
AI-F1




9111  
AI-F1

21/195




9111  
AI-F1

9112  
F2-H12




9112  
F2-H12 AI-C1

PAR Controls HTPB11  
DT



HTPB11  
DT-H12



HTPB11  
DT-H12



12 HTPA008 / HTPB411 PA2 / AD10

1/1/95

Inoculate 200 ul LB + Amp/Kan  
with 20 ul of (+) clones to do  
Small scale inductions (micro)  
incubate at 37°C w/ aeration  
2 hrs.

Add 3 ul of 100 mM IPTG to 20 mM IPTG  
final conc.

Incubate 4 hrs 37°C w/ aeration

Spin 10 min

Resuspension pellet in 5 ul H<sub>2</sub>O

Add 5 ul 2x Buffer

Store -20°C till tomorrow.

Inoculate 5 ml TB + Amp  
with PA2 (+) clones.

- |     |   |    |                      |
|-----|---|----|----------------------|
| (1) | - | 11 | 9N1 / 3128 + PA2 3/4 |
| (2) | - | 8  | 9N1 / 3146 + PA2 1/2 |
| (3) | - | 10 | 9N2 / 3146 + PA2 3/4 |
| (4) | - | 10 | 9N2 / 3146 + PA2 1/2 |
| (7) | - | 1  | HTPB411 + PA2 3/4    |

Incubate 37°C w/ aeration o/n.

2/2/95

Run Protein gels 15% Stacking  
15 ul of Sample + 1 MLE marker.

150V 1 1/2 hrs

Stain 30 min 37°C

Destain 1 hr 37°C


HTPA1508/HTPB411

PA2/PD10

18

HTPA1085 bp + PD10


2/2/95



HTPA1085 bp

Should produce  
a protein:  
32.5 kD

HTPA1085 51 bp + PD10




HTPA1085 18.5 bp

Should produce  
a protein: 27.7 kD

HTPA1085  
18.5 bp + PD10

HTPA1085  
18.5 bp + PD10




looks like HTPA1085  
18.5 bp + PD10 -

does not  
have induction

14

HTPANO8 | HTPB411 PDIO / PAZ

2/2/95



looks like  
all induced  
except  $\odot$  Control.

Do large scale  
prep.

HTPANO8SC1 (SSD) + PDIO

~~Obv~~

Do Boiling Mini Preps. of PAZ Constructs

2ml culture Spin 2min

Remove Supernatant

Resuspend pellet in 75ul STEC +

RNase + Hypoxanthine

Heat 100°C 2min

Spin 10min

Remove pellet

Add 75ul 13% PBG 8000 / 1.64 NaCl

mix well

Spin 10min

Remove Supernatant

Add 75ul 70% EtOH to wash pellet

Spin 5min

Remove Supernatant

Allow Pellet to dry at RT 10min

Resuspend pellets in 150ul TE


Run 2ul on 1% TBE gel with 1 kb

adder

HTPANO81 HTPBY11

PD10/PAZL 15

2/2/95



Mini preps all  
look good


Check for inserts.  
By digesting w/  
Bam / Xba!

DNA	5 $\mu$ l	}
OK	30 $\mu$ l	
H <sub>2</sub> O	21.8 $\mu$ l	
Bam	0.1 $\mu$ l	
Xba	0.1 $\mu$ l	
	30 $\mu$ l	

Incubate 37°C  
overnight

2/3/95

Run 10  $\mu$ l on gel w/ 1 Kb ladder.



All look correct  
- best  
Select 2 40 Sequences  
with internal  
primers.

(1)	1, 2	}	RPO6 FPI6
(2)	1, 2		
(3)	1, 2		
(4)	1, 2		

HTPBY11 + PAZL - FPR3

Carve Cheval Seq. Some looked good so put them



16

HTRAND8/HTRB411 + PAZ/PD10

1/13/95

- Inoculate 3ml TB+Amp w/  
Cultures 1-1, 1-2, 2-1, 2-2, 3-1, 3-2  
4-1, 4-2, & 1/2 HTRB411 + PAZ.  
Make Glycerols.
- Inoculate 200ml LB+Amp with  
HTRB411 + PAZ to do Maxi Prep.
- Inoculate 10ml LB+Amp Kan with  
induced culture of HTRAND8/185bp+PD10  
(#15) - Do large scale induction

1/14/95

- made Glycerol Stocks  
- 80°C Protein Expression Box #1
- HTRAND8 185bp + PD10 #12 -  
Inoculate 300ml LB+Amp/Kan  
with 300ml of O/N culture  
Inoculate 300ml 3hrs - until  $OD_{600} = 0.4-0.6$   
Add 100mM IPTG to 2mM - 6ml  
Inoculate 300ml 4 1/2 hrs w/aeration  
Spin 5K 15 min  
Resuspend Pellet in 30ml of  
2M NaCl 0.1M  
Store at 4°C O/N.
- HTRB411 + PAZ Maxi Prep
- Reagen Maxi  
Spin Culture 6K 20min  
Pour off Supernatant

(pg 27)

HIRANUS / HTPBY11

27

pg 16

2/14/95

HTPB11 FRAZ

Resuspend pellet in 10ml of P1 HRase  
Let sit RT 5min

Add 10ml P2 while mixing

Add 10ml P3 while mixing

Let sit on ice 20min

Spin 20 min 8K

Equilibrate Tip-500 water 10ml QBT

Apply Supernatant to Equilibrated

Column

Wash Column 2x with 30ml

P2

Elate DNA in 15ml QF

Add 0.7x (10.5ml) of Isopropanol

Mix Well

Spin 9K for 25min

Pour off Supernatant

Wash pellet with ice cold 70% Ethanol

(10ml)

Spin 9K 10min

Pour off 70% Ethanol

Allow pellet to dry at RT

Resuspend pellet in 100ul TE

Run Blue on gel

Read OD 260/280 at 1:200 Dilution

abs	abs	bkg abs	260.0 nm	280.0 nm
260.0 nm	280.0 nm	320.0 nm	280.0 nm	260.0 nm
-0.0063	-0.0028	-0.0022	6.9003	0.1449
0.1519	0.0964	0.0232	1.7591	0.5885


1.52 ug/ul + 600ug  
Total

Store 4°C Plasmid Box #2  
Sequence w/ internal primers to confirm sequence

28

HTPB411 &amp; HTPAN08

2/14/95



looks good.  
- See if sequence is good

2/15/95

Digest @ 100 DNA w/ 100 Bam/Klu  
to see if closest "Pop" cut

DNA (250mg)	4
10X #2	3
H <sub>2</sub> O	22.6
Bam	0.2
Klu	0.2
	<hr/> 30

Incubate 37°C

Submit for sequencing w/ internal  
primers. HTPB411 PA2 RA/FP

RPB50A  
RPB50B  
RPB50C  
RPB50D  
RPB50E  
RPB50F  
RPB50G  
RPB50H  
RPB50I  
RPB50J  
RPB50K  
RPB50L  
RPB50M  
RPB50N  
RPB50O  
RPB50P  
RPB50Q  
RPB50R  
RPB50S  
RPB50T  
RPB50U  
RPB50V  
RPB50W  
RPB50X  
RPB50Y  
RPB50Z

RPD1A  
RPD3A  
RPD4A  
RPD5

RPD6A  
RPD7  
RPD8  
RPD9

RPD10  
RPD11  
PA2  
RPD13

FP4  
FP5  
FP10  
FP17

FP18  
FP19  
FP20  
FP21A

FP22A  
FP23C  
RPD50

HTPAN08/HTPBYN

29

Submit for Sequencing w/ internal

2/15/95

HTPAN08 51bp + PAZ

HTPAN08 185bp + PAZ

Empty Plasmids No Stress to submit  
to Protein Expression for baculovirus.

51bp.

RP12	FP14	RP01
FP13	RP05	FP18
RP10	RP06	RP08
FP16	FP17	RP50

185bp

FP16	FP17	RP50
FP14	RP01	
RP05	FP18	
FP08	RP06	

HTPAN08PA5 RP/FP

HTPAN08P18 SRP/FA

HTPAN08 185bp + PD10 #12 large scale

~~Inductions~~

Spin Culture 20 min 8K

Prepare NiSO<sub>4</sub> Column

Wash Apply 2ml Resin to Column

Wash 20ml H<sub>2</sub>OAdd 30ml 0.1M NiSO<sub>4</sub> for ChangeWash 30ml H<sub>2</sub>O

Equilibrate with 30ml 0.1M HCl pH 6

Apply Supernatant - Collect Flow

Wash 45ml pH 8 - Collect pH 8

Wash 45ml pH 6 - Collect pH 6

Elute 5ml pH 5 - Collect pH 5

Strip 45ml pH 2 - Collect pH 2

Run on 15% Acrylamide Stacking gel.

190ul H<sub>2</sub>O

20ul 10% Gm Hcl prep-

75ul 50% TCA

50ul 0.15M NaOAc

mix Well

Spin 10 min


30 HTPB411 & HTPB0508

2/15/95

Remove Supernatant  
Resuspend pellet 10  $\mu$ l 0.2N NaOH  
Add 1  $\mu$ l 2X Dissociation Buffer  
Boil 5 min - lost pH 8  $\rightarrow$  pH 2 = Samples  
mixed with water  
Spin 5 min  
Run 20  $\mu$ l on gel with LMW  
Protein marked  
48 180V 1 hour  
Stain 30 min. 37°C  
DESTAIN overnight.

2/16/95

HTPB08 185 PNO #12



Protein  
looks good  
Reapply 5ml of  
pH 5.0 to Column  
Add 30  $\mu$ l pH 8  
Send to Protein  
Expression &  
have some time  
over column

Inoculate 500  $\mu$ l of 1 LB + Amp/Kan  
with 20  $\mu$ l of HTPB411 + HTPB08

H1P6III + H1P6N08



31

2/16/95

Incubate 37°C w/ aeration  
2 hrs.  
Add 100 mM IPTG to 2 mM  
10 $\mu$ M.  
Incubate 37°C w/ aeration  
over night.

2/17/95

Spin Cultures 2 min  
Remove Supernatant.  
Resuspend culture 200 $\mu$ l H<sub>2</sub>O.  
Add 20 $\mu$ l 2X Densocentr Buffer.  
Heat 100°C 5 min  
Re Spin 5 min.  
Run 10 $\mu$ l on 10% Stacking  
gel.  
Accidentally ran 1 kb Marker instead  
of Rainbow Marker.



H1P6III + P2E60


Run 150V 1/2 hr  
2 min  
DESAIN

16-126  
2/17/95

32

HTP ANDS / HTPB11

2/2/95



HTPB11 + POEG60

Try growing up 1 for induction -  
large scale

2/21/95

Incubate 5ml LB + Amp/Kan  
with HTPB11 + POEG60.  
2-2 + 3-4  
incubate 37°C

Transform - HTPB11 5/6p + PD10 7/13/95  
HTPB11 5/6p + POEG60 7/5/95  
into M15 Chemically Competent Cells

Thaw M15 on ice.  
To 100  $\mu$ l of Cells add ligand  
Incubate on ice 1 hour.

(443)

11.1.1.1 + SV  
Pg 26


33

2/17/95

Mix well  
Spin 10 min  
Remove Supernatant  
Resuspend in 10  $\mu$ l 0.2% NaOH  
Add 10  $\mu$ l 2x dissociation Buffer  
Heat 100°C 5 min  
Run all on gel.  
150V 1 hour.

STAIN 30 min 37°C  
DESTAIN 30 min 37°C

Take Picture.



looks like I have  
protein.  
Looks slightly  
contaminated  
Try re-eluting  
to try and clean  
up prep.  
Need to run PAGE

2/21/95

pink 48 more clones from the  
2/18/95 into 200  $\mu$ l of LB amp/kan



H1PAX08 | H1PBY11

43

pg 32

2/2/95

Heat 48°C 45Sec.  
Place on ice  
Add 400ul of LB  
Incubate 37°C 1 hour.  
Plate 300ul onto LB+Amp 150mm plate  
100ul onto LB+Amp/Kan 100mm plate  
Incubate 37°C O/N.

2/22/95

PICK H1PAX08 ST + PD160 48 AEI-H12  
H1PBY08 ST + PD10 48 AI-D12  
into 200ul LB+Amp/Kan  
Incubate 37°C w/aeration O/N.

H1PBY11 + PD160.  
to 300 ml LB+Amp/Kan. add  
3ml of O/N culture. 2-2/3-4  
Incubate w/aeration 3 hrs until  
OD<sub>600</sub> ~ 0.4 - 0.6  
Add 100mM IPTG to 2mM (10ml)  
Incubate 37°C 4 1/2 hours  
Spin cultures 5K 20min  
Remove supernatant  
Resuspend pellet in 30ml CoMGM/HCP  
pH 8.  
Store 4°C O/N.

2/23/95

Do PCR of H1PAX08 clones.

44

HTRANOS6/HTRB411

2/23/95

HTRANOS6 51bp + PD10  
(50x)

q113	1.5	75
BH46	0.1	5
10x dNTP	3.2	160
10x PCR	3.2	160
H <sub>2</sub> O	21.9	1095
Taq	0.1	5
Culture	2	
30ul/tube		

HTRANOS6 51bp + PD10

q113	1.5	50x
PD10	0.1	75
10x dNTP	3.2	5
10x PCR	3.2	160
H <sub>2</sub> O	21.9	160
Taq	0.1	1095
Culture	2	5
30ul/tube		

PCR program # 10C

Run 10ul on gel w/ 1 Kb ladder



HTRANOS PD10 A1-D12  
A3, A7, B12, C5, C7, D5, D9, D10  
Incubate 2ml LB + Amp/Kan  
w/ for 20ul O/N culture

HTRANOS PD10 E1-H12  
E1, E2, F1, F2, G1, G3, H1, H2  
Incubate 2ml LB + Amp/Kan w/ for  
20ul O/N culture  
Incubate at RT Overnight

HTPBVI/HTPBVII

45

2/23/95

HTPBVII + PDE60. 2-2 + ~~3~~ 3-4  
Spin culture 8K 10 min.  
Transfer supernatant to fresh tube.  
(Crude Extract)

Prepare Column.  
Prepare Column with 2ml NTA Resin  
Wash 30ml H<sub>2</sub>O  
Change 30ml 0.1M Na<sub>2</sub>SO<sub>4</sub>  
Wash 30ml H<sub>2</sub>O  
Equilibrate 30ml 6M GnHCl pH8.

Apply Supernatant to Column.  
Collect as flow.  
Wash 30ml 6M GnHCl pH8  
Collect as pH8.  
Wash 30ml 6M GnHCl pH6  
Collect as pH6.  
Elute protein 6.5ml 6M GnHCl pH5  
Collect as pH5  
Strip Column. 30ml 6M GnHCl pH2  
Collect as pH2.  
Store 4°C O/N.


2/24/95

HTPBVII + PDE60.

400ul H<sub>2</sub>O  
20ul of eluted Protein in 6M GnHCl  
50ul 0.5% NaDOC.  
75ul 50% TCA  
Mix well  
Spin 5 min  
Remove Supernatant  
Resuspend pellet in 0.2N NaOH. -10ul

2/24/95

Add 10ul 2X Dissociation Buffer  
 Heat 100°C 5 min  
 Spin 1 ml  
 Run 15ul on gel with Unclonucel  
 and Rainbow marker  
 12.5% Acrylamide Stacking gel  
 150V 1 hour



STAIN 37°C  
 30 min  
 DESTAIN 30 min  
 37°C

Take Picture  
 looks like isolated  
 protein

Store at 4°C  
 Ask Steve about  
 how to get  
 off of membrane

HTRAND8 5kbp + PDE/PD10

Place tubes at 37°C w/ aeration

For 2 hrs  
 Add 100mM IPTG to 2mM 4ul  
 Incubate 37°C 4 hours  
 Spin 1 ml 2 min  
 Remove supernatant  
 Resuspend pellet in 30ul H<sub>2</sub>O  
 Add 30ul 2X Dissociation Buffer  
 Run 15ul on gel with Rainbow Marker  
 150V 1 hour

HTRANS | HTPBY11


47

2/24/95

Stain 30 min 37°C  
DESTAIN over weekend at RT

2/27/95

Take picture



A2  
C7  
D5  
D9


Induced

Show up to  
make glycyls.  
Pulse 1 to Do  
Case Scale  
Prep

C7 - large Scale  
Inoculate 10ml

LB + Amp/Kan w/  
C7. incubate o/n

w/ aeration. Make  
Glycyls of A2, C7, D5, D7  
Does not look  
like anything  
induced - )



Try running  
again - seriously  
10/11 this time

48

HTPAN08 HTP04

2/27/95

Re Run HTPAN08 51bp + PGE60  
10ul.

Run 10ul each of the 1st + 2nd  
imidazole elution of HTPAN08 185 + PD10  
#12

2/28/95



HTPAN08 185bp + PD10 #12


looks like both fractions  
have protein.

~~HTPAN08 185bp + PD10~~

Amplify 300ul LB +  
amp Kan with  
HTPAN08 51bp + PD10 #7  
to do (quick mass)  
Amplify 37°C 8 generations  
O/N.

Re Run HTPAN08 51bp + PGE60

150V 1 hour  
Stain / DESTAIN



Does not look like  
there was any  
induction  
- true over-crowd  
- shut it up!

E 4, 5, 6, 8  
F 3, 4, 5, 6  
G 4, 6, 7, 8  
H 3, 4, 5, 6, 7, 8

Amplify 300ul LB +

HTPAN08 / HTPBULL

49

Wash 30ml of culture from 4/23  
(log phase) O  
let sit at RT O/N

2/28/95

HTPAN08 51bp + PD10 CF (Glycerol)  
Maxi - See Pex 42 of Pex  
along side PD10 to SV PDECO

HTPAN08 51bp + PD10 CF  
Large Scale Induction  
Induce 300ml LP + amp Kan  
5ml of O/N culture  
Incubate 37°C w/ aeration  
until O.D. 600 = 0.4 - 0.6 - 2 1/2 to  
3 hours  
Add 100mM IPTG to 2mM (10ml)  
Incubate 37°C w/ aeration  
4 1/2 hours.  
Spin culture 5K 15 min  
Pour off supernatant.  
Resuspend pellet in 30ml  
Glycine HCl pH 8  
Store 4°C O/N.

3/1/95

Incubate HTPAN08 51bp + PDECO  
(in 2ml LP + amp Kan 37°C  
w/ aeration)  
Incubate until ~ 2 hours.  
Add 100mM IPTG to 2mM (4ml)  
Incubate 37°C 5 hours  
Spin 1ml culture  
Remove supernatant  
Resuspend pellet in 40ul H<sub>2</sub>O  
Add 40ul 2x Dissociation Buffer  
Alec Hly  
3/1/95

50

H18W H18BY

2/1/95

Heat 100°C 5min

Spin 2min

Run 10ml on gel w/ 2 MW/marker

180V

40min



12.5% gel

Stain O/N -

B12125  
3/2/95~~DESTAIN~~

DESTAIN

Take Picture



Nothing unusual -

Paler



60

New Clono - HTLA Screen

3/7/95

Spin through G25 Column.  
Count Lys.

HTAT33	41	1	899759.00	0.21	1.00
HTACM40	42	1	965819.00	0.20	1.00
HTACB94	43	1	1254146.25	0.20	0.80
HTABG94	44	1	1446598.62	0.20	0.70

Add 100 $\mu$ l of Salmon Sperm DNA  
Heat 100°C 5min  
Quick Chill.Purify prep.  
Add 500 $\mu$ l Hyb Buffer  
(2xSSC / 0.1% Dextran Sulfate, 50% Form.)  
Add prep to hyb buffer  
incubate 48 hr O/N

3/8/95

Wash filters.

Purify hyb  
Purify filters 0.2xSSC / 0.1% SDS  
Wash filters 65°C  
0.2xSSC / 0.1% SDS

Wash 3x

Put on film

- 1-20A - HTAT33

- 21-40B - HTABG94

The others leave in tank at 65°C.  
Not enough cassettes.

3/9/95

Develop film

Place remaining filters on film  
expose O/N at 80°C

No cassettes for HTAT33

KC SV / HIPANOS 81bp + PDE100

61

(pg 42)

(pg 50)

3/1/95

Add 0.7 volumes isopropanol 10.5ml  
 Mix well.  
 Spin 8K 30min  
 Wash pellet 15ml 80% Etanol  
 -20°C.  
 Spin 8K 10min  
 Pour off allow pellet to dry at RT  
 20min.

Resuspend pellet in total of 400µl  
 TE and transfer to eppendorf  
 tube.

Read OD<sub>260</sub>/280 - 1.200 Dilution.

Sample ID	abs	abs	bkg abs	260.0 nm	280.0 nm	
	260.0 nm	280.0 nm	320.0 nm	280.0 nm	260.0 nm	
1 KC SV PDE100	0.1502	0.0995	0.0219	1.6537	0.6047	1.5 µg/µl
2 HIPANOS 81bp PDE100	0.1091	0.0715	0.0144	1.6582	0.6031	1.1 µg/µl

Run 2ul on gel w/ 1kb ladder

Plasmid looks good  
 Strong plasmid #2.  
 LFC



abs	abs	bkg abs	260.0 nm
260.0 nm	280.0 nm	320.0 nm	280.0 nm
0.0942	0.0638	0.0201	1.6947

0.94 µg/µl

3/8/95

Start Culture to do an induction of  
 HIPANOS 81bp + PDE100

Inoculate 30ml LB Amp Kan with  
 C-7  
 inoculate 37°C ON

62 H1PM503 51p PD1D -

3/9/95

inoculate 300ml LB+Amp/Ken  
10ml of O/N culture of  
H1PM503 51p + PD1D 2C3D  
incubate 37°C w/ aeration  
until OD<sub>600</sub> ~ 0.4-0.6  
Add 100ml IPTG to 2mM (10ml)  
incubate 37°C 4 hours.  
Spin 8K 20min  
Remove supernatant  
Resuspend pellet in 30ml 1M Gm HCl  
pH 8  
Store 4°C O/N.

3/10/95

Spin Culture 8K 20 min  
transfer supernatant to fresh tube.  
Ken in gel.

400ul H<sub>2</sub>O  
25ul Protein in Gm HCl  
50ul 0.15% Na DOK  
75ul 50% TCA.

Mix well


Spin 10 min

Remove supernatant

Resuspend pellet in 2ul 0.2N NaOH

Run on 8.5% gel.

Stain Destain



looks good -  
- Load Buffer  
columns

PD 31

H14 Screen / HTA Screen

69

pg 60

3/10/94

Develop film for.

21A-40A HT4CM40  
1-20B HT4CB44

3/13/94

Pick positive clones of

HT4A133

HT4CM40

HT4CB44

into 20ul SM in 96 well dish

Pick 48 from each clone  
incubate samples at RT O/N

Plate HT4A  
HT4CB44 for screening of

Dilute 1/1000 into 20ul into  
6cm<sup>2</sup> 1E372 cells  $OD_{600} = 1.0$

incubate 37°C 15 min

Plate into 150mm NZY plates with

final LB + 0.35% Bacto. Let set

incubate 37°C 5 hours.

Store 4°C O/N

3/14/95

for the 48 clones in SM of

HT4A133

HT4CM40

HT4CB44

To 50ul LB372  $OD_{600} = 1.0$  - add 20ul

of phage

incubate 37°C 15 min

add 15ul NZYM broth

H-TPAD 51bp + PD10

71

pg 62

3/13/95

Prepare NiSO<sub>4</sub> Column.

2 ml Resin

Water 30 ml H<sub>2</sub>O

Strip 30 ml 0.1M HCl pH 2.

Water 30 ml H<sub>2</sub>O

Change 30 ml 0.1M NiSO<sub>4</sub>

Water 30 ml H<sub>2</sub>O

Equilibrate 30 ml 0.1M HCl pH 8.

Apply supernatant - Collect - flow


Water 30 ml pH 8 0.1M HCl - Collect

Water 30 ml 0.1M HCl pH 6 - Collect - pH 6

Strip 30 ml 0.1M HCl pH 5 - Collect

Strip 30 ml 0.1M HCl pH 2 - Collect 3/14/95

Run 20 ml on cpl. with uninduced cultures.



Protein looks good.

HTPAD 51bp + PD10

3/13/95

prepare 125% Preparative cpl  
1.5 mm

72

HTRANDS 80 5162 +PD10

3/15/95

Prep. 80 ul of pH 5 in H<sub>2</sub>O / 50% GMA  
 Add 0.15% NaBO<sub>2</sub>  
 Add 400 ul 0.3N NaOH  
 Add 400 ul 2X Dissociation Buffer  
 Heat 100°C 5 mins  
 Spin 2 min  
 Run on gel, 100V.

Cut off Marker & Port of gel.

Stain DESTAIN



Align w/ corresponding gel.

Cut out from unstained gel.  
 Place in 15ml conical  
 Ready for RT Production

133  
 PM

HT4/HTA screen

75

3/16/95

HT4CB44



HT4CM40

HT4CM40

HT4CB44



HT4CB44

HT4CB44

Wash HTA & HT4 filters  
0.2xSSC / 0.1% SDS. -3x 65°C  
Put on film  
-80°C O/N.

Develop film.

Steve P. Ph

3/17/95

32 Random Prime Probe - HSEEN37, HT4SB02, HSKB009  
4/6/95

Mix by Flicking  
Quick Spin  
Incubate 37°C 10 min

For HSEEN37 use Primer if

Primers	10ul
DNA	4ul
H <sub>2</sub> O	30ul
	44ul

Heat 100°C 5min  
Quick Chill  
Quick Spin

add 10ul 5x dCTP buffer  
5ul 83P dCTP  
1ul Klenow

Incubate 37°C 10min

put through G 25 Column

Count tubes

& Did Not put HT4SB02

HT4SB02 HSEEN37 HSKB009  
through Column

CPM	2SIG
-----	------

2391106.50	0.19 HSKB009
3005574.25	0.19 HT4SB02

1012483.00	0.20 HSEEN37
------------	--------------

Add 100ul Salmon Sperm DNA  
Heat 100°C 5min  
Quick Chill

N. A. Ph.  
4/6/95



Seem H0AAH/H0SEN (H0KBN/H1444/H14CB

103

4/27/95

Inoculate 30 ml TB+Amp  
with HTAB494, S01, S02, S03, S04, S05 & S06  
Incubate 37°C O/N

From plated resus - pick 6 white  
Colonies into ~~each~~ 200 ml TB+Amp  
Incubate 37°C 4 hours  
Do PCR

H0AAH02

FP50	1
M13R	0.05
10x	3.2
10x	3.2
H <sub>2</sub> O	22.4
Taq	0.15
Cult	2
	<u>32</u>

H0SEN37

FP50	1
M13R	0.05
10x	3.2
10x	3.2
H <sub>2</sub> O	22.4
Taq	0.15
Cult	2
	<u>32</u>

H0KBN09

FP50	2
M13R	0.05
10x	3.2
10x	3.2
H <sub>2</sub> O	21.4
Taq	0.15
Cult	2
	<u>32</u>

H144480

FP01	1.2
M13R	0.05
10x	3.2
10x	3.2
H <sub>2</sub> O	22.2
Taq	0.15
Cult	2
	<u>32</u>

H14CB44

FP02	1
M13R	0.05
10x	3.2
10x	3.2
H <sub>2</sub> O	22.4
Taq	0.15
Cult	2
	<u>32</u>

PCR Program



95°C	5 min
95°C	10 sec
55°C	20 sec
72°C	1 min
72°C	7 1/2 min
4°C	Hold

Sta. R. R.  
4/27/95

TNT - HTPA208 51bp Protein Prep

133

5/5/95



1972

5/8/95

Inoculate LB + Amp Kan with  
HTPA208 51bp ATG in PUD10  
into 100 ml.  
Incubate 37°C O/N.

Inoculate 250 ml LB + Amp with  
~~HTPA208~~ H.TABG 94500.  
for Maxi Prep

Inoculate 30 ml TB + Amp with  
HOAAT164 5000

Incubate 37°C O/N w/ aeration

5/9/95

Inoculate 1 l of LB + Amp Kan  
with 50 ml O/N culture  
of HTPA208 51bp ATG in PUD10  
Incubate 37°C until OD<sub>600</sub> = 0.4-0.6  
Add 100 mM IPTG to 2 mM - 20 ml

134

Maxi HTABG94S06. | midi HOAAH62S02

5/9/95

Incubate 37°C 4 1/2 hours

Spin 5K 20 min

Pour off Supernatant

Resuspend pellet in a total of 100ml  
0.1M Gm HCl pH 8

let sit on ice at 4°C

Design Max Prep  
of HTABG94S06  
1:200 Dilution

abs	abs	bkg abs	260.0 nm	280.0 nm
260.0 nm	280.0 nm	320.0 nm	280.0 nm	260.0 nm
0.1548	0.1071	0.0400	1.7092	0.5851

1.55 µg/ml

Run 0.5 µg on gel

HOAAH62S02. Alkylin Lyses #1

Spin Culture

Pour off Supernatant

Resuspend pellet in P1 + RNase  
(from Qiagen) - 2ml

let sit RT 5 min

Add 2ml P2

Mix gently

let sit on ice 10 min

Add 2ml P3

Mix well

let sit on ice 10 min

Spin 20 min 8K

Transfer Supernatant to fresh tube

Add isopropanol 0.7 Volumes - 4.2ml

Mix well

Spin 30 min 8K

HIPANOS 51bp ATG + PD10

135

Pour off Supernatant  
Wash pellet 70% Ethanol  
Spin 10 min 8K  
Allow pellet to Dry at RT O/N

5/9/95

5/10/95

H04A#02502

Resuspend pellet in 1ml TE  
Transfer to 2 microfuge tubes  
Add equal volume 13% PEG/1.6M KAc  
Mix well  
Spin 15 min

1X 70% Ethanol Wash  
Resuspend pellet in a total of 200ul of TE  
Run 1ul on gel -  
Read OD260/280 1:200 Dilution

OD260 OD280 OD260/280

0.0279 0.0168 0.0034 1.8322 0.5458

0.28ug/ul.

HIPANOS 51bp ATG + PD10.

Spin 8K 20 min  
Transfer Supernatant to fresh tube.  
Prepare NiSO<sub>4</sub> Column.

Pour 3ml Resin into Column.  
Wash 20ml H<sub>2</sub>O  
Equilibrate 20ml 10mM Gln HCl pH 8  
Pour on Supernatant (crude extract)  
and Collect Flow through.

Wash Column 60ml 10mM Gln HCl pH 8  
Collect - pH 8  
Wash Column 60ml 10mM Gln HCl pH 6  
Collect - pH 6  
Elute ~~with~~ Protein - 5ml 10mM Gln HCl pH 5  
Collect - pH 5

136

HTPAND3504 61bpA1G1 + PD10

5/10/95

Strip Column 30ml 10M Gm HCl pH 2  
collect pH 2

5/11/95

Run Protein Samples on 12.5% Gel  
To Samples in 10M Gm HCl450  $\mu$ l H<sub>2</sub>O50  $\mu$ l of Sample50  $\mu$ l of 0.15% Na DOK35  $\mu$ l 50% TCA

mix well


Spin 10 min

Remove Supernatant

Resuspend in 10  $\mu$ l of 0.2 N NaOHAdd 10  $\mu$ l of Dissociation Buffer

Heat 5 min 100°C

Spin

Run Samples on 5% SDS PAGE Stacking  
gel. 150 V 1 1/2 Hours

- 1 ~~Crude~~ - Uninduced
- 3 ~~Crude~~ - Crude Extract
- 4 ~~Crude~~ - pH 8
- 5 ~~Crude~~ - pH 6
- 6 ~~Crude~~ - pH 4
- 7 ~~Crude~~ - pH 5
- 2 - Rainbow Marker

Does not look good try increasing loading  
volume and do another  
induction

~~800~~ 800 ~~run~~

# Fragment preps. + Northern Blots

187

5/17/95

Set up Digests of unknowns to give  
to Brent Krieger - He'll do  
Northern.

Clone ID	Conc	5µg DNA	H <sub>2</sub> O	10X P2	XhoI	EcoRI
HSKBND09	2.0 µg/µl	2.5	42.1	5 µl	0.2	0.2
HNBAAZ6	0.6 µg/µl	7.5	37.1		0.2	0.2
HILBY30	0.5 µg/µl	10	34.6			
HTLAT33	PCR Product	20	24.6			
HT4CM40	0.4 µg/µl	11	33.6			
HTLCB44	1.2 µg/µl	4.2	40.4			
HNFAAG4	0.73 µg/µl	6.8	37.8			
HT4AY80	0.34 µg/µl	15.6	29			
HTABG894	0.93 µg/µl	5.2	38.4			

Digest 37°C O/N.

5/12/95

Run 5µl on gel with 1 kb marker  
XhoI EcoRI Digests



- 1 - HSKBND09
- 2 - HNBAAZ6
- 3 - HILBY30
- 4 - HTLAT33
- 5 - HTLCM40
- 6 - HTLCB44
- 7 - HNFAAG4
- 8 - HT4AY80
- 9 - HTABG894

5/12/95

## Fragment preps - Northern Blots

5/12/95

Run on 0.8% LMP Gel with 1 kb ladder.  
 Run 80V 2 1/2 hours



Cut out fragments  
 - place into 1.5 ml  
 microfuge tube  
 Store 4°C over  
 weekend.

Need to digest  
 HILB450 Again

201


H-TRANS 5.1 kb ATG in P.D.V

- INDUCTION
- Resuspend in 40  $\mu$ l COM Gen HILPH8
- Store 4°C over weekend.

5/15/95

Gene Clean fragments  
 Resuspend in 30  $\mu$ l TE

Run 1 kb on gel with 1 kb ladder



- 1 HSKBAX9
- 2 HNBAA216
- 3 HT4AT33
- 4 HTUCM00
- 5 HTUCB44
- 6 HNSFAAGW
- 7 HT4AV80
- 8 HT4BG94
- 9

# Fragments prep - Northern Blots

139

Send to Brent Knäder - Rand Kovacs

5/15/95

Clone ID	Libraries Expressed	Size in Kb	Eco RI/Xho I	Approx. Conc(ng/ul)
HSKBN09	HT4, HSK	-1.4	AMK	-100 15ul
HNBA26	HBJ, HBM, HCA, HNB, HRG, HTO	-0.800		-100 ↓
HILBY30	HIL, HLM	-0.60		Will Send 3ul
HT4AI33	HT4	-0.90		-150 15ul
HT4CM40	HT4	-1.7		-150
HT4CB44	HT4	-2.5		-150
HNFAA64	HNH, HSI, HSK, HTA	-1.80		-150
HT4AY80	HT4, HTX, HT3	-0.85		-50
HTABG94	HTA	-1.7		-150
HT4CI56	HT4	-1.7		-100
HT4AI55	HT4	-1.7		-140
HT4CL32	HT4, HT5, HT3	-1.1		-150
HT4CA46	HT4, HT3, HCE, HGO, HTA, HL3, HMW	-1.7		-250
HMSAF22	HMS, HOS, HHF, HSR HTN	-2.0		-100 15ul
HT3SB70	HNH, HT3, HT4, HT5, HTA	-1.5	CLF	-300 5ul
HT4SB02	HET, HGL, HHF, HSU, HT4, HTE, HTP	-1.3	AMK	-200 10ul

Set up Digest of HILBY30 PCR product

DNA	20ul
YOH	5
H2O	24.6
EcoRI	0.2
Xho I	0.2
	30ul

Digest at 37°C overnight

HIPANC8 5' bpATG + PD10  
 Put over column  
 collect Flow-through  
 pH 8  
 pH 6  
 pH 5




140

HTRANDS 516p AT4 + PD13

(ms136)

Run HTRANDS Samples on 12% Acrylamide gel



- 1 Rainbow Marker
- 2 Crude Extract
- 3 flow through
- 4 pH 8
- 5 pH 6
- 6 pH 5
- 7 pH 2.

Store at 4°C

Carrie Did Gene clean of HTLBY30  
Xho P.I.

Run gel on gel with 1Kb ladder

Give to Brent  
~600 bp  
~100 ng gel.

HTRAND8504 5bp ATG APD/O

141

5/18/95

Reapply pH5 of HTRAND8504 ATG to  
fresh column. (from 3/13/95 - pg 31)

Add 3ml pH5 + 2ml pH8.

Wash 3ml pH8.

Wash 3ml pH8.


Give to protein expression for  
reanalysis over column

5/19-5/22

Computer

5/23/95

Received Column Back  
Elute 2nd Immobilized elution buffer.  
2 times - Run in 12% Gels  
with Marker + HTRAND8504 ATG



Computer

off.

5/23, 5/24  
5/25

5/26/95

pg 28  
Book #245  
p. 10

HTRAND8504 Sub ATG HAD10

141

5/18/95

Relaply pH5 of HTRAND8516pATG to  
fresh Column. L from 3/21/95 - pg 311

Also 3ml pH5 + 3ml pH8.

Wash 3ml pH8.

Wash 3ml pH5.


Give to Protein expression for  
reanalysis over column.

5/19/95

Computer

5/23/95

Received Column Back  
Elute 2ml Immobilized elution Buffer.  
2 times - Run in 12% Gels  
with Marker + HTRAND8516pATG



Computer

off

5/23, 5/24  
5/25

5/26/95

PJ28  
HTRAND8516pATG  
Marker + 10

**This Page is Inserted by IFW Indexing and Scanning  
Operations and is not part of the Official Record**

**BEST AVAILABLE IMAGES**

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

- ☐ **BLACK BORDERS**
- ☐ **IMAGE CUT OFF AT TOP, BOTTOM OR SIDES**
- ☐ **FADED TEXT OR DRAWING**
- ☐ **BLURRED OR ILLEGIBLE TEXT OR DRAWING**
- ☐ **SKEWED/SLANTED IMAGES**
- ☐ **COLOR OR BLACK AND WHITE PHOTOGRAPHS**
- ☐ **GRAY SCALE DOCUMENTS**
- ☐ **LINES OR MARKS ON ORIGINAL DOCUMENT**
- ☒ **REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY**
- ☐ **OTHER:** \_\_\_\_\_

**IMAGES ARE BEST AVAILABLE COPY.**

**As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.**