


Steven M. Ruben  
Appl. No. 10/662,429


BEST AVAILABLE COPY

Human Genome Project  
Lab Notebook # 532

Department Mol. Biol.  
Subject 10/20/94 - 02/02/95  
Name ANN KIM # 8  
Address \_\_\_\_\_

 43-648

**Computation Notebook**  
Dennison Stationery Products Co., Framingham, MA 01701

 75 Sheets  
11 3/4" x 9 1/2"  
4x4 Quad.

0 73333 43648 8

Ruben EXHIBIT #91

HTPB411 & HTPAN08 in PD10

3

pg150 Book 7

10/20/94

# Ligations

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19
HTPAN08 9/30	4	4	4																
10/20				6	6	6													
HTPB411 9/30							3	3	3										
10/20										6	6	6							
PD10 1.1													6	6	6				
PD1024 10/11	2			2			2			2			2			2			
10/14		2			2			2			2			2			2		
10/14			2			2			2			2			2			2	
10x Buffer	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2
H <sub>2</sub> O	11	11	11	9	9	9	12	12	12	9	9	9	9	9	9	15	15	15	17
T4 Ligase	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
	20	20	20	20	20	20	20	20	20	20	20	20	20	20	20	20	20	20	20

Incubate @ 16°C O/N

10/24/94

Q1

Transform into M15 Chemically Competent Cells

Thaw M15 Cells on ice  
 Aliquot 100ul into @ fresh Sterile tubes  
 Add 10ul of ligation to tubes.  
 Incubate on ice 1 hr  
 Heat 42°C 45 sec  
 Place on ice  
 Add 400ul LB  
 Incubate 37°C 1 hr  
 Plate 175ul onto LBt Amp<sup>r</sup> Kan plates  
 Incubate RT over the weekend

HEAVY DUTY NOTEBOOK  
COMPUTATION NOTEBOOK

Department Mar. Biol.  
Subject 10/20/94 - 02/02/95  
Name ANN KIM # 8  
Address \_\_\_\_\_



43-648

### Computation Notebook

Dennison Stationery Products Co., Framingham, MA 01701



75 Sheets  
11 1/4" x 9 1/4"  
4x4 Quad.

0 73333 43648 8

Ruben EXHIBIT 2091  
Ruben v. Wiley et al.  
Interference No. 105,077  
RX 2091

4

HTPB411 &amp; HTPAN08 in PD10

4/94

Plates look OK.

PD10 1.1 + 2.4 looks very good so legations worked well

PD10 vector colonies had many colonies  
 so did Vector + ~~PD10~~ fragment  
 try PCR on colonies.

Pick 200 each into 200  $\mu$ l LB + Amp + Kan  
 in 96 well dish - HTPAN08 & HTPB411

Incubate 37°C 4 hrs w/aeration

PCR

HTPAN08

		200 $\mu$
10x PCR	3.2	640
10x dNTP	3.2	640
2x99	0.2	40
2500	0.2	40
Taq	0.2	40
H <sub>2</sub> O	23	4600
Cult	2	

HTPB411

		200 $\mu$
10x PCR	3.2	640
10x dNTP	3.2	640
2501	0.2	40
2502	0.2	40
Taq	0.2	40
H <sub>2</sub> O	23	4600
Culture	2	

PCR Program #69

95°C 5 min  
 95°C 20 sec  
 55°C 20 sec } 30x  
 72°C 1 min  
 72°C 7 1/2 min  
 4°C Hold.

HTPB411 & HTPA008 in PP10

5

10/24/94



Run 10  $\mu$ l on gel w/ 1 kb ladder

- Nothing worked!

Didn't bother to take all pictures.

HTPA008

had strange size band.



Retry - Make New Primers today later.

- Phenol extract from PCR - then ethanol ppt - Digest O/X.
- Run on LMP gel.

10/26/94

Try using S' Bam primers from  
ligation into pAZ vector.

HTPA008 - 7689

HTPB411 - 2429.

6

HTPAN03 + HTPB11 8' in PD10

10/26/94

DNA

100µl


	HTPAN03			HTPB11	
	①	②		③	④
10xPCR	10.0	10.0	10xPCR	10	10
10x dNTP	10.0	10.0	10x dNTP	10	10
2499	0.5	0.5	2501	0.5	0.5
2500	0.5	—	2502	0.5	—
768A	—	0.5	2429	—	0.5
Taq	0.3	0.3	Taq	0.3	0.3
H <sub>2</sub> O	77.7	77.7	H <sub>2</sub> O	77.7	77.7
DNA	1	1	DNA	1	1
	100	100		100µl	100µl

Run PCR Program #04. Modified

95°C 5min  
 95°C 20sec  
 55°C 20sec  
 72°C 1min  
 72°C 7 1/2 min.  
 4°C hold

25X

Run 5µl on gel with 1Kb ladder



looks good.

Add equal Volume PEG/NaCl  
 Sit on ice 10min  
 Spin 10min  
 Remove supernatant  
 1x70% ethanol wash - 500µl

Spin 5min  
 Remove supernatant  
 let pellet dry slightly  
 Resuspend in 200µl PEG-TG  
 Add equal Volume Phos  
 Vortex

HTPAD008 + H7PB008 in R10

7

10/20/91

Spin 5 min

Transfer upper aqueous layer to fresh tube  
repeat

Add 2 volumes Ethanol (400  $\mu$ l)  
1/10 volume 3M NaAc (20  $\mu$ l)

Mix

Let sit on ice 1/2 hr.

Spin 10 min

Remove supernatant

1x 70% Ethanol wash - 1000  $\mu$ l

Spin 5 min

Remove supernatant

Allow pellet to dry

Resuspend in 200  $\mu$ l TE.


Set up digestions.

DNA	100 $\mu$ l.
10x #2	20 $\mu$ l.
H <sub>2</sub> O	79 $\mu$ l.
XbaI/Bam	0.5/0.5
	<hr/> 200 $\mu$ l.

Incubate 37°C / overnight

10/27/91

Run 10  $\mu$ l on 1% gel with 1 kb ladder.



gel looks good

Add sample Dye. Angel

Run on 0.8% 1MP

Agarose gel


Run 80V 1 1/2 hrs.

Ethanol ppt samples (see above) - 1x 70% wash  
Resuspend in 40  $\mu$ l TE - Now ready for gel

8

HTPA008 &amp; HTPB411 in PD10

10/27/94




for Sample #1 cut out  
1.8 Kb fragment too

Gene Clean  
Add 1000  $\mu$ l NaI  
Heat 55°C 5 min  
Mix well  
Add 1  $\mu$ l glass milk  
Let incubate at  
Room temp 5 min  
w/ occasional mixing  
Spin 7 Sec

3X { Remove Supernatant  
Resuspend pellet in 400  $\mu$ l Wash Buffer  
Spin 7 Sec  
Remove Supernatant  
Spin 7 Sec  
Remove Supernatant  
Resuspend pellet in 20  $\mu$ l TE  
Heat 55°C 1 min  
Spin 7 Sec  
Transfer to fresh tube  
Repeat 20  $\mu$ l TE

Rem 2  $\mu$ l on gel with 1 Kb ladder



Use to Set up digests

Do Not use the  
1.3 Kb fragment from  
#1



HTPANOS + HTPBY11 in pD10

9

10/27/94

# Ligation Reactions.

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	Ke
HTPANOS #1	4	4							4							
HTPANOS #2			4	4						4						
HTPB411 #3					4	4					4					
HTPB411 #4							4	4				4				
PD10 1)																4
PD10 10/19	1	-	1	-	1	-	1	-					1			1
PD10 10/12	-	1	-	1	-	1	-	1						1		
10x Buffer	2	2	2	2	2	2	2	2	2	2	2	-	-	-	-	-
H2O	12	12	12	12	12	12	12	12	13	13	13	13	16	16	17	12
T4 Ligase	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
	20	20	20	20	20	20	20	20	20	20	20	20	20	20	20	20

Incubate at 16°C O/N.

10/28/94

## Transform M15 Chemically Competent cells.

Thaw Cells on Ice  
 Aliquot 100ul of cells into fresh tubes  
 #16 may get 450ul.  
 Add 10ul of ligation to tubes  
 Incubate on ice 1 hr.  
 Heat 42°C 25 sec  
 Put on ice  
 Add 400ul LB  
 Incubate 37°C 1 hr  
 Plate onto LB + Amp + Kan plates, 150mm.  
 plate 1-16 100ul  
 1-8 250ul also  
 plate #6 (100ul) mixed w/ #7

10

HTPANOS

poll in

iii PDID

10/28/94

Incubate  $37^{\circ}\text{C}$  O/N.

Asked Paul Moore

O/K.

place at  $4^{\circ}\text{C}$

Tomorrow.

10/31/94

Panel did NOT take plates out  
- there were colonies. So try to PCR

Dick 24 of each 1-8 & PCR

14  
TRANOS

545-8  
7TPB4

	1-2	50	3-4	50	5-6	7-8
10x PCK	3.2	160	3.2	160	3.2	160
10x dNTP	3.2	160	3.2	160	3.2	160
2499	0.2	10	2499 0.2	10	2501 0.2	10
2500	0.2	10	7689 0.2	10	2502 0.2	10
Taq	0.2	10	0.2	10	0.2	10
H <sub>2</sub> O	23	1150	23	1150	23	1150
<del>DNA</del>	B					
Cult	<u>2</u>	—	<u>2</u>		<u>2</u>	<u>2</u>

✓ CC Proxy 69.  
✓ PH Proxy let up.

④ make.  
Try by own.

Bob

13

# HTPB111 & HTPAN08 m. PD10

13

(Pg 12)

10/31/94

Make new insert again  
PCR w/ New and old primers.

	(1)	(2)	(3)	(4)
HTPAN08				
3' Xba 2499	0.2	—	0.2	—
3' Xba 2656	—	0.2	—	0.2
5' Bam 2500	0.2	—	—	0.2
5' Bam 7689	—	0.2	0.2	—
10x PCR	10	10	10	10
10x dNTP	10	10	10	10
Taq	0.2	0.2	0.2	0.2
DNA (100ng/ul)	1	1	1	1
H <sub>2</sub> O	78.4	78.4	78.4	78.4
	100	100	100	100

	(5)	(6)	(7)	(8)
HTPB111				
3' Xba 2501	0.2	—	0.2	—
3' Xba 2653	—	0.2	—	0.2
5' Bam 2502	0.2	—	—	0.2
5' Bam 2429	—	0.2	0.2	—
10x PCR	10	10	10	10
10x dNTP	10	10	10	10
Taq	0.2	0.2	0.2	0.2
DNA (100ng/ul)	1	1	1	1
H <sub>2</sub> O	78.4	78.4	78.4	78.4
	100	100	100	100

PCR Prog # 619

95°C 5min  
 95°C 20 sec  
 58°C 20 sec  
 72°C 1min  
 72°C 7 1/2 min  
 4°C Hold

25x


14

HTPAND08 &amp; HTPB411

in PD10

1/1/94

Run 5ul on gel with 1Kb ladder.



Add equal Volume PEG/NaO  
Spin 15 min  
70% Ethanol Wash.  
Resuspend 100ul.  
Digest O/N.


DNA	50ul
H <sub>2</sub> O	38
10x#2	10
Xba/Bam	1/1
	100ul

incubate reaction  
37°C Overnight

3

11/2/94

Prep. Precipitate Digest  
Add 1/10 Vol 3M NaAcetate  
2 Vol 100% Ethanol  
Spin 15 min  
70% Ethanol Wash.  
Allow pellet to dry  
Resuspend in 45ul of TE  
Add 5ul 10x Sample Buffer  
Run on 0.8% LMP gel.  
HTPAND08 HTPB411



HTP4008 & HTP4111 8 in PMS

15

11/2/94

Con  
Cut out fragment

Gene clean

Add 1000ul Pa I

Heat 65°C 2 min

Mix well

Add 7ul Glass milk

Incubate at Room Temp 5 min

With occasional mixing

Spin 7 sec

Remove Supernatant

Resuspend pellet in 400ul Wash

Buffer

Spin 7 sec

Remove Supernatant

Spin 7 sec

Remove Supernatant

Resuspend pellet in 20ul TE


Heat 65°C 1 min

Spin 10 sec

Transfer to fresh tube

Repeat with 20ul TE

Run 2ul of eluted fragment on  
1% TAE gel with 1 kb ladder



With fragments  
from 10/27  
and the New Fragment  
Set up ligations

16

HTRANSACTORY 11/2/94 into PD10

11/2/94

## ligation Reactions

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
HTRANS 2499+2500 1/2	3														
2499+7689		3													
11/2 2499+2500			3												
2656+7689				3											
2499+7689					3										
2656+2500						3									
HTRANS 2501+2502 1/2							3								
2501+2429								3							
11/2 2502+2502									3						
2653+2429										3					
2501+2429											3				
2653+2502												3			
PD10 Bam/Hin	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
10x Lig Buffer	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2
H <sub>2</sub> O	13	13	13	13	13	13	13	13	13	13	13	13	13	13	13
T4 Ligase	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
PD10 .1 kb															3

Incubate 16°C O/N

11/3/94

Transform M15 cells

Use 1-8 from 10/27 ligation  
as well as 1-15 from 11/2

To 100 µl of cells add 10 µl of ligation

Let sit on ice 1 hr  
Heat 42°C for 45 sec  
Add 100 µl LB

Incubate 1 hr at 37°C

Plate onto LB + Amp + Kan  
plates

HTPAW08 + HTPB411 in pD10

17

11/5/94

150 mm plates

use 150 ul 1-8 & 1-12  
use 300 ul 1-8 & 1-12

Incubate 37°C O/N.

11/4/94

Pick 1-8 10/27 16/28/29  
into 200 ul LB + Kan + Amp plates  
Pick 1-12 11/2 into  
200 ul LB + Amp + Kan

Incubate 37°C w/ aeration  
120 PCR

use internal primers for confirm

HTPB411 RPI0 + FPI5  
HTPAW08 RPI2 + FPI8

HTPB411

		200x
2400	0.25	50
2.7945	1	200
10x PCR	3.2	1400
10x dNTP	3.2	1400
Taq	0.2	60
H <sub>2</sub> O	2.15	2430
Culture	2	

HTPAW08

		200x
16638	0.25	50
3407	0.25	80
10x PCR	3.2	240
10x dNTP	3.2	240
Taq	23.1	460
H <sub>2</sub> O	0.2	40
Culture	2	

PCR Program HGA

95°C	5 min
95°C	20 sec
55°C	20 sec
72°C	1 min
72°C	7 1/2 min
40°C	hold

Stir Exsist  
-20°C pres  
incubated

Control Plasmid  
Control H<sub>2</sub>O

18

HTPAN08 + HTPB411 in PD10

11/7/94

Rem 10  $\mu$ l of Reaction on gel work (163)  
 loaded

Nothing set up again!

Talked w/ Steve -

Need to remake primers,  
 Sites on old ones wrong

Submit both HTPAN08 + HTPB411  
 5' Bam a 3' ~~At~~ Xba

11/8/94

3' Xba HTPB411 - Oligo Rec'd.  
 (60°C. O/X)

11/10/94

Rec'd 5' Bam New HTPB411  
 5' Bam New HTPAN08

11/11/94

Rec'd 3' Xba New HTPAN08

11/12/94

Dried Oligos

Rem PCR Reactions w/ new Primers

39



HTPAN08 / HTPB411 in PD10

39

pg 18

11/14/94

HTPAN08S04

HTPB411S15


DNA (10ng/ul)	1
3' Xba New	0.6
5' Bam New	0.7
10x dNTP	10
10x PCR	10
H <sub>2</sub> O	72.1
Taq	0.3
	100

DNA (10ng/ul)	1
3' Xba New	0.6
5' Bam New	0.7
10x dNTP	10
10x PCR	10
H <sub>2</sub> O	71.4
Taq	0.3
	100

PCR Program #69 modified

95°C 5min  
 95°C 20sec  
 55°C 20sec } 25X  
 72°C 1min  
 72°C 7 1/2 min  
 4°C Hold

11/14 Run 5ul on gel w/ kb ladder



Add equal Volume PEG-NaCl  
 Spin 10min  
 Pour off Supernatant  
 Wash 1ml 70% ethanol  
 Spin 5min  
 Remove Supernatant  
 Resuspend in 100ul TE  
 Digestions:

DNA	2.5
H <sub>2</sub> O	1.9
10X #2	5
Xba / Bam	0.5 / 0.5
	30ul

↑ 200bp Digest

Incubate RT over full Monday 11/14


40 HTPA008 / HTPB411 in PDI0

11/14/94

Remain

Run digest on gel - see pg 39.  
with 1 kb ladder

Run Digest on 0.8% LMP Gel  
Run 80V 2hrs



Cut out fragments  
take picture

Use Qiaquick Gel Extraction  
Kit

- Add 900ul QXI Buffer
- incubate 60°C 10 min
- Place in Spin Column
- Spin 60 sec - Remove Flowthrough
- Add remaining gel/QXI
- Spin 60 sec
- Remove Flowthrough

Add 750ul Buffer PE

Spin 60 sec

Remove flowthrough

Add 750ul Buffer PE

Spin 60 sec

Remove flowthrough

Spin 60 sec


Transfer Column to fresh tube

Add 50ul TE Heated 60°C

Spin 60 sec

Run Gel on gel with 1 kb ladder

looks good.



HTPANC8 / HTPB411 in PD10

41

11/14/94

Set-up ligations

	①	②	③	④	⑤
HTPANC8	4	—	—	—	—
HTPB411	—	4	—	—	—
PD10 1:1	—	—	4	—	—
PD10 2:4	1	1	1	1	2
10X Buffer	2	2	2	2	2
14 ligase	1	1	1	1	1
H <sub>2</sub> O	12	12	12	16	17

Set up 2 sets - 1 for 16°C overnight  
A-1 set for RT 1 hr.

for A - 1 hr at RT

100 µl ligation

100 µl H15 Chem Competent cells

let sit on ice 1 hr

heat 42°C 45 sec

Place on ice

Add 400 µl LB

Heat 37°C 1 hr

Plate onto LBt Amp<sup>r</sup> Kan plates

100 mm - 50 µl

150 mm - 200 µl

Incubate 37°C O/N

11/15/94

Plates look OK

very few colonies on Vector alone

No colonies on lig rxn alone

11/15/94

Pick Colonies from 1, 2 & 3 into  
LB + Amp + Kan.  
Incubate 37°C 4 hrs.

Run PCR Reagents

HTPAN08504	
3' Xba New	0.2
5' Bam New	2
10x dNTP	3.2
10x PCR	3.2
H <sub>2</sub> O	20
Temp	0.2
Culture	2
	32

HTPB411S05	
3' Xba New	0.2
5' Bam New	2.5
10x dNTP	3.2
10x PCR	3.2
H <sub>2</sub> O	20.7
Temp	0.2
Culture	2
	32

① Control

PCR &amp; insert -

② Control LB + Amp + Kan

with vector alone Colonies

take both sets of primers &amp; make sure

No contamination by fragments

Run PCR Program 69.

95°C	5 min	30x
95°C	20 sec	
55°C	20 sec	
72°C	1 min	
72°C	7 1/2 min	
4°C	Hold.	

Run Dial on gel with 1 kb ladder  
in 1% Agarose gel in TAE


(pg 47)

HTPAN08/HTPB411 m PD10.

47

(pg 42)

11/15/94



HTPAN08S04

HTPAN08S04  
looks good.  
HTPB411S05

Nothing so Re  
plate ligations  
from 16°C O/N.



HTPAN08S04 -

Incubate 5ml TB+  
Amp + Kan with 15ul  
of culture 2-11  
Incubate 37°C O/N.

With ligation Rxns left at 16°C O/N

Thaw m15 chemically competent cells once

To 100ul of thawed cells add

10ul of ligation Rxn

45 (tube) add 10ul PD10 ~~reaction~~ planned

DNA as (+) control.

Incubate on ice 1 hr.

Heat 42°C 45 sec

Place on ice

Add 400ul LB

Incubate 37°C 1 hr.

Plate onto LB + Amp + Kan plates

- 100 mm plates - 50ul

- 150 mm plates - 250ul.

Incubate 37°C O/N.

11/10/94

Boiling Min pups

Spin Outlines 2min

Remove Supernatant

Resuspend pellet in 750ul STEH

Rnase + Lipzyme

Boil 45 sec

Spin 10 min

Remove Pellet

Add equal Volume (750ul) of  
13% PEG / 1.6M NaCl

Mix Well

Spin 10 min

Remove Supernatant


Add 1000ul 70% Etanol to wash  
pellet

Spin 5 min

Remove ~~pellet~~ Supernatant

Allow pellet to dry slightly at RT

Resuspend in 200ul TE

Run 2ul on gel with 1kb ladder  
and pD10

plasmid looks  
good.  
looks like a lot of  
chromosome.

Digest with Bam 1/kb to see if  
the insert will "Pop" out  
N800 bp

pg 53

HTPAW08 + HTPB411 is RPO10.


58

(pg 48)

11/18/94

DNA	5ul
H <sub>2</sub> O	17.1
10x#2	2.5
Kam	0.2
Xba	0.2
	<hr/> 25ul

Incubate 37°C 4hrs  
Run on gel with 1 kb ladder  
Digest pD10 - Should appear  
1.1 kb fragment



looks like all are  
correct  
Submit for sequencing  
with pD10 5'  
primer

FABPDD1-10 RPO1

Setup digestion for HTPB411

Fragment	10
Vector	1
T4 ligase	1
10x Buffer	2
	<hr/> 20ul

Incubate RT 5 hrs

Thaw M15 cells on ice

To 100ul of M15 Chemically Competent cells

add 10ul of ligation reaction - ~~add~~ Use up all of ligation

Incubate on ice 1 hr

Heat 42°C 45 sec

Place on ice

Add 400ul LB

Incubate 37°C 1 hr

Plate 250ul into LB + Amp + Kan

plates - 150 mm

HTPAN08 + ~~HTPCE9~~ / 16 + HTPB4/11 in PD10


11/17/94

Some Colonies grew.  
 Pick into 200  $\mu$ l LB + Amp + Kan  
 in 96 well dish.  
 Incubate 37°C 4 hrs  
 PCR -

3' Xba	0.2
5' Bam	2.5
10x PCR	3.2
10x dNTP	3.2
Taq	0.2
H <sub>2</sub> O	20.7
Culture	2
	32.

PCR Program #69  
 use HTPB4/15/15  
 3' Xba New / 5' Bam New  
 Fragment as  
 positive Control

Run 10  $\mu$ l on gel with 1 kb ladder.



looks like  
 1 positive !!!

Incubate  
 50 ml TP + Amp  
 Kan  
 with #20

Incubate 37°C  
 O/N

11/18/94

Dragon Modified Media  
 Spin Culture 4.5K 15 min  
 Pour off supernatant  
 Resuspended in Tris PI



HTPA008 & HTPB411 in PD10

55

11/18/98

Add 7ml P2

Mix gently

Add 7ml P3

Mix

Set on ice 20 min

Spin 4.5 K 15 min

Equilibrate TIP-100 with 5ml

QBT

Apply Supernatant to tube through

Rem wipe

allow to flow through

2X Wash Column 10ml @ QC Buffer

Elute DNA 10ml QF

Add 7ml isopropanol (0.7 volumes)

Mix well

Spin 8K 30 min

Remove Supernatant

Wash 10ml 70% Ethanol

Spin 8K 15 min

Pour off Supernatant carefully


Let pellet dry at RT

Resuspend in 30ul TE

Run gel on gel with 1 Kb ladder

Dilute 1:200 in H<sub>2</sub>O

Read D<sub>260</sub>/D<sub>280</sub>



Sample ID	abs 260.0 nm	abs 280.0 nm	abs 320.0 nm	260.0 nm / 280.0 nm	260.0 nm / 320.0 nm
1	0.0557	0.0350	0.0031	1.6777	0.5708

assayed

HTPB411  
PD10

~ 160 ng total

Digest with Bam / Xba to pop out 2.4 Kb fragment

DNA 4

H<sub>2</sub>O 13.16

10X 2

Bam 0.2

Xba 0.2

Incubate 37°C

2 hrs


56

HTPAN08 + HTPB411

in PD10

11/18/94

Run 5ul on gel w/ 1Kb ladder  
 a 1ul of undigested.



11/28/94

Made New PA2 Primers 5' end  
 for HTPAN08.

2782 5' Bam 51bp New 7.1pmul/ul  
 2793 6' Bam 185bp New 8pmul/ul  
 Use 3 A<sub>1</sub> from previous PA2 cloning.  
 7690

PCR.

5' Bam	3
3' A <sub>1</sub>	0.5
POX PCR	10
10x dNTP	10
H <sub>2</sub> O	75.2
Taq	0.3
DNA (10ngul)	1
	<hr/> 100

[POX]

PCR Program # 69.


95°C	5 min	} 25 X
95°C	20 sec	
55°C	20 sec	
72°C	1 min	
72°C	7 1/2 min	
11°C	HOLD	

HTPANOS + HTPBY11

57

11/21/94

Rem 5ul on gel with 1Kb ladder



Add equal volume  
PEG / SacI - mix

Spin 10 min

Remove Supernatant

Wash pellet 1ml 70°C

Ethanol

Resuspend in 300ul TE


Rem 2ul on gel with  
1Kb ladder

Digest fragment

Bam / ~~Asp~~

DNA	20
H <sub>2</sub> O	23
10X#2	5
Bam	1
<del>Asp</del> xba	1
50ul	

Digest 0/1N  
store 20°C



11/28/94

Digested w/ xba by mistake.  
Rel Digest with

Bam / ~~Asp~~ Asp.

DNA	20
H <sub>2</sub> O	23
10X#2	5
Bam	1
Asp	1
50ul	

incubate 0/1N 37°C.

## HTRANS + PAZ lig

11/30/94

Isolate on 0.8% LMP Gel &  
 open clean  
 (see pg 102)

Set up ligations

185 bp Bam/Asp	10	—	—
51 bp Bam/Asp	—	10	—
PAZ Bam/Asp	1	1	1
10x Buffer	2	2	2
T4 Ligase	1	1	1
H <sub>2</sub> O	16	16	16
	20	20	20

Incubate 16°C O/N.

12/1/94

Transform into M15 cells

Thaw M15 Chem. Competent Cells  
 To 10  $\mu$ l of ligation add 100  $\mu$ l  
 of cells  
 Incubate on ice 1 hr  
 heat 42°C 45 sec  
 Place on Ice  
 add 400  $\mu$ l LB  
 incubate 37°C 1 hr.  
 plate 100  $\mu$ l on 100 mm plates  
 300  $\mu$ l on 150 mm plates

12/2/94

M15 cells contaminated?

(pg 107)

62

PD10. B/HIE / HTPANOS frag / DNase P frag

11/30/94


Del isolate fragment on  
0.8% agarose gel

Also Red. Lam.

- 1 HTPANOS 185bp Bam/Aap
- 2 HTPANOS 51/185bp Bam/Aap
- 3 DNase P Bam/Xba
- 4 PD10 Bam/HII

Wash 1 kb ladder

Run 80V 1 1/2 hrs  
Cut out bands & take picture



Gene Clean fragments  
Add 1ml Na.II

Heat 55°C 5 min

Add 5ul Glass milk  
incubate RT 2 min

Spin 10 sec

Remove Supernatant

Resuspend pellet in

50ul Wash Buffer

Spin 10 sec

Remove Supernatant

Spin 10 sec

Remove all supernatant possible

Resuspend pellet in 20ul TE

HTPANOS 51bp + 185bp DNase P

Resuspend pellet in 40ul TE PD10 Bam/HII


incubate 55°C 2 min

Spin 10 sec

Transfer Supernatant to fresh tube

Re-plate

Run 2ul on 1% gel with 1 kb ladder



store at -20°C

HTPANOS - see pg 58  
DNase P - same


HPCC91, HE20142

65

pg60

12/5/94

- ① HE20142S15
- ② S14
- ③ 1 kb ladder
- ④ HPCC91S14
- ⑤ S15



Decided to use a 3' fragment for  
screening instead  
Inoculate 150 ml TB + Amp.  
HE20142S14  
HPCC91S14

12/6/94

Diagnose max prep of HE20142S14  
HPCC91S14

Spin cultures 4.5 K 15 min  
Pour off supernatant  
Resuspend pellet in 10 ml P1 + RNase  
let sit at RT 5 min  
Add 10 ml P2 mix gently  
let sit at RT 5 min  
Add 10 ml Cold P3 & incubate on  
ice 20 min  
Spin 4.5 K 15 min  
Transfer to Cofence 10 ml QBT  
Apply supernatant to Cofence  
Wash 2x 30 ml QC  
Elute 15 ml QF  
PPT w/ 10.5 ml isopropanol  
Spin 9K 30 min  
Pour off supernatant  
Wash 10 ml 70% ethanol  
Spin 9K 10 min  
Pour off supernatant

Next Ph  
12/15/94

HTPAW08504 51bp/185bp in PA2

67

pg 58

12/2/94

The lig Rn alone & M15 cells alone  
produced colonies -  
Don't use M15 - use DH5α.

Thaw Chemically Competent DH5α  
cells on ice

12/3/94

10 100ul Cells add 10ul of  
ligation

incubate on ice 1 hr

Heat 42°C 45 sec

Put on ice

Add 400ul 1B

incubate 37°C 1 hr

Plate 200ul onto 150 mm LB+Amp

incubate 37°C O/N

12/3/94

Take plates out of 37°C

Put at 4°C

12/6/94

Pick into 200ul LB+Amp

incubate 37°C 2/N w/antibiotic

PCR To test for inserts

12/7/94

51bp		40X
5' Bam SI	0.2	8
3' Asp	0.2	8
10x dNTP	3.2	128
10x PCR	3.2	128
Taq	0.2	85
H <sub>2</sub> O	21	910
Cult	2	—

185bp		40X
5' Bam 185	0.2	4
3' Asp	0.2	4
10x dNTP	3.2	64
10x PCR	3.2	64
Taq	0.2	4
H <sub>2</sub> O	21	420
Cult	2	—

68 HTPAN08504 51bp/185bp in PAZ

12/31/92


PCR prog # 66

95°C	5min	} 30x
95°C	20sec	
55°C	20sec	
72°C	1min	
72°C	7 1/2 min	
4°C	Hold	

Run Out on 1% gel with 1 kb ladder.

HTPAN08504 51bp

HTPAN08504 185bp



Pick @ HTPAN08504 185bp into 5ml  
TB+Amp  
Incubate 37°C O/N w/aeration

12/8/94

Boiling min prep

Spin 2ml culture 5min  
~~Remove~~ Remove Supernatant  
Resuspend pellet 750ul STE + RNase +  
lysozyme.

(pg 84)




HTP4W08 PAZ

89

pg 68

D/8/94

Boil 1 min  
Spin 10 min  
Remove pellet  
Add 750  $\mu$ l PEG/NaCl  
Mix well  
Spin 10 min  
Remove supernatant  
Add 1000  $\mu$ l 70% Ethanol  
Wash pellet  
Spin 5 min  
Remove supernatant  
Allow pellet to dry at RT 10 min.  
Resuspend pellet 200  $\mu$ l  $\phi$ TE  
Run 2  $\mu$ l on gel with 1 kb ladder +  $\lambda$  Hae



Substrate for

Try digesting with Bam/Asp

			19x
Asp	DNA	4	429.4
	H <sub>2</sub> O	22.6	22.6
	10x <sup>+</sup> 2	3	3.6
	Bam	0.2	3.6
	Asp	0.2	3.6

Incubate 37C O/A

16. APR 94  
12/2/94

90

HTPAND08 PA2 / PD10

12/9/94

Does not look like anything digested  
 Submit for Sequencing  
 RP06  
 RP19

11/11/95

Remake inserts w/ new primers  
 for BlnA PA2 & PD10

PA2 has Bam ~~HT~~ HT site with  
 Kozak seq.

51bp + 185 bp fragments

PCR ③

PD10

HTPAND08 51bp

9113  
 2742  
 10xPCR  
 10xDNTP  
 1µl  
 H<sub>2</sub>O  
 DNA

20  
 2  
 50  
 50  
 2  
 375  
 1

④

PD10

HTPAND08 185bp

9114 15  
 2742 2  
 50  
 50  
 2  
 382  
 1

①

PA2

HTPAND08 51bp

9111 18  
 2742 2  
 50  
 50  
 2  
 377  
 1

②

PA2

HTPAND08 185bp

9112 15  
 2742 2  
 50  
 50  
 2  
 382  
 1

DNA-HTPAND0804 plasmid DNA - 10 µg/µl.  
 PCR using Program # 66

③ Control  
 H<sub>2</sub>O only

PA2 125

PQE60 Nco / Bam Vector

97

12/19/94


Digest PQE60 = 0.54  $\mu$ g.

DNA	4 <del>10</del> $\mu$ l
10x #4	5
H <sub>2</sub> O	40
Bam	0.5
Nco	0.5
	<hr/> 50


Incubate 37°C O/N

12/20/94

Run 2  $\mu$ l on gel with 1 kb ladder



looks good.  
Add 5  $\mu$ l 10x SB  
Run on 0.8% TBE gel  
80V 1 1/2 hr  
Cut out fragment.  
Take picture



Gene Clean.

Add 900  $\mu$ l NaI  
Heat 55°C 5 min  
Mix well.  
Add 7  $\mu$ l Glass beads  
Mix well  
Incubate at RT 5 min  
with 0.1% Mucosin  
Spin 10 sec  
Remove Supernatant  
Resuspend pellet 500  $\mu$ l  
Wash Buffer.  
Spin 10 sec

3A

W

HTPAN08/HTPB411 in PGE60

99

12/20/14

PCR

HTPAN08504 3' NcoI / 3' Bgl II  
 HTPB411S15 5' BspHI / BamHI  
 for ligation into PGE60

HTPAN		HTPB4	
5' NcoI	20	5' BspHI (5bp)	20
3' Bgl II	4	3' Bam	4
10x dNTP	30	10x dNTP	30
10x PCR	30	10x PCR	30
H <sub>2</sub> O	234	H <sub>2</sub> O	234
Taq	1	Taq	1
DNA (100ng/ul)	1	DNA (100ng/ul)	1
	300		300

Set up rxn w/o DNA (30ul vol total)


PCR. B10x 30ul rxns.

Prog 66.

Run 5ul on gel with 1kb ladder

5  
 5  
 39ul  
 0.5  
 0.5  
 50ul

95°C	5min	30x
95°C	20sec	
65°C	20sec	
72°C	1min	
72°C	7 1/2 min	
4°C	hold	




Add equal Vol PEG/NaCl to ppt. - mix well  
 Spin 10min  
 Remove Supernatant  
 Wash 1x 70% Ethanol (100ul)  
 Spin 5min  
 Remove Supernatant.  
 Allow pellet to dry at RT 10min  
 Resuspend in 100ul TE

100

HTPAN08 HTPB411 in PQE100

12/20/94

Run 50  $\mu$ l on 0.8% LMP gel  
80 V 1 1/2 hrs

Cut out fragments

Gene clean  
Resuspended - Total 40  $\mu$ l TERun 2  $\mu$ l on gel with  
1 kb ladder

looks good.

Set up digestions:

DNA	30 $\mu$ l
H <sub>2</sub> O	13
10x#4	5
Enz/Enz	0/1
	30 $\mu$ l

For	HTPAN08	-	Nco	/	Bgl	II
	HTPB411		Bsp	HI	/	Bam

Incubate 37°C O/N.


12/22/94

Add 5  $\mu$ l 10x SB  
Run on 0.8% LMP gel with  
1 kb ladder  
80 V 1 1/2 hrs

HTPB008 | HTPB411 in PQE60

101


12/22/9



Gene Clean

Resuspend in 40 ul Total

Run 2 ul on gel with 1 Kb ladder



Set up ligation

	①	②	③	④	⑤	⑥
HTPB008 5' Nco / 3' Bgl II	4	—	4	—	—	—
HTPB411 5' 5'p H3 / 3' Bam	—	4	—	4	—	—
10x	2	2	2	2	2	2
H <sub>2</sub> O	12	12	13	13	16	1
T4	1	1	1	1	1	1
PQE 60 Nco / Bam	—	1	—	—	—	—
PQE 60 Nco / Bgl II	1	—	1	—	—	—

Incubate 4°C

1/3/95

Transform M15 supH cells  
Chemically competent (long)

Thaw cells on ice  
To 100 ul of cells add 10 ul of ligation  
To 1 tube add 10 ul PQE60 only  
To 1 tube add 10 ul of H<sub>2</sub>O only

1/3/95

Incubate on ice 1 hr  
 Heat  $42^{\circ}\text{C}$  1 min  
 Let sit on ice 2 min  
 Add 400ul LB  
 Incubate  $37^{\circ}\text{C}$  1 hr.  
 Plate onto 2 Bt Amp + Kan plates.

1A - 200ul, 300ul  
 1B - 100ul, 400ul  
 2A - 200ul, 300ul  
 2B - 100ul, 400ul  
 3 - 300ul  
 4 - 300ul  
 5 - 300ul  
 6 - 300ul ligation Rxn.  
 7 - PDE60 10ng 100ul  
 8 - M15 cells alone 300ul

Spread Evenly with Beads  
 Incubate  $37^{\circ}\text{C}$  O/N.

1/4/95

All plates grew -  
 Even M15 cells alone and ligation Rxn  
 alone  
 Cells are Contaminated!

Told Lang - said he would test cells  
 & make New.

Try Redoing ligation with  
 HTPANOS 504 185 bp fragment.  
 (Pg 107)

HTPAN08504 185bp Fragment

103

10/28/94

PCR HTPAN08504 185bp fragment  
for PQE10 Nco I Bgl II

DNA	1 $\mu$ l
H <sub>2</sub> O	237
10x PCR	30
100x dNTP	30
1 $\mu$ g	1
5' Nco I 185bp	20
3' Bgl II	4
	300 $\mu$ l

100x H<sub>2</sub>O only  
do @ Prep.

Split 10 tubes

PCR Program #101.

95°C	5 min
95°C	20 sec
55°C	20 sec
72°C	1 min
72°C	7.5 min
4°C	hold.

30x

done 4°C

Run  
PCR (10  $\mu$ l each) 1 Kb ladder

1/3/94



Add equal Volume 13% PEG/  
1.2M NaCl  
incubate on ice 10 min  
mixing well  
spin 10 min  
pour off supernatant.



104

HTRANO8504

185bp Fragment

1/3/95

Wash pellet in 700ul 70% Ethanol.

Spin 5 min.

Pour off supernatant

Air dry pellet to dry at RT

Resuspend pellet in 200ul TE

Now Set up Digestion


DNA	35
H <sub>2</sub> O	8
10X#4	5
Nco	1
Bgl II	1
	<hr/> 50ul

Incubate 37°C O/N.

1/4/95

Run EMP Gel - 0.8% in TAE.

80V 1 1/2 hrs

Cut out fragment from gel and place in tube  
Take picture of gel.

2 Tubes

1 Tube use 10M NaI from Bio 101

Other Tube use 10M NaI made here

Add 1000ul NaI  
Heat 55°C 5 min

pg 107

~~HTLAD31~~ 1986 bp Fragment

105

HTLAD31 PER Protein

1/2/95

PER - Protein

Digest HTLAD31 = 0.25 µg/µl

- (1) Xho / Eco RI  
(2) Sac / Kpn

(1)  
DNA 6  
H<sub>2</sub>O 20.6  
10X #2 3  
Xho 0.2  
Eco RI 0.2  
30 µl

(2)  
DNA 6  
H<sub>2</sub>O 20.6  
10X #4 3  
Kpn 0.2  
Sac 0.2  
30 µl

Incubate 37°C O/N

Submit for sequencing

1/4/95

Run 10 µl in 1% TAE gel  
with 1 Kb ladder

1- Xho / Eco RI ~ 4.7 Kb

2- Kpn / Sac ~ 800 bp + ~ 600 bp

Insert looks ~ 1.6 Kb

make primers -

RPD6 & PPD5

1/9/95 Received RPO6 & RPO4.

Submit for sequencing  
with other RPO2, RPO3  
and Reverse & Forward

1/11/95 Received FPO5.

Submit for seq RPO6 & FPO5

1/13/95

large 3' untranslated region  
~ 1.5 Kb long

about 1.2 Kb is untranslated

HTPAV08 SOL 185bp Fragment

107

pg 104

+ pg 102

1/4/95

Mix well to make sure all of gel  
is dissolved.

Add 8ul Glass milk + mix.

Incubate at RT 2min w/ occasional  
mixing

Spin 10 sec

Pour off Supernatant

Resuspend pellet in 40ul Wash Buffer

3 Spin 10 sec

Pour off Supernatant

Spin 10 sec

Remove as much of Supernatant as  
possible

Resuspend pellet 30ul TE

Heat 55°C 2min

Spin 10 sec

Transfer Supernatant to fresh Tube

Resuspend pellet in 20ul TE

Heat 55°C 2min

Spin 10 min sec

Transfer to tube

Run 2ul on gel with 1 kb ladder

1- using Gene Clean Kit  
6A NaI

2- using Home made 6A NaI

Both look good...  
Ready for ligations

1/5/95

Set up ligation for

HTPB08 & HTPB411 in PDE60

1/5/95

	①	②	③	④	⑤	⑥	⑦	⑧	⑨
HTPANOS 5' S' Aco / 3' Bgl II	6			6	6				
HTPANOS 5' S' Aco / 3' Bgl II		6			6				
HTPBYS 5' S' Bsp HI / 3' Bam HI			6			6			
10X Buffer	2	2	2	2	2	2	2	2	2
PGE60 Nco / Bgl II	2	2						2	
PGE60 Nco / Bam HI			2						
H <sub>2</sub> O	9	9	9	11	11	11	15	15	17
T4 DNA Ligase	1	1	1	1	1	1	1	1	1

Make Cocktail

		9x
10X Buffer	2	18
H <sub>2</sub> O	9	81
T4 DNA lig	1	9
	12	

Store ligations at 4°C till  
fused M15 cells are found

1/9/95

STEVE said to transform into XL-1  
Blue plasmid is correct - ~~it~~ can  
transform into M15 cells.

Thaw XL-1 Blue Chemically Competent  
Cells on ice  
To 100 µl of ligation add 100 µl  
thawed cells  
Incubate on ice 1 hr - as ② controls  
Add 100 µl cells only  
Heat Shock 42°C 45 sec

P115

Maxi Prep HTPAN08504/HMEAA88

111

1/4/95

Inoculate 200ml TB + Amp  
from frozen stocks of  
HTPAN08504  
HMEAA88A.

Cultivate 37°C w/ aeration  
overnight.

1/5/95

Diagen Maxi Prep

Spin Cultures 4.5K 15min

Pour off supernatant  
Resuspend pellet in 10ml P1 +  
RNase

Cultivate RT 5min

Add 10ml P2

Cultivate RT 10min

Add 10ml P3 - ice cold

Cultivate on ice 20min

Spin 4.5K 15min

Equilibrate Tip -500 with 10ml  
QBT

Apply Supernatant to tip through  
Kim wipe

Allow to flow through

Wash Column 3x

30ml QC

Elute DNA - 15ml QF

Add 10.5ml (0.7 times volume)

Mix well

Spin 8K 30min

Pour off supernatant

Wash pellet 15ml Cold 70%  
ethanol


112

Maxi Prep HIRAND8504, HMEAA88

1/5/95

Spin 8K 15min  
 Pour off supernatant.  
 Allow off pellet to dry -  
 Dry o/n at Room Temp -  
 with parafilm w/ holes covering.

1/6/95

Resuspend pellet in 400  $\mu$ l TERun 1  $\mu$ l on gel with 1 Kb ladder  
 & 2 Hind III MarkerDilute 1:100 in H<sub>2</sub>O  
 Read OD<sub>260</sub>/280

Sample ID	abs	abs	bkg abs	260.0 nm	280.0 nm	
	260.0 nm	280.0 nm	320.0 nm	280.0 nm	260.0 nm	
1 HIRAND85	0.0402	0.0319	0.0231	1.9374	0.5161	0.2 $\mu$ g/ $\mu$ l
2 HIRAND85	0.2051	0.1314	0.0315	1.7389	0.5754	1.03 $\mu$ g/ $\mu$ l

looks like HMEAA88 is all chromosomal

1/9/95

Steve said to try transforming with  
 Boring pup DNA to get a  
 clone

~~HTPB411~~ ~~HTPB411~~ HTPAN08 HTPB411 P115

(Pg 108)

1/9/95

Set on ice 2 min  
Add 400ul LB  
Incubate 37°C 1 1/2 hrs  
Plate 200ul onto LB + Amp  
150mm plates

Transform ligations from 12/22  
too.

Incubate 37°C O/N  
Plate M15 Cells on LB + Amp  
LB + Kan  
LB + Amp/Kan

1/10/95

Plates look good.  
No contamination

M15 Cells

Grew on LB + Kan  
on LB + Amp + Kan  
Not on LB + Amp.

∴ LB + Amp/Kan plates don't have  
Amp or enough Amp so contamination  
can be seen.

Can use M15 cells - Not contaminated

Picks - ~~Q1~~ Colonies onto LB + Amp in  
96 well dish.

HTPB411

HTPAN08.

12/22 2A - 36  
2B - 44  
1/5 3 - 25

12/22 2 48  
1/5 1 48  
1/5 2 48

pick PRELID Vector + PA2 Vector as controls



116

HTPAN08/HTPB411 in PDE60

1/10/95

Incubate 37°C w/aeration 5hrs  
 PCR using internal primer  
 as well as the 2 external primers

HTPB411 + PDE60		HTPAN08 51bp + PDE60		HTPAN08 185bp + PDE60	
FR20	0.5	FR14	0.2	FR4	0.2
2887	0.2	2888	0.2	2867	0.2
2886	0.2	2885	0.2	2865	0.2
10x dNTP	3.2	10x dNTP	3.2	10x dNTP	3.2
10x PCR	3.2	10x PCR	3.2	10x PCR	3.2
H <sub>2</sub> O	18.7	H <sub>2</sub> O	17	H <sub>2</sub> O	28.8
Taq	0.2	Taq	0.2	Taq	0.2
Culture	2	Culture	2	Cult	2.0
	32		32		32

PCR - (+) control Plasmid DNA  
 (-) control LB + Amp alone

PCR Program

95°C 5min  
 95°C 20sec  
 55°C 20sec } 30x  
 72°C 1min  
 72°C 7 1/2 min  
 4°C Hold

Set up Rxns  
 Let Cultures  
 grow 37°C O/N  
 w/aeration

1/11/95

Run 10 µl of PCR Rxn on a 1.3% Agarose  
 gel in TAE with 1 kb ladder

89/117

HTPB411 in PAZ / HTPANOS + HTPB411 in PGE60 119





pg 110

pg 116

1/11/95




① HTPB411 + PAZ A-E  
HTPB411 + PGE60 H

② HTPB411 + PGE60 A-E  
HTPB411 + PAZ F-H



③ HTPANOS 51bp + PGE60

④ HTPANOS 185bp + PGE60



120

HTPB411 PAZ / HTPB411 + TRAN in PDE60

1/10/95

Inoculate ~~28~~ TB + Amp 5ml  $\phi$   
do Baling pups.

① HTPB411 + PAZ

1, 5, 7, 8, 12, 16, 17, 18, 22, 23, 24, 29, 30, 31, 32, 41, 42, 43

② HTPB411 + PDE60

2, 6, 7, 19, 23, 24, 29, 38, 41, 45, 48, 57, 58, 66, 69, 73, 77, 79

③ HTPA085.6p + PDE60

1, 2, 3, 5, 6, 7, 8, 9, 10, 11, 14, 21, 23, 24, 27, 31, 33, 34

④ HTPA08185.6p + PDE60

5, 6, 7, 8, 11, 13, 15, 16, 17, 20, 21, 22, 23, 24, 25, 26, 28, 32

Incubate 37°C w/ aeration overnight

1/12/95

Boiling min Pups

Spin 2ml of culture 2min

Remove Supernatant

Resuspend pellet in 750  $\mu$ l STE +  
RNase + Lysozyme

Boil 1min

Spin 10 min

Remove Pellet

Add equal Volume - 750  $\mu$ l of 13% PEG 8000  
1.0 M NaCl

Mix well

Spin 10 min

Remove Supernatant

Wash Pellet 1000  $\mu$ l 70% Ethanol



Mix

HIPB411 PAZ / HIPB411 + HIPANOS POE60. 121

1/12/95

Spin 5 min.  
Remove ethanol  
Allow pellets to dry at RT for 30 min  
Resuspend pellets in 100  $\mu$ l TE

Run 2  $\mu$ l of minipreps on 1% TAE  
gel with 1 kb ladder + PAZ & POE60



mix up in labeling - #1 is really #3  
#2 is really #4  
#3 is really #1  
#4 is really #2

Set up digests

For POE60.

EcoRI / HII

55X

DNA

4

10X

2

H<sub>2</sub>O

13.6

EcoRI

0.2

HII

0.2

For PAZ

Bam / Xba

DNA

4

10X

2

H<sub>2</sub>O

13.6

Bam

0.2

Xba

0.2

Incubate 37°C O/N

Digest

HIPANOS504

HIPB411S15

PAZ + POE60

122


HTPB4H P22 | HTPAN08 + HTPB411 PDE60

1/13/95

1,2,4 - Eco / HII

3, Bam / Xba digested

Run 10ul on gel with 1 kb ladder



looks like

4- HTPAN08 185bp + PDE60

1-17 correct - Sequence 4.

2. HTPB411 in PDE60

1-4,6,8-16 looks good Seq 4.

14,3 Need to be digested w/correct enzymes

Submit 4x2 (4 of each)  
 #4 with PDE60 & P16  
 #2 P23  
 for sequencing

1/16/95

Cleanup DNA of 4-1, 4-2, 4-3, 4-4  
 2-1, 2-2, 2-3, 2-4

2x Phenol

2x 0.5% SDS

Elute into PPT - Wash

Resuspend in 80 ul TE - use to transform 15 cells

pg 227

HTPAND8 PAZ/PD10

125

pg 10

95°C 5min  
 95°C 20sec  
 55°C 20sec } 30V  
 72°C 1min  
 72°C 7.5min  
 74°C Hold

Run 5µl on gel with 116 ladder

1/12/95

3 - HTPAND8 51 bp for PAZ  
 4 - HTPAND8 185 bp for PAZ  
 1 - HTPAND8 51 bp for PD10  
 2 - HTPAND8 185 bp for PD10



1/12/95

2

Precipitate Add equal Volume  
 3% PEG/NaCl

Spin 10 min

Wash pellet 1000µl 70% Ethanol

Spin 5 min

Remove Supernatant

Set up digestion

DNA (PCR Fragment) 20

10x #2 Buffer 4

H<sub>2</sub>O 15

Bam 0.5

Xba 0.5


40µl

Incubate 37°C overnight

1/13/95

Run All of digestion of 0.8% Agarose

Cut out gel fragments  
Take picture




- 1 - HTPAN0851bp } PD10
- 2 - HTPAN08185bp } PD10
- 3 - HTPAN0851bp } PAZ
- 4 - HTPAN08185bp } PAZ

Gene Clean  
Add 800  $\mu$ l NaI  
Heat 55°C 5 min  
Add 8  $\mu$ l Glass Milk  
Mix well  
Incubate at RT 5 min  
w/ occasional mixing  
Spin 10800  
Remove Supernatant

\* Add 500  $\mu$ l Wash Buffer  
Resuspend pellet  
Spin 10800  
Remove Supernatant

Spin 10800  
Remove Supernatant  
Resuspend pellet in 20  $\mu$ l TE  
Heat 55°C 1 hr  
Spin 10800  
Transfer Supernatant to fresh tube

Run on gel with 1 kb ladder and 1  $\mu$ l  
of PAZ and 1  $\mu$ l PD10 B/x



- 1 HTPAN08 51bp PAZ
- 2 HTPAN08 185bp PAZ
- 3 PAZ Bam Hba
- 4 HTPAN08 51bp PD10
- 5 HTPAN08 185bp PD10
- 6 PD10 Bam Hba

pg 122

HTPANO8/HTPB41 PAZ/PD10 127

Set-up ligations

1/13/95

	1	2	3	4	5	6	7	8	9	10	11	12	13
HTPANO8 51 PAZ	6					6							
HTPANO8 185 PAZ		6					6						
HTPANO8 210 PD10			6					6					
HTPANO8 185 PD10				6					6				
HTPB41 PD10					6					6			
10X Buffer	2	2	2	2	2	2	2	2	2	2	2	2	2
H <sub>2</sub> O	9	9	9	9	9	11	11	11	11	11	15	15	17
T4 Ligase 100ul	1	1	1	1	1	1	1	1	1	1	1	1	1
PAZ Bam/xba	2	2									2		
PD10 Bam/xba			2	2	2							2	

Set up ligations on ice

Incubate 16°C over the weekend

1/14/95

Then M15 Chemically Competent cells

To 100ul of Cells add 10ul of  
ligation (from 1/13/95) and 1ul  
of DNA (from pg 122)

DI Tube of M15 cells only / 1 of 10ng PD10

Incubate on ice 1 hr

Heat 42°C 1:45 sec

Place on ice

Add 400ul LB

Incubate 37°C 1 1/2 hrs

Plate onto LB + Amp 150mm plates

- ligations - 1-5 plate 100ul + 30ul

- plasmid DNA - plate 50ul

- ligation controls 1-3 plate 30ul

- M15 cells plate 30ul

- M15 cells + PD10 plate 100ul

Incubate 37°C O/N @



## H1PB411 + H1TPAN08

1/16/95

H1PB411 + PA2

1/10/95 Transformalins

Pick 96 Colonies into LB + Amp

H1TPAN08 51bp + PDE60

1/10/95 Transformalins

48 from 1/5/95 ligation

48 from 2/22/94 ligation

into 200ul LB + Amp

Incubate 37°C w/ aeration o/n.

1/17/95

Transformations worked well -

No Colonies seen on  
fragment only plates2-125 colonies seen on  
Vector alone plateNo Colonies on - MIS cells alone  
- leg. R. alone.

Pick colonies into 200ul LB + Amp + Kan

Plate 1) ① H1TPAN08 51bp + PA2 - 48 A17D12

② H1TPAN08 185bp + PA2 - 48 E12H12

12) : ③ H1TPAN08 51bp + PD10 - 48 A17D12

④ H1TPAN08 185bp + PD10 - 48 E12H12

13) : ⑤ H1PB411 + PA2 PD10 - 16 A17B4

PA2 PD10 4

PA2 4

14) 2-1 (12) A1-A12

2-2 (3) H1PB411 + PDE60 B1-B3

2-B (12) in MIS cells C1-C12

2-4 (12) in MIS cells D1-D12

# HTPB411 + HTPANOS

129

1/17/95

E1-E12 4-1 (12)  
F1-F12 4-2 (12)  
G1-G12 4-3 (12)  
H1-H12 4-4 (12) } HTPANOS 185bp +  
PQEG60 M15 cells.

Incubate at 37°C w/ airtight 4 hrs.

PCR - HTPB411 + PA2 + HTPANOS 51bp + PQEG60

HTPANOS 51bp + PQEG60.

		100X
2888	2	200
2409 FPI6.	0.15	15
10x dNTP	3.2	320
10x PCR	3.2	320
H <sub>2</sub> O	2125	2125
Temp	0.2	20
Cult	2	
	32ul	30ul/tube

HTPB411 + PA2

100X


2741	0.2	20
2796 RPO6.	1.2	120
10x dNTP	3.2	320
10x PCR	3.2	320
H <sub>2</sub> O	22	2200
Temp	0.2	20
Cult	2	
	32ul	30ul/tube

PCR Prog. Col.

95°C 5min  
95°C 20sec  
55°C 20sec  
72°C 1min  
72°C 7 1/2 min  
4°C Hold

Run 10ul on gel with 1 kb ladder.

Control use 1 B Kestrel only



HTPB411 + PA2  
2741 (3' xba-) + 2796 (RPO6)

1-12  
1-12

71012  
2410  
784



-A12  
-B3  
-C12  
-D12

30



H7PB411 + H7PAN08

1/17/95

H7PAN08 51bp + PDE60 (2888 5' BspH1 + 7409 (FP16) S)



H7PB411 + PA2  
2741 (3' Xba) + 7796 (R66)



H7PAN08 51bp + PDE60  
2888 (5' BspH1) +  
7409 (FP16)

Inoculate 5ml TB + Amp w/ln.  
① clones 1-9 from H7PB411 + PA2  
H7PAN08 51bp + PDE60

Incubate 37°C w/aeration O/N

PCR Plates [1], [2], [3], [4]

# HTPB411 & HTPAN08

131

## Plate 11

11/17/95

### ① HTPAN08 51bp + PAZ

		55x
9111	1.7	66
2742	0.1	5.5
10x dNTP	3.2	174
10x PCR	3.2	174
H <sub>2</sub> O	22.1	1215.5
Tag	0.2	11
Cult.	2	
	32ul	30ul/tube

### ② HTPAN08 185bp + PAZ

		55x
9112	1.7	66
2742	0.1	5.5
10x dNTP	3.2	174
10x PCR	3.2	174
H <sub>2</sub> O	22.1	1215.5
Tag	0.2	11
Cult.	2	
	32ul	30ul/tube

control  
HTPB411

③ Control PAZ Vector + LB Amp/Km

## Plate 12

### ③ HTPAN08 51bp + PD10

		55x
9113	1.4	77
2742	0.1	5.5
10x dNTP	3.2	174
10x PCR	3.2	174
H <sub>2</sub> O	21.9	1204.5
Tag	0.2	11
Cult.	2	
	32	30ul/tube

control  
HTPB411

### ④ HTPAN08 185bp + PD10

		55x
9114	1.2	66
2742	0.1	5.5
10x dNTP	3.2	174
10x PCR	3.2	174
H <sub>2</sub> O	22.1	1215.5
Tag	0.2	11
Cult.	2	
	32	30ul/tube

control  
HTPB411

⑤ Control PD10 Vector + LB Amp/Km

## Plate 13

### ⑤ HTPB411 + PD10

		24x
2752	1.2	28.8
2741	0.2	4.8
10x dNTP	3.2	76.8
10x PCR	3.2	76.8
H <sub>2</sub> O	22	528
Tag	0.2	4.8
Cult.	2	
	32ul	30ul/tube

control  
PD10 Vector  
+ LB Amp/Km

control  
HTPB411

1/17/95

plate 4

① HTPBULL + PDE60

		SDX
28057	3	150
28056	0.2	110
10x dNA	3.2	160
10x PCR	3.2	160
H <sub>2</sub> O	20.2	1010
Temp	0.2	10
Cult	2	
	32	32 / tube

② HTPANOB185 + PDE60

		SDX
29007	0.2	10
28057	0.2	10
10x dNA	3.2	160
10x PCR	3.2	160
H <sub>2</sub> O	23	1150
Temp	0.2	10
Cult	2	
	32	32 / Tube

③ control HTPBULL

④ control HTPANOB

⑤ control LB Amp / Km

Per Prog 66-

Run 10ul on gel with 1k ladder



95°C 5min  
 95°C 20sec  
 55°C 20sec  
 72°C 1min  
 72°C 7 1/2 min  
 4°C hold

BOX

Run by plate

Plate 1

Plate 2



HTPB411 & HTPAN

183

plate [3] + plate [2]

plate [2]

11/17/95





plate [2]

plate [3]







plate [3] + plate [4]



plate [4]



11/17/95

Plate [4]

Controls:



Controls		
① 9111 + 2742	HTPAN08	51bp PAZ
② 9112 + 2742	HTPAN08	185bp PAZ
③ 9113 + 2742	↓	51bp PD10
④ 9114 + 2742		185bp PD10
⑤ 2752 + 2741	HTPB411	PD10
⑥ 2887 + 2886	HTPB411	PQ E60
⑦ 2967 + 2865	HTPAN08	185bp PQ E60

Inoculate 5ml LB + Kan/Amp  
with 2 colonies ea.

With 2 clones each of  
HTPB411 + PQE60 - 2-1, 2-2, 2-3, 2-4 (A+B)  
HTPAN08 185bp + PQE60 - 4-1, 4-2, 4-3, 4-4 (A+B)

Inoculate 5ml TB + Amp  
with 9 clones of each.

HTPAN08	51bp + PAZ	- 9
HTPAN08	185bp + PAZ	- 9
HTPAN08	51bp + PD10	- 3
HTPAN08	185bp + PD10	- 9
<del>HTPB411</del>	<del>PQE60</del>	
HTPB411	PD10	- 9

Incubate 37°C O/D w/ Auralum

HIPB411 & HIPANOS8

135

1/18/95

HIPB411 + PQE60.

HIPANOS8 185bp + PQE60.

To 5ml LB + Amp/Kan add 300  $\mu$ l  
of ON cultures.  
Incubate 37°C w/ aeration till  
OD<sub>600</sub> = 0.4 - 0.6. - 2 hrs  
Add 100mM IPTG to 2mM + 110  $\mu$ l.  
Incubate 37°C 4 hrs  
Spin 750  $\mu$ l of culture - induced  
and 1 of uninduced of HIPANOS8  
HIPB411.

Remove supernatant

Resuspend pellet in 100  $\mu$ l H<sub>2</sub>O

Add 100  $\mu$ l 2x dissociation Buffer.

2x  
Dissociation Buffer: 0.25M Tris pH 6.9  
4% SDS  
20% Glycerol.  
10%  $\beta$ -Mercapto Ethanol.  
0.2% Bromophenol Blue.

Mix well

Heat 100°C 5min

Spin 5min

Put on ice

Run 10  $\mu$ l on 12.5% Acrylamide  
stacking gel with Rainbow  
marker

Run at 150V in 1x Running  
Buffer for 1 1/2 hrs. - till dye  
front is at bottom of gel.

Stain - 1/2 hr at 37°C

DeStain overnight at RT -

Boiling Miniprep.

A45)  
A45)



1/18/95

HTPB411-LPAZ  
 HTPAN08 51bp PAZ  
 HTPAN08 185bp PAZ  
 HTPAN08 51bp PD10  
 HTPAN08 185bp PD10  
 HTPAN08 51bp PDE } Boiling mixture

Some Cultures did Not grow

Open 2ml Culture

Remove Supernatant

Resuspend pellet in 750ul STET  
 + RNaseB / Lysozyme

Boil 1 min

Spin 10 min

Remove Pellet

Add 750ul 13% PEG 8000 / 16HNaCl

Mix Well

Spin 10 min

Remove Supernatant

Wash pellet 1000ul 40% Etanol


Spin 5 min

Remove Supernatant

Allow Pellet to dry at RT 30 min

Resuspend pellet in 150ul TE

Run 2ul on gel with 1 Kb ladder



A - HTPB411 + PAZ

B - HTPAN08 51 + PAZ

C - HTPAN08 185 + PAZ

D - HTPAN08 51 + PDE

E - HTPAN08 185 + PD10

F - HTPAN08 51 + PD10

Setup digests

H7PB411 & H7P0N08.

137

for PAZ & PD10  
Constructs  
Digest: Bam/Xba

		40X
DNA	5	—
10X#2	3	120
H <sub>2</sub> O	21.6	874
Bam	0.2	8
Xba	0.2	8
	30ul	25ul/tube

For PGE Constructs  
Digest: EcoRI/HindIII

		10X
DNA	5	—
10X#2	3	30
H <sub>2</sub> O	21.6	216
EcoRI	0.2	2
HindIII	0.2	2
	30ul	25ul/tube

Incubate all reactions at 37°C O/N.

Digestions of Minipreps:

Run 10ul of Digestion on 1% TAE  
Agarose gel with 1 Kb ladder



A H7PB411 + PAZ  
B H7P0N08 51bp + PAZ  
C 185bp + PAZ  
D 51bp + PD10  
E 185bp + PD10  
F 51bp + PGE100

looks like some digested correctly

clean-up and send in to be sequenced with internal primers

1/19/95

Dilute HTPB411 + ~~PAZ~~ PAZ 1, 2, 3, 4  
 HTPAN08 51bp + PAZ 1, 2, 3, 4  
 HTPAN08 185bp + PAZ 1, 3, 4, 5

1:200 - use 10  $\mu$ l to transform into DH5 $\alpha$ .

Thaw DH5 $\alpha$  Chemically Competent cells  
 on ice  
 To 100  $\mu$ l of thawed cells Add 10  $\mu$ l  
 of Diluted DNA.  
 incubate on ice 1 hr  
 Heat 42°C - 45 sec  
 place on ice  
 Add 400  $\mu$ l LB  
 incubate 37°C 1 hr  
 plate 300  $\mu$ l onto 150 mm LB + Amp  
 plates  
 incubate 37°C O/N.

PCR HTPB411 + PAZ  
 HTPAN08 51bp + PAZ  
 HTPAN08 185bp + PAZ  
 with a PAZ specific T7 promoter  
 primer. and the 3' end primer

HTPAN	
T7	0.2
2742	0.2
10X	10
10X	10
H <sub>2</sub> O	78.2
Taq	0.4
DNA	1 $\mu$ l of Diluted
	100 $\mu$ l

HTPB411	
T7	0.2
2741	0.2
10X	10
10X	10
H <sub>2</sub> O	78.2
Taq	0.4
DNA	1 of diluted
	100 $\mu$ l

4 Tubes of each

HTPAN08 / HTPBY11

139

Run Program #66.

1/19/95

95°C	5min	} 30x
95°C	20sec	
55°C	20sec	
72°C	1min	
72°C	7 1/2 min	
4°C	Hold	

Run 5ul of PCR on 1% TAE gel with 1Kb ladder



HTPB411 Did Not  
PCR - Need to  
get more  
clones.

HTPAN08 look  
good.

HTPAN08 - precipitate with equal Vol  
13% PEG / 16M NaCl  
mix well

Spin 10min

Pour off supernatant

Wash pellet 1ml

70% ethanol


Spin 5min

Pour off supernatant

Allow pellet to dry at RT 15min.

Resuspend in 100ul TE

Run 2ul on 1% gel with 1Kb ladder



Don + Know about  
upper band 2Kb - but  
band looks good  
Do TNT 1/20.

1/19/95

Clean up Borling pups to submit  
for sequencing

Add TE to GS all samples to 400ul.  
2X Phenol / SEVAG (1:1)

Extract with Equal Volume  
2X SEVAG Extract with Equal  
Volume

Add 1/10 vol (40ul) 3M NaOAcate pH 5.5  
2 vol (800ul) 100% ethanol  
Mix well

Sit on ice 10min

Spin 15min

Remove supernatant

Wash pellet 100ul 70% ethanol


Spin 5min

Remove supernatant

allow pellet to dry at RT 5min

Resuspend pellet in 100ul TE

Run Sinton gel with 1kb ladder



(A) 1-4 HTPAN08 516bp +  
PAZ

(B) 5-8 HTPAN08 185bp +  
PAZ

(C) 9-12 HTPAN08 51bp +  
PD10

(D) 13-16 HTPAN08 185bp +  
PD10

(E) 17-19 HTPAN08 51bp +  
PQED

Submit for Sequencing  
FP16 + RPO6.

IRIS Names:

(A) FAS51BP1PAZRP06/FP

(B) FAS185B1PAZRP06/FP

(C) FAS51BP1PD10RPO6/FP

(D) FAS185B1PD10RPO6/FP

(E) FAS51BP1PQEDRP06/FP

FP16

HTPANO8 (HTPB11)

141

1/19/95  
1/24/95

circulate 5ml (BTBmgs) (Kew)  
with HTPB411 + PQE60


HTPANO8 135 + PQE constructs

pg 135

Destained gel

HTPB11 + PQE60

HTPANO8 135 + PQE60



Rem more of sample  
on gel - 10%

looks like #4 induced

stain - Destain

1/20/95



Can't really  
tell  
- try inducing more

PQ6 / 100  
dp / 100  
06 / 100  
rd / 100  
100 / 100  
100 / 100

1/20/95

DO TNT with PCR &  
T<sub>7</sub> fragments (pg 139)

Rabbit Reticulocyte Lysate	12.5
T <sub>7</sub> Buffer	1
T <sub>7</sub> Polymerase	0.5
Amino Acid Omet	0.5
<sup>35</sup> S methionine	2
RNA <sub>in</sub>	0.5
PCR Product	4
H <sub>2</sub> O	4
	<u>25 ul</u>

Incubate 30°C 2 hours

Heat 90°C 5 min

Quick Spin

To 20 ul of 2X Association Buffer  
Add 5 ul of TNTRun 15 ul on 15% Acrylamide gel  
with C<sub>14</sub> labeled Rainbow marker

Run gel 150V 1 1/2 hrs

a until Dry front reaches to bottom

Cut off stacking gel

Fix gel in Fixative

10% H<sub>2</sub>AC, 30% MeOH

for 20 min at 37°C

Run gel fix

Amplify in 30 ml Amplify

for 20 min at 37°C

Dry gel 1 1/2 hrs at 80°C

Place on film

leave at RT over weekend

1/22/95

Incubate 150 ul T<sub>7</sub> Amp / T<sub>7</sub> Amp Km  
with cultures 100°C 0.5 hr



(pg 149)

IL-6 PQECs / PD10

143

pg 70

1/19/95



IL-6 PD10

looks like correct clones

IL-6 PQE					IL-6 PD10				
1-A2	6-B6	11-D12	14-G5		1-A2	6-C4	11-E5	16-G1	
2-A3	7-C1	12-E5	17-G8		2-A7	7-C5	12-E11	17-H2	
3-A7	8-C3	13-E8	18-H5		3-B1	8-C8	13-F3	18-H3	
4-A10	9-C9	14-E11	19-H8		4-B4	9-C9	14-F9	19-H5	
5-A11	10-D10	15-F1	20-F10	7	5-B11	10-D7	15-G8	20-H7	

1/23/95

Inoculate 200ul LB + Amp/Kan  
with cultures 1-20 of each  
Incubate o/n at 37°C w/ aeration to DO  
Mini induction

1/24/95

To 200ul of fresh LB + Amp/Kan  
add 50ul of fresh cultures



HTPB4.11 + HTPAN 308

149

pg 42

Develop 1/23/95

1/23/95

Qiagen Maxi Preps of:

HUEAA88

HTPAN08504 51bp + PAZ #1

+ PAZ #3

+ PAZ #4

185bp + PAZ #1

+ PAZ #3

+ PAZ #4

+ PAZ #5

+ PQE 60

+ PQE 60

4-1

4-2

HTPB411515 + PAZ 1



Make Glycerol stocks of all

800ul Bacterial Stock + 800ul 30% Glycerol

pg42

HTPB411 + HTPAN08

Develop TNT.



Qiagen Maxi Preps of:

HMEAA88

HTPAN08S04 51bp + PAZ #1

+ PAZ #3

+ PAZ #4

185bp + PAZ #1

+ PAZ #3

+ PAZ #4

+ PAZ #5

+ PQE60 4-1

+ PQE60 4-2

HTPB411S15 + PAZ 1

Make Glycerol stocks of all

800ul Bacterial Stock + 800ul 30% Glycerol


HTPB411 + HTPAN08

py42

1/23/95

HTPAN08504  
51bp

HTPAN08504  
188bp



TNT 4/20/95  
PCR of PA2 clones  
HTPAN08504 51bp  
HTPAN08504 188bp  
T2 + 3' Xba primers

15% Gel  
1/31/95  
ON -80C

HTPB411S15 + PA2


+PQE60 4-2  
+PQE60

0.1 0.1 0.1 0.1 0.1 0.1 0.1 0.1 0.1 0.1

PAJ42

HTPB411 + HTPANJ08

1/23/95




HTPB411S15 + PAZ + PÄE60 4-2

Make Glycerol stocks of all  
800ul Bacterial Stock + 800ul 30% Glycerol.

1/23/95

Spin Culture 4.5 K 20 min  
 Pour off Supernatant  
 Resuspend pellet in 10 ml P1 + RNase  
 let sit RT 5 min  
 Add 10 ml P2 - mix by swirling  
~~the~~ looked for cell lyses  
 Add 10 ml P3 -  
 Incubate on ice 30 min  
 Spin 8 K 30 min  
 Prepare Column  
 apply 10 ml QBT  
 allow to flow through  
 Apply Supernatant to Column through  
~~the~~ Kim wipe  
 Allow to flow through  
 Wash Column 30 ml Wash Buffer  
 Elute DNA 15 ml TE  
 Add isopropanol 0.7 vol (10.5 ml)  
 Mix well  
 Spin 8 K 30 min  
 Pour off Supernatant  
 Wash pellet 10 ml 70% Ethanol (-20°C)  
 Spin 8 K 15 min  
 Pour off Supernatant  
 Allow pellet to dry at RT for 10 min  
 Resuspend pellet in 400 µl TE pH 7.6  
 Run 1 µl on gel with 1 kb ladder  
 PA2 & PA2eo

looks good



HTPAN08 / HTPB44

151

1/23/9

Read OD<sub>260</sub>/280  
Dilute 1:200 in H<sub>2</sub>O

Sample ID	abs	abs	bkg abs	260.0 nm	280.0 nm		
	260.0 nm	280.0 nm	320.0 nm	260.0 nm	260.0 nm		400ul
1 HMEABE	0.0466	0.0293	0.0014	1.6230	0.6162	0.4 ug/ml	188 ug
2 HTPAN08 PAZ 10.1271	0.0721	0.0721	-0.0021	1.7411	0.5743	1.3 ug/ml	520 ug
3 HTPAN08 PAZ 30.1188	0.0663	0.0663	-0.0002	1.7859	0.5599	1.2 ug/ml	480 ug
4 HTPAN08 PAZ 40.0983	0.0547	0.0547	-0.0023	1.7855	0.5664	1.0 ug/ml	400 ug
5 HTPAN08 PAZ 60.0919	0.0533	0.0533	-0.0003	1.7201	0.5814	0.9 ug/ml	368 ug
6 " 30.0033	0.0001	0.0001	-0.0053	1.5858	0.6306	X	
7 " 40.1070	0.0630	0.0630	-0.0008	1.6905	0.5915	1.1 ug/ml	440 ug
8 " 50.1304	0.0771	0.0771	0.0035	1.7244	0.5799	1.3 ug/ml	520 ug
9 HTPAN08 PAZ 40.0967	0.0605	0.0605	0.0045	1.6471	0.6071	0.9 ug/ml	358 ug
10 HTPAN08 PAZ 10.0481	0.0287	0.0287	0.0002	1.6798	0.5953	0.48 ug/ml	192 ug
11 HTPAN08 PAZ 10.0502	0.0281	0.0281	-0.0017	1.7434	0.5736	0.50 ug/ml	200 ug

Dilute DNA to 200-250 ng/ml  
Submit for sequencing  
for HTPAN08 - FR16 & R226  
for HMEABE - R201, R202, R204, R205 & R206  
HTPB44

Inoculate 300 ml LB + Amp + Kan  
with 5 ml HTPAN08 504 185 bp + P6560  
O/N culture (41-1)

Inoculate at 37°C w/ aeration  
until OD<sub>600</sub> ~ 0.4-0.6

Add IPTG to 2 mM (1 ml of 100 mM IPTG)  
Inoculate at 37°C w/ aeration  
5 hours

Spin culture 8K 20 min  
Resuspend in 50 ml 6M GdnHCl pH 8  
Store 4°C till tomorrow

Barbara

pg 2  
1/23/9

## HTRANOS 7 HTPBY 11

1/24/95

Inoculate 200  $\mu$ l LBT Amp Kan  
w/  $\oplus$  clones of  $\frac{7}{11}$ 

HTRANOS 51bp + PD10 - 3

HTRANOS 185bp + PD10 - 19

HTRANOS 51bp + PDE60 - 17

} 50  $\mu$ l

Incubate w/ aeration 1 hour - 37°C

Add 100mM IPTG to 2mM  $\rightarrow$  25  $\mu$ l

Incubate 37°C w/ aeration 5 hrs

HTRANOS 51bp + PDE60 - None (w/ 10)

Spin Culture

Add 20  $\mu$ l H<sub>2</sub>O to Resuspend pelletAdd 20  $\mu$ l 2X Dissociation Buffer


Heat 100°C 5 min

Run on 10% Acrylamide gel w/

100 V 1 hr w/ Weight Marker

Stain 50 min

De stain 30 min



Not enough gels to run remaining  
samples, but does not look  
like anything induced

looks like from sequences,  
there is a problem with  
3'bp end -  
Reverse primer

Ag 5  
9  
book  
#27

**This Page is Inserted by IFW Indexing and Scanning  
Operations and is not part of the Official Record**

**BEST AVAILABLE IMAGES**

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

- ☐ BLACK BORDERS
- ☐ IMAGE CUT OFF AT TOP, BOTTOM OR SIDES
- ☐ FADED TEXT OR DRAWING
- ☐ BLURRED OR ILLEGIBLE TEXT OR DRAWING
- ☐ SKEWED/SLANTED IMAGES
- ☐ COLOR OR BLACK AND WHITE PHOTOGRAPHS
- ☐ GRAY SCALE DOCUMENTS
- ☐ LINES OR MARKS ON ORIGINAL DOCUMENT
- ☒ REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY
- ☐ OTHER: \_\_\_\_\_

**IMAGES ARE BEST AVAILABLE COPY.**

**As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.**