



Steven M. Ruben
Appl. No. 10/662,429

Human Genome Project
Lab Notebook # 552

| | |
|--|--|
| Department | <u>Mol. Biol.</u> |
| Subject | <u>10/20/94 - 02/02/95</u> |
| Name | <u>ANN KIM</u> # <u>8</u> |
| Address | <u></u> |
|  43-648 | |
| Computation Notebook | |
| Dennison Stationery Products Co., Framingham, MA 01701 | |
|  | 75 Sheets 11 3/4" x 9 1/2" 4x4 Quad. |
| 0 73333 43648 8 | |

BEST AVAILABLE COPY

Ruben EXHIBIT #91

HTPB411 & HTPAN08 in PD10

3

pg150 Book 7

10/20/94

Ligations

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16 | 17 | 18 | 19 |
|------------------|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|
| HTPAN08 9/30 | 4 | 4 | 4 | | | | | | | | | | | | | | | | |
| 10/20 | | | | 6 | 6 | 6 | | | | | | | | | | | | | |
| HTPB411 9/30 | | | | | | | 3 | 3 | 3 | | | | | | | | | | |
| 10/20 | | | | | | | | | | 6 | 6 | 6 | | | | | | | |
| PD10 1.1 | | | | | | | | | | | | | 6 | 6 | 6 | | | | |
| PD1024 10/11 | 2 | | | 2 | | | 2 | | | 2 | | | 2 | | | 2 | | | |
| 10/14 | | 2 | | | 2 | | | 2 | | | 2 | | | 2 | | | 2 | | |
| 10/18 | | | 2 | | | 2 | | | 2 | | | 2 | | | 2 | | | 2 | |
| 10x Buffer | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 |
| H ₂ O | 11 | 11 | 11 | 9 | 9 | 9 | 12 | 12 | 12 | 9 | 9 | 9 | 9 | 9 | 9 | 15 | 15 | 15 | 17 |
| T4 Ligase | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 |
| | 20 | 20 | 20 | 20 | 20 | 20 | 20 | 20 | 20 | 20 | 20 | 20 | 20 | 20 | 20 | 20 | 20 | 20 | 20 |

Incubate @ 16°C O/N

10/24/94

Q1

Transform into M15 Chemically Competent Cells

Thaw M15 Cells on ice
 Aliquot 100ul into @ fresh Sterile tubes
 Add 10ul of ligation to tubes.
 Incubate on ice 1 hr
 Heat 42°C 45 sec
 Place on ice
 Add 400ul LB
 Incubate 37°C 1 hr
 Plate 175ul onto LBt Amp Kan plates
 Incubate RT over the weekend

HEWLETT-PACKARD PRODUCTS CO.
COMPUTATION NOTEBOOK # 535

Department Mol. Biol.
Subject 10/20/94 - 02/02/95
Name ANN KIM # 8
Address _____



43-648

Computation Notebook

Dennison Stationery Products Co., Framingham, MA 01701



75 Sheets
11 7/8" x 9 1/2"
4x4 Quad.

0 73333 43648 8

Ruben EXHIBIT 2091
Ruben v. Wiley et al.
Interference No. 105,077
RX 2091

4

HTPB411 + HTPAN08 in PD10

4/94

Plates look OK.

PD10 1.1 + 2.4 looks very good so ligations worked well

PD10 vector colonies had many colonies
 So did Vector + ~~PCR~~ fragment
 try PCR on colonies.

Pick 200 each into 200 μ l LB + Amp + Kan
 in 96 well dish - HTPAN08 + HTPB411

Incubate 37°C 4 hrs w/ aeration

PCR

HTPAN08

| | | 200 μ |
|------------------|-----|-----------|
| 10X PCR | 3.2 | 640 |
| 10X dNTP | 3.2 | 640 |
| Tag | 0.2 | 40 |
| 2500 | 0.2 | 40 |
| Tag | 0.2 | 40 |
| H ₂ O | 23 | 4600 |
| Cult | 2 | |

HTPB411

| | | 200 μ |
|------------------|-----|-----------|
| 10X PCR | 3.2 | 640 |
| 10X dNTP | 3.2 | 640 |
| 2501 | 0.2 | 40 |
| 2502 | 0.2 | 40 |
| Tag | 0.2 | 40 |
| H ₂ O | 23 | 4600 |
| Culture | 2 | |

PCR Program #69

95°C 5 min
 95°C 20 sec
 55°C 20 sec } 30x
 72°C 1 min
 72°C 7 1/2 min
 4°C Hold.

HTPB411 & HTPAN08 in PP10

5

10/24/94


Run 10 μ l on gel w/ 1 kb ladder

- Nothing worked!

Didn't bother to take all pictures.

HTPAN08

had strange size band.



Retry - Make New Primers today later.

- Phenol extract from PCR - then
ethanol ppt - Digest O/N.

- Run on LMP gel.

10/26/94

Try using S' Bam primers from
ligation into pAZ vector.

HTPAN08 - 7689

HTPB411 - 2429.

6

HTPAN03 + HTPB11 8' in PD10

10/26/94

DNA

100µl

| | HTPAN03 | | | HTPB11 | |
|------------------|---------|------|------------------|--------|-------|
| | ① | ② | | ③ | ④ |
| 10xPCR | 10.0 | 10.0 | 10xPCR | 10 | 10 |
| 10x dNTP | 10.0 | 10.0 | 10x dNTP | 10 | 10 |
| 2499 | 0.5 | 0.5 | 2501 | 0.5 | 0.5 |
| 2500 | 0.5 | — | 2502 | 0.5 | — |
| 368A | — | 0.5 | 2429 | — | 0.5 |
| Taq | 0.3 | 0.3 | Taq | 0.3 | 0.3 |
| H ₂ O | 77.7 | 77.7 | H ₂ O | 77.7 | 77.7 |
| DNA | 1 | 1 | DNA | 1 | 1 |
| | 100 | 100 | | 100µl | 100µl |

Run PCR Program #101. Modified

95°C 5min
 95°C 20sec
 55°C 20sec
 72°C 1min
 72°C 7 1/2 min.
 4°C hold

25X

Run 5µl on gel with 1Kb ladder



looks good.

Add equal Volume PEG/NaCl
 Sit on ice 10min
 Spin 10min
 Remove supernatant
 1x70% ethanol wash - 500µl

Spin 5min
 Remove supernatant
 let pellet dry slightly
 Resuspend in 200µl TE
 Add equal Volume Phos
 Vortex

HTPAD08 + HTPB08 in R10

7

10/20/91

Spin 5 min

Transfer upper Aqueous layer to fresh tube
repeat

Add 2 volumes Ethanol (400 μ l)
1/10 volume 3M NaAc (20 μ l)

Mix

Let sit on ice 1/2 hr

Spin 10 min

Remove supernatant

1x 70% Ethanol wash - 1000 μ l

Spin 5 min

Remove supernatant

Allow pellet to dry

Resuspend in 200 μ l TE

Set up digestions

| | |
|------------------|-------------------|
| DNA | 100 μ l |
| 10x #2 | 20 μ l |
| H ₂ O | 79 μ l |
| XbaI/Bam | 0.5/0.5 |
| | <hr/> 200 μ l |

Incubate 37°C / overnight

10/27/91

Run 10 μ l on 1% gel with 1 kb ladder



gel looks good

Add sample Dye. Angel

Run on 0.8% 1 MP

Agarose gel


Run 80V 1 1/2 hrs

Ethanol ppt samples (see above) - 1x 70% wash
Resuspend in 40 μ l TE - Now ready for gel

8

HTPA008 & HTPB411 in PD10

10/27/94




for Sample #1 cut out
1.8 Kb fragment too

Gene Clean
Add 1000 μ l NaI
Heat 55°C 5 min
Mix well
Add 1 μ l glass milk
Let incubate at
Room temp 5 min
w/ occasional mixing
Spin 7 Sec

3X { Remove Supernatant
Resuspend pellet in 400 μ l Wash Buffer
Spin 7 Sec
Remove Supernatant
Spin 7 Sec
Remove Supernatant
Resuspend pellet in 20 μ l TE
Heat 55°C 1 min
Spin 7 Sec
Transfer to fresh tube
Repeat 20 μ l TE

Rem 2 μ l on gel with 1 Kb ladder



Use to Set up digests

Do Not use the
1.3 Kb fragment from
#1.

HTPAN08 + HTPB411 in pD10

9

10/27/94

Ligation Reactions.

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | 15 | Ke |
|------------|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|
| HTPAN08 #1 | 4 | 4 | | | | | | | 4 | | | | | | | |
| HTPAN08 #2 | | | 4 | 4 | | | | | | 4 | | | | | | |
| HTPB411 #3 | | | | | 4 | 4 | | | | | 4 | | | | | |
| HTPB411 #4 | | | | | | | 4 | 4 | | | | 4 | | | | |
| PD10 1) | | | | | | | | | | | | | | | | 4 |
| PD10 10/19 | 1 | - | 1 | - | 1 | - | 1 | - | | | | | 1 | | | 1 |
| PD10 10/12 | - | 1 | - | 1 | - | 1 | - | 1 | | | | | | 1 | | |
| 10x Buffer | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | - | - | - | - | - |
| H2O | 12 | 12 | 12 | 12 | 12 | 12 | 12 | 12 | 13 | 13 | 13 | 13 | 16 | 16 | 17 | 12 |
| T4 Ligase | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 |
| | 20 | 20 | 20 | 20 | 20 | 20 | 20 | 20 | 20 | 20 | 20 | 20 | 20 | 20 | 20 | 20 |

Incubate at 16°C O/N.

10/28/94

Transform M/S Chemically Competent cells.

Thaw Cells on Ice
 Aliquot 100ul of cells into fresh tubes
 #16 may get 450ul.
 Add 10ul of ligation to tubes
 Incubate on ice 1 hr.
 Heat 42°C 245 sec
 Put on ice
 Add 400ul LB
 Incubate 37°C 1 hr
 Plate onto LB + Amp + Kan plates, 150mm.
 plate 1-16 100ul
 1-8 250ul also
 plate #6 (100ul) mixed w/ #7

10 HTPANOS & HTPBY11 in PD1D

10/28/94

Incubate 37°C O/N.
 Asked Paul Moore to place at 4°C
 Tomorrow.

10/31/94

Paul did NOT take plates out
 - there were colonies. So try to PCR

Pick 2400 each 1-8 & PCR

| | 1+4 HTPANOS | | | | 545-8 HTPB ^y | | | |
|------------------|----------------|------|----------|------|----------------------------|------|----------|------|
| | 1-2 | 50 | 3-4 | 50 | 5-6 | 180 | 7-8 | 160 |
| 10x PCR | 3.2 | 160 | 3.2 | 160 | 3.2 | 180 | 3.2 | 160 |
| 10x dNTP | 3.2 | 160 | 3.2 | 160 | 3.2 | 160 | 3.2 | 160 |
| 2499 | 0.2 | 10 | 2499 0.2 | 10 | 2501 0.2 | 10 | 2501 0.2 | 10 |
| 2500 | 0.2 | 10 | 7089 0.2 | 10 | 2502 0.2 | 10 | 2499 0.2 | 10 |
| Taq | 0.2 | 10 | 0.2 | 10 | 0.2 | 10 | 0.2 | 10 |
| H ₂ O | 23 | 1150 | 23 | 1150 | 23 | 1150 | 23 | 1150 |
| 8 | 2 | — | 2 | — | 2 | — | 2 | — |
| Cult | 2 | — | 2 | — | 2 | — | 2 | — |

PCR Prog 69.
 Nothing set up.

Re-make
 Try by aw.

PCR 2400

13
 100

HTPB111 & HTPAN08 m. PD10

13

(Pg 12)

10/31/94

Make new insert again
PCR w/ New and old primers.

| | (1) | (2) | (3) | (4) |
|------------------|------|------|------|------|
| HTPAN08 | | | | |
| 3' Xba 2499 | 0.2 | — | 0.2 | — |
| 3' Xba 2656 | — | 0.2 | — | 0.2 |
| 5' Bam 2500 | 0.2 | — | — | 0.2 |
| 5' Bam 7689 | — | 0.2 | 0.2 | — |
| 10x PCR | 10 | 10 | 10 | 10 |
| 10x dNTP | 10 | 10 | 10 | 10 |
| Taq | 0.2 | 0.2 | 0.2 | 0.2 |
| DNA (100ng/ul) | 1 | 1 | 1 | 1 |
| H ₂ O | 78.4 | 78.4 | 78.4 | 78.4 |
| | 100 | 100 | 100 | 100 |

| | (5) | (6) | (7) | (8) |
|------------------|------|------|------|------|
| HTPB111 | | | | |
| 3' Xba 2501 | 0.2 | — | 0.2 | — |
| 3' Xba 2653 | — | 0.2 | — | 0.2 |
| 5' Bam 2502 | 0.2 | — | — | 0.2 |
| 5' Bam 2429 | — | 0.2 | 0.2 | — |
| 10x PCR | 10 | 10 | 10 | 10 |
| 10x dNTP | 10 | 10 | 10 | 10 |
| Taq | 0.2 | 0.2 | 0.2 | 0.2 |
| DNA (100ng/ul) | 1 | 1 | 1 | 1 |
| H ₂ O | 78.4 | 78.4 | 78.4 | 78.4 |
| | 100 | 100 | 100 | 100 |

PCR Prog # 619

95°C 5min
 95°C 20 sec
 58°C 20 sec
 72°C 1min
 72°C 7 1/2 min
 4°C Hold

25x

14

HTPANDOS & HTPBY11

in PD10

1/1/94

Run 5ul on gel with 1Kb ladder.



Add equal Volume PEG/NaO
Spin 15 min
70% Ethanol Wash.
Resuspend 100ul.
Digest D/N.


| | |
|------------------|-------|
| DNA | 50ul |
| H ₂ O | 38 |
| 10x#2 | 10 |
| Xba/Bam | 1/1 |
| | 100ul |

incubate reaction
37°C Overnight

(3)

11/2/94

Prep Precipitate Digest
Add 1/10 Vol 3M NaAcetate
2 Vol 100% Ethanol
Spin 15 min
70% Ethanol Wash.
Allow pellet to dry
Resuspend in 45ul of TE
Add 5ul 10x Sample Buffer
Run on 0.8% LMP gel.
HTPANDOS HTPBY11




HTP4008 & HTP4011 8 in PMA

15

11/2/94

Con
Cut out fragment
Gene clean
Add 1000ul Pa.I
Heat 65°C 2 min
Mix well
Add 7ul Glass milk
Incubate at Room Temp 5 min
With occasional mixing
Spin 7 sec
Remove Supernatant
Resuspend pellet in 400ul Wash
Buffer
3X Spin 7 sec
Remove Supernatant
Spin 7 sec
Remove Supernatant
Resuspend pellet in 20ul TE
Heat 65°C 1 min
Spin 10 sec
Transfer to fresh tube
Repeat with 20ul TE
Run 2ul of eluted fragment on
1% TAE gel with 1 kb ladder



With fragments
from 10/27
and the New Fragment
Set up ligations

16

HTRANDS AT 10/27 & 11/2 into PD10

11/2/94

ligation Reactions

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | 15 |
|------------------------|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|
| HTRANDS 2499+2500 9/27 | 3 | | | | | | | | | | | | | | |
| 2499+7689 | | 3 | | | | | | | | | | | | | |
| 11/2 2499+2500 | | | 3 | | | | | | | | | | | | |
| 2656+7689 | | | | 3 | | | | | | | | | | | |
| 2499+7689 | | | | | 3 | | | | | | | | | | |
| 2656+2500 | | | | | | 3 | | | | | | | | | |
| HTRANDS 2501+2502 9/27 | | | | | | | 3 | | | | | | | | |
| 2501+2429 | | | | | | | | 3 | | | | | | | |
| 11/2 2502+2502 | | | | | | | | | 3 | | | | | | |
| 2653+2429 | | | | | | | | | | 3 | | | | | |
| 2501+2429 | | | | | | | | | | | 3 | | | | |
| 2653+2502 | | | | | | | | | | | | 3 | | | |
| PD10 Bam/Hha | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 |
| IDX Lig Buffer | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 |
| H ₂ O | 13 | 13 | 13 | 13 | 13 | 13 | 13 | 13 | 13 | 13 | 13 | 13 | 13 | 13 | 13 |
| T4 Ligase | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 |
| PD10 .1 kb | | | | | | | | | | | | | | | 3 |

Incubate 16°C O/N

11/3/94

Transform M15 cells

use 1-8 from 10/27 ligation
as well as 1-15 from 11/2

To 100 µl of cells add 10 µl of ligation

Let sit on ice 1 hr
Heat 42°C for 45 sec
Add 400 µl LB

Incubate 1 hr at 37°C

plate onto LB + Amp + Kan
plates

HTPAW08 + HTPB411 in pD10

17

11/5/94

150 mm plates

use 150 ul 1-8 & 1-12
use 300 ul 1-8 & 1-12

Incubate 37°C O/N

11/4/94

Pick 1-8 10/27 16/28/29
into 200 ul LB + Kan + Amp plates
Pick 1-12 11/2 into
200 ul LB + Amp + Kan

Incubate 37°C w/ aeration
20 PCR

use internal primers for confirm

HTPB411 RP10 + FP15
HTPAW08 RP12 + FP18

HTPB411

| | | 200x |
|------------------|------|------|
| 2400 | 0.25 | 50 |
| 2.7945 | 1 | 200 |
| 10x PCR | 3.2 | 1400 |
| 10x dNTP | 3.2 | 1400 |
| Taq | 0.2 | 60 |
| H ₂ O | 2.15 | 2430 |
| Culture | 2 | |

HTPAW08

| | | 200x |
|------------------|------|------|
| 16638 | 0.25 | 50 |
| 3407 | 0.25 | 80 |
| 10x PCR | 3.2 | 240 |
| 10x dNTP | 3.2 | 240 |
| Taq | 23.1 | 460 |
| H ₂ O | 0.2 | 40 |
| Culture | 2 | |

PCR Program H69

95°C 5 min
95°C 20 sec
55°C 20 sec
72°C 1 min
72°C 7 1/2 min
40°C hold

Store Exs at
-20°C over
weekend

Ⓢ Control Plasmid
Ⓢ Control H₂O

18

HTPAN08 + HTPBY11 in PD10

11/7/94

Rem 10 μ l of Reaction on gel with 16s
 loaded

Nothing set up again!

Talked w/ Steve -

Need to remake primers,
 Sites on old ones wrong

Submit both HTPAN08 + HTPBY11
 5' Bam a 3' ~~to~~ Xba

11/8/94

3' Xba HTPBY11 - Oligo Rec'd.
 (60°C. O/X)

11/10/94

Rec'd 5' Bam New HTPBY11
 5' Bam New HTPAN08

11/11/94

Rec'd 3' Xba New HTPAN08

11/12/94

Dried Oligos

Rem PCR Reactions w/ new Primers

pg 39

HTPAN08 / HTPB411 in PD10

39

pg 18

11/15/94

HTPAN08S04

HTPB411S15

| | |
|------------------|------|
| DNA (10ng/ul) | 1 |
| 3' Xba New | 0.6 |
| 5' Bam New | 6 |
| 10x dNTP | 10 |
| 10x PCR | 10 |
| H ₂ O | 72.1 |
| Taq | 0.3 |
| | 100 |

| | |
|------------------|------|
| DNA (10ng/ul) | 1 |
| 3' Xba New | 0.6 |
| 5' Bam New | 6.7 |
| 10x dNTP | 10 |
| 10x PCR | 10 |
| H ₂ O | 71.4 |
| Taq | 0.3 |
| | 100 |

PCR Program #69 modified

114 Run 5ul on gel w/ kb ladder

95°C 5min
95°C 20sec
55°C 20sec } 25X
72°C 1min
72°C 7 1/2 min
4°C Hold



Add equal Volume PEG-Mall

Spin 10min

Pour off Supernatant

Wash 1ml 70% ethanol

Spin 5min

Remove Supernatant

Resuspend in 100ul TE

Digestions:

| | |
|------------------|-----------|
| DNA | 2.5 |
| H ₂ O | 1.9 |
| 10X #2 | 5 |
| Xba / Bam | 0.5 / 0.5 |
| | 30ul |

↑ 2ul of Digest

Incubate RT over full Monday 11/14


40 HTPA008 / HTPB411 in PDI

11/14/94

Remain

Run digest on gel - see pg 39.
with 1 kb ladder

Rem Digest on 0.8% LMP Gel
Run 80V 2hrs



Cut out fragments
take picture

Use Qiaquick Gel Extraction
Kit

- Add 900ul QXI Buffer
- incubate 60°C 10 min
- Place in Spin Column
- Spin 60 sec - Remove Flowthrough
- Add remaining gel / QXI
- Spin 60 sec
- Remove Flowthrough

Add 750ul Buffer PE

Spin 60 sec

Remove flowthrough

Add 750ul Buffer PE

Spin 60 sec

Remove flowthrough

Spin 60 sec

Transfer Column to fresh tube

Add 50ul TE Heated 60°C

Spin 60 sec

Run Gel on gel with 1 kb ladder



looks good

HTPANC8 / HTPB411 in PD10

41

11/14/94

Set-up ligations

| | ① | ② | ③ | ④ | ⑤ |
|------------------|----|----|----|----|----|
| HTPANC8 | 4 | — | — | — | — |
| HTPB411 | — | 4 | — | — | — |
| PD10 1:1 | — | — | 4 | — | — |
| PD10 2:4 | 1 | 1 | 1 | 1 | 2 |
| 10X Buffer | 2 | 2 | 2 | 2 | 2 |
| 14 ligase | 1 | 1 | 1 | 1 | 1 |
| H ₂ O | 12 | 12 | 12 | 16 | 17 |

Set up 2 sets - 1 for 16°C overnight
A-1 set for RT 1 hr.

for A - 1 hr at RT

10 µl ligation

100 µl H15 Chem Competent cells

let sit on ice 1 hr

heat 42°C 45 sec

Place on ice

Add 400 µl LB

Heat 37°C 1 hr

Plate onto LBt Amp^r Kan plates

100 mm - 50 µl

or 150 mm - 200 µl

Incubate 37°C O/N

11/15/94

Plates look OK

very few colonies on Vector alone

No colonies on lig rxn alone

42

HTPAN08 / HTPB411 in PD10

11/15/94

Pick Colonies from 1, 2 & 3 into
LB + Amp + Kan
Incubate 37°C 4 hrs.

Run PCR Reactions

| HTPAN08504 | |
|------------------|-----------|
| 3' Xba New | 0.2 |
| 5' Bam New | 2 |
| 10x dNTP | 3.2 |
| 10x PCR | 3.2 |
| H ₂ O | 20 |
| Temp | 0.2 |
| Culture | 2 |
| | <u>32</u> |

| HTPB411S05 | |
|------------------|-----------|
| 3' Xba New | 0.2 |
| 5' Bam New | 2.5 |
| 10x dNTP | 3.2 |
| 10x PCR | 3.2 |
| H ₂ O | 20.7 |
| Temp | 0.2 |
| Cell time | 2 |
| | <u>32</u> |

① Control

PCR & insert -

② Control LB + Amp + Kan

with vector alone Colonies

take both sets of primers & make sure

No contamination by fragments

Run PCR Program 69.

| | | |
|------|-----------|-----|
| 95°C | 5 min | 30x |
| 95°C | 20 sec | |
| 55°C | 20 sec | |
| 72°C | 1 min | |
| 72°C | 7 1/2 min | |
| 4°C | Hold. | |

Run Dial on gel with 1 kb ladder
in 1% Agarose gel in TAE


(pg 47)

HTPAN08/HTPB411 in PD10.

47

(pg 42)

11/15/94



HTPAN08S04

HTPAN08S04
looks good.
HTPB411S15

Nothing so Re
plate ligations
from 16°C O/N



HTPAN08S04 -
Inoculate 5ml TB+
Amp + Kan with 15ul
of culture 2-11
Incubate 37°C O/N.

With ligation Rxns left at 16°C O/N

Thaw m15 chemically competent cells once
To 100ul of thawed cells add
10ul of ligation Rxn
45 (tube) add 10ul PD10 ~~reaction~~ planned
DNA as (+) control
Incubate on ice 1 hr.
Heat 42°C 45 sec
Place on ice
Add 400ul LB
Incubate 37°C 1 hr -
Plate onto LB + Amp + Kan plates
- 100 mm plates - 50ul
- 150 mm plates - 250ul
Incubate 37°C O/N

11/10/94

Boiling Mm pups

Spin Outlines 2min

Remove Supernatant

Resuspended pellet in 750ul STEH

RNase + Lipzyme

Boil 45 sec

Spin 10 min

Remove Pellet

Add equal Volume (750ul) of
13% PEG / 1.6M NaCl

Mix Well

Spin 10 min

Remove Supernatant

Add 1000ul 70% Etanol to wash
pellet

Spin 5 min

Remove ~~pellet~~ Supernatant

Allow pellet to dry slightly at RT.

Resuspended in 200ul TE

Run 2ul on gel with 1kb ladder
and pD10

plasmid looks
good.
looks like a lot of
chromosome.

Digest with Bam 1Xba to see if
the insert will "Pop" out
N300 bp

pg 53

HTPAW08 + HTPB411 is RPO10.


58

(pg 48)

11/18/94

| | |
|------------------|------------|
| DNA | 5ul |
| H ₂ O | 17.1 |
| IX #2 | 2.5 |
| Kan | 0.2 |
| Xba | 0.2 |
| | <hr/> 25ul |

Incubate 37°C 4hrs
Run on gel with 1 kb ladder
Digest pD10 - Should appear
1.1 kb fragment



looks like all are
correct
Submit for sequencing
with RPO10 5'
primer

FABPDD1-10 RPO1

Setup ligation for HTPB411

| | |
|------------|------------|
| Fragment | 10 |
| Vector | 1 |
| T4 ligase | 1 |
| 10X Buffer | 2 |
| | <hr/> 20ul |

Incubate RT 5 hrs

Thaw M15 cells on ice

To 100ul of M15 Chemically Competent cells

add 10ul of ligation reaction - ~~add~~ Use up all of ligation

Incubate on ice 1 hr

Heat 42°C 45 sec

Place on ice

Add 400ul LB

Incubate 37°C 1 hr

Plate 250ul into LB + Amp + Kan

Plates - 150 mm

HTPAN08 + ~~HTPCC9~~ / 16 + TPB4/11 in PD10


11/17/94

Some Colonies grew
 Pick into 200 μ l LB + Amp + Kan
 in 96 well dish.
 Incubate 37°C 4 hrs
 PCR -

| | |
|------------------|------|
| 3' Xba | 0.2 |
| 5' Bam | 2.5 |
| 10x PCR | 3.2 |
| 10x dNTP | 3.2 |
| Taq | 0.2 |
| H ₂ O | 20.7 |
| Culture | 2 |
| | 32. |

PCR Program #69
 use HTPB4/15/15
 3' Xba New / 5' Bam New
 Fragment as
 positive Control

Run 10 μ l on gel with 1 kb ladder.



looks like
 1 product !!!

Inoculate
 50 ml TP + Amp
 Kan
 with #20

Incubate 37°C
 O/N

11/18/94

Dragon Modified Media
 Spin Culture 4.5K 15 min
 Pour off supernatant
 Resuspended in Tris PI

HTPA008 & HTPB411 in PD10

55

11/18/98

Add 7ml P2

Mix gently

Add 7ml P3

Mix

Set on ice 20 min

Spin 4.5 K 15 min

Equilibrate TIP-100 with 5ml

QBT

Apply Supernatant to tube through

Rem wipe

allow to flow through

2X Wash Column 10ml @ QC Buffer

Elute DNA 10ml QF

Add 7ml isopropanol (0.7 volumes)

Mix well

Spin 8K 30 min

Remove Supernatant

Wash 10ml 70% Ethanol

Spin 8K 15 min

Pour off Supernatant carefully

Let pellet dry at RT

Resuspend in 30ul TE

Run gel on gel with 1 Kb ladder

Dilute 1:200 in H₂O

Read D₂₆₀/D₂₈₀



| Sample ID | abs 260.0 nm | abs 280.0 nm | abs 320.0 nm | 260.0 nm / 280.0 nm | 260.0 nm / 320.0 nm |
|-----------|--------------|--------------|--------------|---------------------|---------------------|
| 1 | 0.0557 | 0.0350 | 0.0031 | 1.6777 | 0.5708 |

as expected

HTPB411
PD10

~ 167 ng total

Digest with Bam / Xba to pop out 2.4 Kb fragment

DNA 4
H₂O 13.6
10X 2
Bam 0.2
Xba 0.2

Incubate 37°C 2 hrs

56

HTPAN08 + HTPB411

in PD10

11/18/94

Run 5ul on gel w/fin 1Kb ladder
 a 1ul of undigested.



11/28/94

Made New PA2 Primers 5' end
 for HTPAN08.

2782 5' Bam 51bp New 7.1pmul/ul
 2793 5' Bam 185bp New 8pmul/ul
 Use 3 A₂ from previous PA2 cloning
 7690

PCR.

| | |
|-------------------|-----------|
| 5' Bam | 3 |
| 3' A ₂ | 0.5 |
| POX PCR | 10 |
| 10x dNTP | 10 |
| H ₂ O | 75.2 |
| Taq | 0.3 |
| DNA (10ng/ul) | 1 |
| | <hr/> 100 |

[POX]

PCR Program # 69.

| | | |
|------|-----------|-------|
| 95°C | 5 min | } 25X |
| 95°C | 20 sec | |
| 55°C | 20 sec | |
| 72°C | 1 min | |
| 72°C | 7 1/2 min | |
| 11°C | HOLD | |

HTPANOS + HTPBY11

57

11/21/94

Rem 5ul on gel with 1Kb ladder



Add equal volume
PEG/NaCl-mix

Spin 10 min

Remove Supernatant

Wash pellet 1ml 70%
Etanol

Resuspend in 300ul TE

Rem 2ul on gel with
1Kb ladder

Digest fragment

Bam / ~~Hpa~~

| | |
|--------------------|------|
| DNA | 20 |
| H ₂ O | 23 |
| 10X#2 | 5 |
| Bam | 1 |
| Asp xba | 1 |
| | 50ul |

Digest 0/1N
store 20°C



11/28/94

Digested w/ xba by mistake.
Rel Digest with

Bam / ~~Hpa~~ Asp.

| | |
|------------------|------|
| DNA | 20 |
| H ₂ O | 23 |
| 10X#2 | 5 |
| Bam | 1 |
| Asp | 1 |
| | 50ul |

incubate 0/1N 37°C.

HTRANS + PAZ lig

11/30/94

Isolate on 0.8% LMP Gel &
 open clean
 (see pg 102)

Set up ligations

| | | | |
|------------------|----|----|----|
| 185 bp Bam/Asp | 10 | — | — |
| 51 bp Bam/Asp | — | 10 | — |
| PAZ Bam/Asp | 1 | 1 | 1 |
| 10x Buffer | 2 | 2 | 2 |
| T4 Ligase | 1 | 1 | 1 |
| H ₂ O | 16 | 16 | 16 |
| | 20 | 20 | 20 |

Incubate 16°C O/N.

12/1/94

Transform into M15 cells

Thaw M15 Chem. Competent Cells
 To 10 μ l of ligation add 100 μ l
 of cells
 Incubate on ice 1 hr
 Heat 42°C 45 sec
 Place on Ice
 Add 400 μ l LB
 Incubate 37°C 1 hr.
 Plate 100 μ l on 100 mm plates
 300 μ l on 150 mm plates

12/2/94

M15 cells Contaminated?

(pg 107)

62

PD10 B/HIE / HTPANOS frag / DNase P frag

11/30/94


Del isolate fragment on
0.8% agarose gel

Also Red Lam.

- 1 HTPANOS 185bp Bam/Aap
- 2 HTPANOS 51/185bp Bam/Aap
- 3 DNase P Bam/Xba
- 4 PD10 Bam/Hae

Wash 1 kb ladder

Run 80V 1 1/2 hrs
Cut out bands & take picture



Gene Clean fragments

Add 1ml NaCl

Heat 55°C 5 min

Add 5ul Glass milk

incubate RT 2 min

Spin 10 sec

Remove supernatant

Resuspend pellet in

50ul wash Buffer

Spin 10 sec

Remove supernatant

Spin 10 sec

Remove all supernatant possible

Resuspend pellet in 20ul TE

HTPANOS 51bp + 185bp DNase P

Resuspend pellet in 40ul TE PD10 Bam/Hae

incubate 55°C 2 min

Spin 10 sec

Transfer supernatant to fresh tube

Re-plate

Run 2ul on 1% gel with 1 kb ladder



store at -20°C

HTPANOS - see pg 58
DNase P - same



HPCC91, HE20142

65

pg60

12/5/94

- ① HE20142S15
- ② S14
- ③ 1 kb ladder
- ④ HPCC91S14
- ⑤ S15



Decided to use a 3' fragment for
screening instead
Inoculate 150 ml TB + Amp.
HE20142S14
HPCC91S14

12/6/94

Diagn. Max prep of HE20142S14
HPCC91S14

Spin cultures 4.5 K 15 min
Pour off supernatant
Resuspend pellet in 10 ml P1 + RNase
let sit at RT 5 min
Add 10 ml P2 mix gently
let sit at RT 5 min
Add 10 ml Cold P3 & incubate on
ice 20 min
Spin 4.5 K 15 min
Eqw. l. brgt. Colunox 10 ml QBT
Apply supernatant to column
Wash 2x 30 ml QC
Elute 15 ml QF
PPT w/ 10.5 ml isopropanol
Spin 9K 30 min
Pour off supernatant
Wash 10 ml 70% ethanol
Spin 9K 10 min
Pour off supernatant

Next Ph
12/15/94

HTPAWOS 504 51bp/185bp in PA2

67

pg 58

12/2/94

The lig Rn alone & M15 cells alone
produced colonies -
Don't use m15 - use DH5α.

Thaw Chemically Competent DH5α
cells on ice

12/3/94

10 100ul Cells add 10ul of
ligation

incubate on ice 30 min

Heat 42°C 45 sec

Put on ice

Add 400ul 1B

incubate 37°C 1 hr

Plate 200ul onto 150 mm LB+Amp

incubate 37°C O/N

12/3/94

Take plates out of 37°C

Put at 4°C

12/6/94

Pick into 200ul LB+Amp

incubate 37°C 2/N w/antibiotic

PCR To test for inserts

12/7/94

| 51bp | | 40X | 185bp | | 40X |
|-----------------------|-----|-----|-----------------------|-----|-----|
| 5' Bam SI | 0.2 | 8 | 5' Bam 185 | 0.2 | 4 |
| 3' Asp | 0.2 | 8 | 3' Asp | 0.2 | 4 |
| 10x dH ₂ O | 3.2 | 120 | 10x dH ₂ O | 3.2 | 64 |
| 10x PCR | 3.2 | 128 | 10x PCR | 3.2 | 64 |
| Taq | 0.2 | 85 | Taq | 0.2 | 4 |
| H ₂ O | 21 | 840 | H ₂ O | 21 | 420 |
| Cult | 2 | — | Cult | 2 | — |

68 HTPAN08504 51bp/185bp in PAZ

12/31/92


PCR prog # 66

| | | |
|------|-----------|-------|
| 95°C | 5min | } 30x |
| 95°C | 20sec | |
| 55°C | 20sec | |
| 72°C | 1min | |
| 72°C | 7 1/2 min | |
| 4°C | Hold | |

Run 10ul on 1% gel with 1kb ladder.

HTPAN08504 51bp

HTPAN08504 185bp



Pick @ HTPAN08504 185bp into 5ml
TB+Amp
Incubate 37°C O/N w/aeration

12/8/94

Boiling min prep

Spin 2ml culture 5min

~~Remove~~ Remove Supernatant

Resuspend pellet 750ul STE + RNase +
lysozyme.

(pg 84)

HTPHW08 PA2

89

pg 68

D/8/94

Boil 1 min
 Spin 10 min
 Remove pellet
 Add 750 μ l PEG/NaCl
 Mix well
 Spin 10 min
 Remove supernatant
 Add 1000 μ l 70% Ethanol
 Wash pellet
 Spin 5 min
 Remove supernatant
 Allow pellet to dry at RT 10 min.
 Resuspend pellet 200 μ l ϕ TE
 Run 2 μ l on gel with 1 kb ladder + λ HVI



Substrate for

Try digesting with Bam/Asp

| | | | |
|-----|--------------------|------|-------|
| | | | 19x |
| amp | DNA | 4 | 429.4 |
| | H ₂ O | 22.6 | 57 |
| | 10x ⁺ 2 | 3 | 3.6 |
| | Bam | 0.2 | 3.6 |
| | Asp | 0.2 | 3.6 |

Incubate 37C O/A

16. APR 94
 12/2/94

90

HTPAND08 PA2 / PD10

12/9/94

Does not look like anything digested
 Submit for Sequencing
 RP06
 RP19

1/11/95

Remake inserts w/ new primers
 for BlnA PA2 & PD10

PA2 has Bam ~~HT~~ HT site with
 Kozak seq.

51bp + 185 bp fragments

PCR ③

PD10

HTPAND08 51bp

| | |
|------------------|-----|
| 9113 | 20 |
| 2742 | 2 |
| 10xPCR | 50 |
| 10xDNTP | 50 |
| 1µl | 2 |
| H ₂ O | 375 |
| DNA | 1 |

④

PD10

HTPAND08 185bp

| | |
|------|-----|
| 9114 | 15 |
| 2742 | 2 |
| | 50 |
| | 50 |
| | 2 |
| | 382 |
| | 1 |

①

PA2

HTPAND08 51bp

| | |
|------|-----|
| 9111 | 18 |
| 2742 | 2 |
| | 50 |
| | 50 |
| | 2 |
| | 377 |
| | 1 |

②

PA2

HTPAND08 185bp

| | |
|------|-----|
| 9112 | 15 |
| 2742 | 2 |
| | 50 |
| | 50 |
| | 2 |
| | 382 |
| | 1 |

DNA-HTPAND0804 plasmid DNA - 10 µg/µl.
 PCR using Program # 66

③ Control
 H₂O only

④ 125

PQE60 Nco / Bam Vector

97

12/19/94

Digest PQE60 = 0.54 μ g.

| | |
|------------------|-------------------------|
| DNA | 4 10 μ l |
| 10x #4 | 5 |
| H ₂ O | 40 |
| Bam | 0.5 |
| Nco | 0.5 |
| | <hr/> 50 |

Incubate 37°C O/N

12/20/94

Run 2 μ l on gel with 1 kb ladder



looks good.
Add 5 μ l 10x SB
Run on 0.8% TBE gel
80V 1 1/2 hr
Cut out fragment.
Take picture



Gene Clean.

Add 900 μ l NaI
Heat 55°C 5 min
Mix well.
Add 7 μ l Glass beads
Mix well
Incubate at RT 5 min
with occ. mixing
Spin 10 sec
Remove Supernatant
Resuspend pellet 500 μ l
Wash Buffer
Spin 10 sec

W

HTPAN08/HTPB411 in PGE60

99

12/20/9

PCR

HTPAN08504 3' NcoI / 3' Bgl II
 HTPB411S15 5' BspHI / BamHI
 for ligation into PGE60

| HTPAN | | HTPB4 | |
|------------------|-----|------------------|-----|
| 5' NcoI | 20 | 5' BspHI (5bp) | 20 |
| 3' Bgl II | 4 | 3' BamHI | 4 |
| 10x dNTP | 30 | 10x dNTP | 30 |
| 10x PCR | 30 | 10x PCR | 30 |
| H ₂ O | 234 | H ₂ O | 234 |
| Taq | 1 | Taq | 1 |
| DNA (100ng/ul) | 1 | DNA (100ng/ul) | 1 |
| | 300 | | 300 |

Set up rxn w/o DNA (30ul vol total)

PCR: B10x 30ul rxns.

Prog 66.

Run 5ul on gel with 1kb ladder

5
5
39ul
0.5
0.5
50ul

95°C 5min
 95°C 20sec
 65°C 20sec } 30x
 72°C 1min
 72°C 7 1/2 min
 4°C hold




Add equal vol PEG/NaCl to ppt. - mix well
 Spin 10min
 Remove Supernatant
 Wash 1x 70% Ethanol (100ul)
 Spin 5min
 Remove Supernatant
 Allow pellet to dry at RT 10min
 Resuspend in 100ul TE

100

HTPAN08 HTPB411 in PQE60

12/20/94

Run 50 μ l on 0.8% LMP gel
80 V 1 1/2 hrs

Cut out fragments

Gene clean
Resuspended - Total 40 μ l TERun 2 μ l on gel with
1 kb ladder

looks good.

Set up digestions:

| | |
|------------------|------------|
| DNA | 30 μ l |
| H ₂ O | 13 |
| 10x#4 | 5 |
| Enz/En | 1/1 |
| | 30 μ l |

| | | | | | |
|-----|---------|---|--------|---|--------|
| For | HTPAN08 | - | Nco | / | Bgl II |
| | HTPB411 | | Bsp HI | / | Bam |

Incubate 37°C O/N


12/22/94

Add 5 μ l 10x SB-
Run on 0.8% LMP gel with
1 kb ladder
80 V 1 1/2 hrs

HTP8008 | HTPB411 in PQEC60

101


12/22/9



Gene Clean

Resuspend dried elute
in 40 μ l Total

Run 2 μ l on gel with
1 Kb ladder



Set up ligation

| | ① | ② | ③ | ④ | ⑤ | ⑥ |
|--------------------------|----|----|----|----|----|---|
| HEPABOSS045'Nco/S'Bg/II | 4 | — | 4 | — | — | — |
| HTPB411S15 5'Bg/H3/S'Nco | — | 4 | — | 4 | — | — |
| 10x | 2 | 2 | 2 | 2 | 2 | 2 |
| H ₂ O | 12 | 12 | 13 | 13 | 16 | 1 |
| T4 | 1 | 1 | 1 | 1 | 1 | 1 |
| PQE 60 Nco/Bam | — | 1 | — | # | — | — |
| PQE 60 Nco/Bgl/II | 1 | — | + | — | — | — |

Incubate 4°C

1/3/95

Transform M15 supH cells
Chemically competent (long)

Thaw cells on ice
To 100 μ l of cells add 10 μ l of ligation
To 1 tube add 10 μ l PQE60 only
To 1 tube add 10 μ l of H₂O only

1/3/95

Incubate on ice 1 hr
 Heat 42°C 1 min
 Let sit on ice 2 min
 Add 400ul LB
 Incubate 37°C 1 hr.
 Plate onto 2 Bt Amp + Kan plates.

1A - 200ul, 300ul
 1B - 100ul, 400ul
 2A - 200ul, 300ul
 2B - 100ul, 400ul
 3 - 300ul
 4 - 300ul
 5 - 300ul
 6 - 300ul ligation Rxn.
 7 - PDE60 10ng 100ul
 8 - M15 cells alone 300ul

Spread Evenly with Beads
 Incubate 37°C O/N.

1/4/95

All plates grew -
 Even M15 cells alone and ligation Rxn
 alone
 Cells are Contaminated!

Told Lang - Said he would test cells
 & make New.

Try Redoing ligation with
 HTPANOS 504 185 bp fragment.
 (Pg 107)

HTPAN08504 185bp Fragment

103

10/28/94

PCR HTPAN08504 185bp fragment
for PQE10 Nco I Bgl II

| | |
|------------------|-------------|
| DNA | 1 μ l |
| H ₂ O | 237 |
| 10x PCR | 30 |
| 10x dNTP | 30 |
| 1 μ g | 1 |
| 5' Nco I 185bp | 20 |
| 3' Bgl II | 4 |
| | 300 μ l |

10x PCR H₂O only
do @ Prep.

Split 10 tubes

PCR Program #101

| | |
|------|---------|
| 95°C | 5 min |
| 95°C | 20 sec |
| 55°C | 20 sec |
| 72°C | 1 min |
| 72°C | 7.5 min |
| 72°C | hold. |

30x

done 4°C

Run
PCR (10 μ l each) 1 Kb ladder

1/3/94



Add equal Volume 13% PEG
1.4M NaCl
incubate on ice 10 min
mixing well
spin 10 min
pour off supernatant.

104

HTRANO8504

185bp Fragment

1/3/95

Wash pellet in 700ul 70% Ethanol.

Spin 5 min.

Pour off supernatant

Air dry pellet to dry at RT

Resuspend pellet in 200ul TE

Now Set up Digestion


| | |
|------------------|------------|
| DNA | 35 |
| H ₂ O | 8 |
| 10X#4 | 5 |
| Nco | 1 |
| Bgl II | 1 |
| | <hr/> 50ul |

Incubate 37°C O/N.

1/4/95

Run EMP Gel - 0.8% in TAE.

80V 1 1/2 hrs

Cut out fragment from gel and place in tube
Take picture of gel.

2 Tubes

1 Tube use 10M NaI from Bio 101

Other Tube use 10M NaI made here

Add 1000ul NaI
Heat 55°C 5 min

pg 107

~~HTLAD31~~ 185 bp Fragment

105

HTLAD31 PER Protein

1/2/95

PER - Protein

Digest HTLAD31 = 0.25 µg/µl

- (1) Xho / Eco RI
- (2) Sac / Kpn

| | (1) |
|------------------|--------------|
| DNA | 6 |
| H ₂ O | 20.6 |
| 10X #2 | 3 |
| Xho | 0.2 |
| Eco RI | 0.2 |
| | <u>30 µl</u> |

| | (2) |
|------------------|--------------|
| DNA | 6 |
| H ₂ O | 20.6 |
| 10X #4 | 3 |
| Kpn | 0.2 |
| Sac | 0.2 |
| | <u>30 µl</u> |

Incubate 37°C O/N

Submit for sequencing -

1/4/95

Run 10 µl in 1% TAE gel
with 1 Kb ladder

1- Xho / Eco RI ~ 4.7 Kb

2- Kpn / Sac ~ 800 bp + ~ 600 bp

Insert looks ~ 1.6 Kb

make primers -

RPD6 & PPD5

1/9/95 Received RP06 & RP04.

Submit for sequencing
with other RP02, RP03
and Reverse & Forward

1/11/95 Received FP05.

Submit for seq RP06 & FP05

1/13/95

large 3' untranslated region

~ 1.5 Kb long

about 1.2 Kb is untranslated

HTPAV08 SOL 185bp Fragment

107

pg 104

+ pg 102

1/4/95

Mix well to make sure all of gel
is dissolved.

Add 8ul Glass milk + mix.

Incubate at RT 2min w/ occasional
mixing

Spin 10 sec

Pour off Supernatant

Resuspend pellet in 40ul Wash Buffer

3 Spin 10 sec

Pour off Supernatant

Spin 10 sec

Remove as much of Supernatant as
possible

Resuspend pellet 30ul TE

Heat 55°C 2min

Spin 10 sec

Transfer Supernatant to fresh Tube

Resuspend pellet in 20ul TE

Heat 55°C 2min

Spin 10 min sec

Transfer to Tube

Run 2ul on gel with 1 kb ladder

1- using Gene Clean Kit
6A NaI

2- using Home made 6A NaI

Both look good...
Ready for ligations

1/5/95

Set up ligation for

HTPAV08 & HTPB411 in PDE60

1/5/95

| | ① | ② | ③ | ④ | ⑤ | ⑥ | ⑦ | ⑧ | ⑨ |
|---------------------------------|---|---|---|----|----|----|----|----|----|
| HTPANOS04 5bp 5'Nco/3'BglII | 6 | — | — | 6 | — | — | — | — | — |
| HTPANOS04 18bp 5'Nco/3'BglII | — | 6 | — | — | 6 | — | — | — | — |
| HTPBYS15 5' 5' BspHI / 3' BamHI | — | — | 6 | — | — | 6 | — | — | — |
| 10X Buffer | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 |
| PGE60 Nco/BglII | 2 | 2 | — | — | — | — | — | — | — |
| PGE60 Nco/BamHI | — | — | 2 | — | — | — | — | 2 | — |
| H ₂ O | 9 | 9 | 9 | 11 | 11 | 11 | 15 | 15 | 17 |
| T4 DNA Ligase | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 |

Make Cocktail

| | | |
|------------------|----|----|
| | | 9x |
| 10X Buffer | 2 | 18 |
| H ₂ O | 9 | 81 |
| T4 DNA lig | 1 | 9 |
| | 12 | |

Store ligations at 4°C till
fused M15 cells are found

1/9/95

STEVE said to transform into XL-1
Blue plasmid is correct - ~~it~~ can
transform into M15 cells.

Thaw XL-1 Blue Chemically Competent
Cells on ice
To 100 µl of ligation add 100 µl
thawed cells
Incubate on ice 1 hr - as ② controls
Use 100 µl cells only
Heat Shock 42°C 45 sec

P115

Maxi Prep HTPAN08504/HMEAA88

111

11/4/95

Inoculate 200ml TB + Amp
from frozen stocks of
HTPAN08504
HMEAA88A.

Cultivate 37°C w/ aeration
overnight.

11/5/95

Diagen Maxi Prep

Spin Cultures 4.5K 15min

Pour off Supernatant

Resuspend pellet in 10ml P1 +
RNase

Cultivate RT 5min

Add 10ml P2

Cultivate RT 10min

Add 10ml P3 - ice cold

Cultivate on ice 20min

Spin 4.5K 15min

Equilibrate Tip -500 with 10ml
DBT

Apply Supernatant to tip through
Kim wipe

Allow to flow through

Wash Column 3x

30ml QC

Elute DNA - 1.5ml (IF

Add 10.5ml (all 0.7 times & volumes)

Mix well

Spin 8K 30min

Pour off Supernatant

Wash pellet 15ml Cold 70%
ethanol

112

Maxi Prep HIRAND8504, HMEAA88

1/5/95

Spin 8K 15min
 Pour off supernatant
 Allow off pellet to dry -
 Dry o/n at Room Temp -
 with parafilm w/ holes covering

1/6/95

Resuspend pellet in 400 μ l TE
 Run 1/2 on gel with 1 Kb ladder
 & 2 Hind III marker

Dilute 1:100 in H₂O
 Read OD₂₆₀/280



| Sample ID | abs | abs | bg abs | 260.0 nm | 280.0 nm | |
|--------------|----------|----------|----------|----------|----------|-----------------|
| | 260.0 nm | 280.0 nm | 320.0 nm | 280.0 nm | 260.0 nm | |
| 1 HMEAA88 | 0.0402 | 0.0319 | 0.0231 | 1.9374 | 0.5181 | 0.2 μ g/ml |
| 2 HIRAND8504 | 0.2051 | 0.1314 | 0.0315 | 1.7389 | 0.5754 | 1.03 μ g/ml |

looks like HMEAA88 is all chromosomal

1/9/95

Steve said to try transforming with
 Boring prep DNA to get a
 clone

~~HTPB411~~ ~~HTPB411~~ HTPAN08 HTPB411 P115

(Pg 108)

1/9/95

Set on ice 2 min
Add 400ul LB
Incubate 37°C 1 1/2 hrs
Plate 200ul onto LB + Amp
150mm plates

Transform ligations from 12/22

Incubate 37°C O/N
Plate M15 Cells on LB + Amp
LB + Kan
LB + Amp/Kan

1/10/95

Plates look good.
No contamination

M15 Cells

Grew on LB + Kan
on LB + Amp + Kan
Not on LB + Amp.

LB + Amp/Kan plates don't have
Amp or enough Amp so contamination
can be seen.

Can use M15 Cells - Not Contaminated

Pick - colonies onto LB + Amp in
96 well dish.

HTPB411

HTPAN08

12/22 2A - 36
2B - 44
1/5 3 - 25

12/22 2 48
1/5 1 48
1/5 2 48

pick PQUEO Vector + PA2 Vector as controls

116

HTPANOS/HTPB411 in PQEG60

1/10/95

Incubate 37°C w/aeration 5hrs
 PCR using internal primer
 as well as the 2 external primers

| HTPB411 + PQEG60 | | HTPANOS 51bp + PQEG60 | | HTPANOS 185bp + PQEG60 | |
|------------------|------|-----------------------|-----|------------------------|------|
| FP20 | 0.5 | FR14 | 0.2 | FR4 | 0.2 |
| 2887 | 0.2 | 2888 | 0.2 | 2867 | 0.2 |
| 2886 | 0.2 | 2885 | 0.2 | 2865 | 0.2 |
| 10x dNTP | 3.2 | 10x dNTP | 3.2 | 10x dNTP | 3.2 |
| 10x PCR | 3.2 | 10x PCR | 3.2 | 10x PCR | 3.2 |
| H ₂ O | 18.7 | H ₂ O | 17 | H ₂ O | 28.8 |
| Taq | 0.2 | Taq | 0.2 | Taq | 0.2 |
| Culture | 2 | Culture | 2 | Cult | 2.0 |
| | 32 | | 32 | | 32 |

PCR - (+) control Plasmid DNA
 (-) control LB + Amp alone

PCR Prog 166

95°C 5min

95°C 20sec

55°C 20sec

72°C 1min

72°C 7 1/2 min

4°C Hold

30x

Set up Rxns

Let Cultures

grow 37°C O/N

w/aeration

1/11/95

Run 10 µl of PCR Rxn on a 1.3% Agarose
 gel in TAE with 1kb ladder

pg 117

HTPB411 in PAZ / HTPANOS + HTPB411 in PCE60 119





pg 110

pg 116

1/11/95



① HTPB411 + PAZ A-E
HTPB411 + PCE60 H

② HTPB411 + PCE60 A-E
HTPB411 + PAZ F-H



③ HTPANOS 51bp + PCE60

④ HTPANOS 185bp + PCE60



120

HTPB411 PAZ / HTPB411 + TRAN in PDE60

1/10/95

Inoculate ~~28~~ TB + Amp 5ml ϕ
do Boiling preps.

① HTPB411 + PAZ

1, 5, 7, 8, 12, 16, 17, 18, 22, 23, 24, 29, 30, 31, 32, 41, 42, 43

② HTPB411 + PDE60

2, 6, 7, 19, 23, 24, 29, 38, 41, 45, 48, 57, 58, 60, 69, 73, 77, 79

③ HTPA085.6p + PDE60

1, 2, 3, 5, 6, 7, 8, 9, 10, 11, 14, 21, 23, 24, 27, 31, 33, 34

④ HTPA08185.6p + PDE60

5, 6, 7, 8, 11, 13, 15, 16, 17, 20, 21, 22, 23, 24, 25, 26, 28, 32

Incubate 37°C w/ aeration overnight

1/12/95

Boiling min. Preps

Spin 2ml of culture 2min

Remove Supernatant

Resuspend pellet in 750 μ l STET +
RNase + Lysozyme

Boil 1min

Spin 10 min

Remove Pellet

Add equal Volume - 750 μ l of 13% PEG 8000
1.0 M NaCl

Mix well

Spin 10 min

Remove Supernatant

Wash Pellet 1000 μ l 70% EtOH


Mix

HIPB411 PAZ / HIPB411 + HTRANOS POE60. 121

1/12/95

Spin 5 min.
Remove ethanol
Allow pellets to dry at RT for 30 min
Resuspend pellets in 100 μ l TE

Run 2 μ l of minipreps on 1% TAE
gel with 1 kb ladder + PAZ & POE60



mix up in labeling -
 (1) is really #3
 (2) is really #4
 (3) is really #1
 (4) is really #2

Set up digests

For POE60.

| EcoRI / HII | | |
|------------------|------|-----|
| | | 55X |
| DNA | 4 | |
| 10X | 2 | 110 |
| H ₂ O | 13.6 | 748 |
| EcoRI | 0.2 | 11 |
| HII | 0.2 | 11 |

For PAZ

| Bam / Xba | | |
|------------------|------|-----|
| | | 20X |
| DNA | 4 | |
| 10X | 2 | 40 |
| H ₂ O | 13.6 | 272 |
| Bam | 0.2 | 4 |
| Xba | 0.2 | 4 |

Incubate 37°C O/N
 Digest HTRANOS504 HIPB411S15
 PAZ & POE60

122


HTPB4H P22 / HTPAN08 + HTPB411 PDE60

1/13/95

1,2,4 - Eco / HII

3, Bam / Xba digested

Run 10ul on gel with 1 kb ladder



looks like

4- HTPAN08 185bp + PDE60

1-17 correct - Sequence 4.

2- HTPB411 in PDE60

1-4,6,8-16 looks good Seq 4.

14,3 Need to be digested w/correct enzymes

Submit 4x2 (4 of each)
#4 with PDE60 & P16
#2 P23
for sequencing

1/16/95

Cleanup DNA of 4-1, 4-2, 4-3, 4-4
2-1, 2-2, 2-3, 2-4

2x Phenol

2x Qiagen

Elute into PPT - Wash

Resuspend in 80 ul TE - use to transform 15 cells

pg 227

HTPAND8 PAZ/PD10

125


pg 10

95°C 5min
 95°C 20sec
 55°C 20sec } 30x
 72°C 1min
 72°C 75min
 74°C Hold

Run 5µl on gel with 1L6 ladder

11/12/95

3 - HTPAND8 51 bp for PAZ
 4 - HTPAND8 185 bp for PAZ
 1 - HTPAND8 51 bp for PD10
 2 - HTPAND8 185 bp for PD10



1/20

2

Precipitate Add equal Volume
 3% PEG/NaCl

Spin 10 min

Wash pellet 1000µl 70% Ethanol

Spin 5 min

Remove Supernatant

Set up digestion

DNA (PCR Fragment) 20

10x #2 Buffer 4

H₂O 15

Bam 0.5

Xba 0.5

40µl

Incubate 37°C overnight

1/13/95

Run All of digestion of 0.8% Agarose

Cut out gel fragments
Take picture



- 1 - HTPAN08 51bp } PD10
- 2 - HTPAN08 185bp } PD10
- 3 - HTPAN08 51bp } PAZ
- 4 - HTPAN08 185bp } PAZ

Gene Clean
Add 800 μ l NaI
Heat 55°C 5 min
Add 8 μ l Glass Milk
Mix well
Incubate at RT 5 min
w/ occasional mixing
Spin 10 sec
Remove Supernatant

- 3x { Add 500 μ l Wash Buffer
Resuspend pellet
Spin 10 sec
Remove Supernatant

- 2x { Spin 10 sec
Remove Supernatant
Resuspend pellet in 20 μ l TE
Heat 55°C 1 min
Spin 10 sec
Transfer Supernatant to fresh tube

Run on gel with 1 kb ladder and 1 μ l of PAZ and 1 μ l PD10 B/x



- 1 HTPAN08 51bp PAZ
- 2 HTPAN08 185bp PAZ
- 3 PAZ Run N/A
- 4 HTPAN08 51bp PD10
- 5 HTPAN08 185bp PD10
- 6 PD10 Run N/A

pg 122

HTPANO8/HTPB41 PAZ/PD10 127

Set-up ligations

1/13/95

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 |
|------------------|---|---|---|---|---|----|----|----|----|----|----|----|----|
| HTPANO8 51 PAZ | 6 | | | | | 6 | | | | | | | |
| HTPANO8 185 PAZ | | 6 | | | | | 6 | | | | | | |
| HTPANO8 210 PD10 | | | 6 | | | | | 6 | | | | | |
| HTPANO8 185 PD10 | | | | 6 | | | | | 6 | | | | |
| HTPB41 PD10 | | | | | 6 | | | | | 6 | | | |
| 10x Buffer | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 |
| H ₂ O | 9 | 9 | 9 | 9 | 9 | 11 | 11 | 11 | 11 | 11 | 15 | 15 | 17 |
| T4 Ligase 10ul | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 |
| PAZ Bam/xba | 2 | 2 | | | | | | | | | 2 | | |
| PD10 Bam/xba | | | 2 | 2 | 2 | | | | | | | 2 | |

Set up ligations on ice

Incubate 16°C over the weekend

1/14/95

Then M15 Chemically Competent cells

To 100ul of Cells add 10ul of
ligation (from 1/13/95) and 1ul
of DNA (from pg 122)

DI Tube of M15 cells only / 1 of 10ng PD10

Incubate on ice 1 hr

Heat 42°C 1:45 sec

Place on ice

Add 400ul LB

Incubate 37°C 1 1/2 hrs

Plate into LB + Amp 150mm plates

- ligations - 1-5 plate 100+300ul

- plasmid DNA - plate 50ul

- ligation controls 10-13 plate 300ul

- M15 cells plate 300ul

- M15 cells + PD10 plate 100ul

Incubate 37°C O/N @

1/16/95

H1PB411 + PA2

1/10/95 Transformalins

Pick 96 Colonies into LB + Amp

H1TPAN08 51bp + PDE10

1/10/95 Transformalins

48 from 1/5/95 lysis

48 from 2/22/94 lysis

into 200ul LB + Amp

Incubate 37°C w/ aeration o/n.

1/17/95

Transformations worked well -

No Colonies seen on
fragment only plates2-125 colonies seen on
Vector alone plateNo Colonies on - MIS cells alone
- leg. R. alone

Pick colonies into 200ul LB + Amp + Kan

Plate 1) ① H1TPAN08 51bp + PA2 - 48 A1 → D12

② H1TPAN08 185bp + PA2 - 48 E1 → H12

12) : ③ H1TPAN08 51bp + PD10 - 48 A1 → D12

④ H1TPAN08 185bp + PD10 - 48 E1 → H12

13) : ⑤ H1PB411 + PA2 PD10 - 16 A1 → B4

PA2 PD10 4

PA2 4

14) 2-1 ⑥ H1PB411 + PDE10 A1 → A12

2-2 ⑦ H1PB411 + PDE10 B1 → B3

2-3 ⑧ in MIS cells C1 → C12

2-4 ⑨ in MIS cells D1 → D12

HTPB411 + HTPANOS

129

1/17/95

E1-E12 4-1 (12)
 F1-F12 4-2 (12)
 G1-G12 4-3 (12)
 H1-H12 4-4 (12) } HTPANOS 185bp +
 PQEG60 M.15 cells.

Incubate at 37°C w/ airtight 4 hrs.

PCR - HTPB411 + PA2 + HTPANOS 51bp + PQEG60


| HTPANOS 51bp + PQEG60. | | |
|------------------------|------|-----------|
| | | 100X |
| 2888 | 2 | 200 |
| 2409 FPI6. | 0.15 | 15 |
| 10x dNTP | 3.2 | 320 |
| 10x PCR | 3.2 | 320 |
| H ₂ O | 2125 | 2125 |
| Temp | 0.2 | 20 |
| Cult. | 2 | |
| | 32ul | 30ul/Tube |

| HTPB411 + PA2 100X | | |
|--------------------|------|-----------|
| 2741 | 0.2 | 20 |
| 2796 RPO. | 1.2 | 120 |
| 10x dNTP | 3.2 | 320 |
| 10x PCR | 3.2 | 320 |
| H ₂ O | 22 | 2200 |
| Temp | 0.2 | 20 |
| Cult. | 2 | |
| | 32ul | 30ul/Tube |

PCR Prog. Col.
 95°C 5min
 95°C 20sec
 55°C 20sec
 72°C 1min
 72°C 7 1/2 min
 4°C Hold

Run 10ul on gel with 1 kb ladder.

Control use 1 B Kestrel only



HTPB411 + PA2
 2741 (3' x 10') + 2796 (RPO)



1-12
 1-12
 71012
 2410
 784
 -A12
 -B3
 -C12
 -D12

30


HTRB411 + HTPAN08

1/17/95

HTRAN08 51bp + PDE60 (2888 5' BspH1 + 7409 (FPII) S)



HTRB411 + PA2
2741 (3' Xba) + 7796 (RPA6)



HTRAN08 51bp + PDE60
2888 (5' BspH1) +
7409 (FPII)

Inoculate 5ml TB + Amp w/in.
① clones 1-9 from HTRB411 + PA2
HTRAN08 51bp + PDE60

Incubate 37°C w/aeration O/N

PCR Plates [1], [2], [3], [4]

HTPB411 & HTPAN08

131

Plate 11

11/17/95

① HTPAN08 51bp + PAZ

② HTPAN08 185bp + PAZ

| | | |
|------------------|------|-----------|
| | | 55x |
| 9111 | 1.7 | 66 |
| 2742 | 0.1 | 5.5 |
| 10x dNTP | 3.2 | 174 |
| 10x PCR | 3.2 | 174 |
| H ₂ O | 22.1 | 1215.5 |
| Tag | 0.2 | 11 |
| Cult. | 2 | |
| | 32ul | 30ul/tube |

| | | |
|------------------|------|-----------|
| | | 55x |
| 9112 | 1.7 | 66 |
| 2742 | 0.1 | 5.5 |
| 10x dNTP | 3.2 | 174 |
| 10x PCR | 3.2 | 174 |
| H ₂ O | 22.1 | 1215.5 |
| Tag | 0.2 | 11 |
| Cult. | 2 | |
| | 32ul | 30ul/tube |

control
HTPB411

③ Control PAZ Vector + LB Amp/Km

Plate 12

③ HTPAN08 51bp + PD10

④ HTPAN08 185bp + PD10

| | | |
|------------------|------|-----------|
| | | 55x |
| 9113 | 1.4 | 77 |
| 2742 | 0.1 | 5.5 |
| 10x dNTP | 3.2 | 174 |
| 10x PCR | 3.2 | 174 |
| H ₂ O | 21.9 | 1204.5 |
| Tag | 0.2 | 11 |
| Cult. | 2 | |
| | 32ul | 30ul/tube |

| | | |
|------------------|------|-----------|
| | | 55x |
| 9114 | 1.2 | 66 |
| 2742 | 0.1 | 5.5 |
| 10x dNTP | 3.2 | 174 |
| 10x PCR | 3.2 | 174 |
| H ₂ O | 22.1 | 1215.5 |
| Tag | 0.2 | 11 |
| Cult. | 2 | |
| | 32ul | 30ul/tube |

control
HTPB411

⑤ Control PD10 Vector + LB Amp/Km

Plate 13

⑤ HTPB411 + PD10

| | | |
|------------------|------|-----------|
| | | 24x |
| 2752 | 1.2 | 28.8 |
| 2741 | 0.2 | 4.8 |
| 10x dNTP | 3.2 | 76.8 |
| 10x PCR | 3.2 | 76.8 |
| H ₂ O | 22 | 528 |
| Tag | 0.2 | 4.8 |
| Cult. | 2 | |
| | 32ul | 30ul/tube |

⑥ Control PD10 Vector + LB Amp/Km

⑦ Control HTPB411

1/17/95

plate 4

① H1PBULL + PDE60

| | | SDX |
|------------------|------|----------|
| 28057 | 3 | 150 |
| 28056 | 0.2 | 10 |
| 10x dNTP | 3.2 | 160 |
| 10x PCR | 3.2 | 160 |
| H ₂ O | 20.2 | 1010 |
| Temp | 0.2 | 10 |
| Cult | 2 | |
| | 32 | 32 tubes |

② H1PANOB185 + PDE60

| | | SDX |
|------------------|-----|----------|
| 29027 | 0.2 | 10 |
| 28057 | 0.2 | 10 |
| 10x dNTP | 3.2 | 160 |
| 10x PCR | 3.2 | 160 |
| H ₂ O | 23 | 1150 |
| Temp | 0.2 | 10 |
| Cult | 2 | |
| | 32 | 32 tubes |

③ control H1PBULL

④ control H1PBULL H1PANOB

⑤ control LB Amp / Km

Per Prog 66-



Run 10ul on gel with 1k ladder

95°C 5min
 95°C 20sec
 55°C 20sec
 72°C 1min
 72°C 7 1/2 min
 4°C hold

BOX Run by plate

Plate 1

Plate 2



HTPB411 & HTPAN

183

plate [3] + plate [2]

plate [2]

11/13/95





plate [2]

plate [3]





plate [3] + plate [4]

plate [4]





HTPAN08 - HTPB411

11/17/95

Plate [4]

Controls



| Controls | | |
|---------------|---------|--------------|
| ① 9111 + 2742 | HTPAN08 | 51bp PAZ |
| ② 9112 + 2742 | HTPAN08 | 185bp PAZ |
| ③ 9113 + 2742 | ↓ | 51bp PD10 |
| ④ 9114 + 2742 | | 185bp PD10 |
| ⑤ 2752 + 2741 | HTPB411 | PD10 |
| ⑥ 2887 + 2886 | HTPB411 | PQ E60 |
| ⑦ 2967 + 2865 | HTPAN08 | 185bp PQ E60 |

Inoculate 5ml LB + Kan/Amp
with 2 colonies ea.

With 2 clones each of
HTPB411 + PQE60 - 2-1, 2-2, 2-3, 2-4 (A+B)
HTPAN08 185bp + PQE60 - 4-1, 4-2, 4-3, 4-4 (A+B)

Inoculate 5ml TB + Amp
with 9 clones of each.

| | | |
|--------------------|------------------|-----|
| HTPAN08 | 51bp + PAZ | - 9 |
| HTPAN08 | 185bp + PAZ | - 9 |
| HTPAN08 | 51bp + PD10 | - 3 |
| HTPAN08 | 185bp + PD10 | - 9 |
| HTPB411 | PQE60 | |
| HTPB411 | PD10 | - 9 |

Incubate 37°C O/D w/ Auralin

HIPB411 & HIPANOS8

135

1/18/95

HIPB411 + PQE60.

HIPANOS8 185bp + PQE60.

To 5ml LB + Amp/Kan add 300 μ l
of ON cultures.
Incubate 37°C w/ aeration till
OD₆₀₀ = 0.4 - 0.6. - 2 hrs
Add 100mM IPTG to 2mM + 110 μ l.
Incubate 37°C 4 hrs
Spin 750 μ l of culture - induced
and 1 of uninduced of HIPANOS8
HIPB411.

Remove supernatant
Resuspend pellet in 100 μ l H₂O
Add 100 μ l 2x dissociation Buffer.

2x
Dissociation Buffer: 0.25M Tris pH 6.9
4% SDS
20% Glycerol.
10% β -Mercapto Ethanol.
0.2% Bromophenol Blue.

Mix well
Heat 100°C 5 min
Spin 5 min

Put on ice
Run 10 μ l on 12.5% Acrylamide
stacking gel with Rainbow
marker

Run at 150V in 1x Running
Buffer for 1 1/2 hrs. - till dye
front is at bottom of gel.

Stain - 1/2 hr at 37°C
DeStain overnight at RT -

Boiling Miniprep.

A45)
A45)

1/18/95

HTPB411-LPAZ
 HTPAN08 51bp PAZ
 HTPAN08 185bp PAZ
 HTPAN08 51bp PD10
 HTPAN08 185bp PD10
 HTPAN08 51bp PDE } Boiling mixture

Some Cultures did Not grow.
 Open 2ml Culture
 Remove Supernatant.
 Resuspend pellet in 750ul STET
 + RNase/ Lysogyme.
 Boil. 1 min
 Spin 10 min
 Remove Pellet
 Add 750ul 13% PEG 8000 / 16HNaCl
 Mix Well
 Spin 10 min
 Remove Supernatant
 Wash pellet 1000ul 40% Ethanol
 Spin 5 min
 Remove Supernatant
 Allow Pellet to dry at RT 30 min.
 Resuspend pellet in 150ul TE
 Run 2ul on gel with 1 Kb ladder



A - HTPB411 + PAZ
 B - HTPAN08 51 + PAZ
 C - HTPAN08 185 + PAZ
 D - HTPAN08 51 + PDE
 E - HTPAN08 185 + PD10
 F - HTPAN08 51 + PD10

Setup digests

H7PB411 & H7PON08

137

for PAZ & PD10
Constructs
Digest: Bam/Xba

| | | |
|------------------|------|-----------|
| | | 40X |
| DNA | 5 | — |
| 10X#2 | 3 | 120 |
| H ₂ O | 21.6 | 874 |
| Bam | 0.2 | 8 |
| Xba | 0.2 | 8 |
| | 30ul | 25ul/tube |

For PGE Constructs
Digest: EcoRI/HindIII

| | | |
|------------------|------|-----------|
| | | 10X |
| DNA | 5 | — |
| 10X#2 | 3 | 30 |
| H ₂ O | 21.6 | 216 |
| EcoRI | 0.2 | 2 |
| HindIII | 0.2 | 2 |
| | 30ul | 25ul/tube |

Incubate all reactions at 37°C O/N.

Digestions of Minipreps:

Run 10ul of Digestion on 1% TAE
Agarose gel with 1kb ladder



A H7PB411 + PAZ
B H7PON08 51bp + PAZ
C 185bp + PAZ
D 51bp + PD10
E 185bp + PD10
F 51bp + PGE60

Looks like some digested correctly

clean-up and send in to be sequenced with internal primers

1/19/95

Dilute HTPB411 + ~~PAZ~~ PAZ 1, 2, 3, 4
 HTPAN08 51bp + PAZ 1, 2, 3, 4
 HTPAN08 185bp + PAZ 1, 3, 4, 5

1:200 - use 10 μ l to transform into DH5 α .

Thaw DH5 α Chemically Competent cells
 on ice
 To 100 μ l of thawed cells Add 10 μ l
 of Diluted DNA.
 incubate on ice 1 hr.
 Heat 42 $^{\circ}$ C - 45 sec
 place on ice
 Add 400 μ l LB
 incubate 37 $^{\circ}$ C 1 hr.
 plate 300 μ l onto 150 mm LB + Amp
 plates
 incubate 37 $^{\circ}$ C O/N.

PCR HTPB411 + PAZ
 HTPAN08 51bp + PAZ
 HTPAN08 185bp + PAZ
 with a PAZ specific T7 promoter
 primer. and the 3' end primer

| HTPAN | |
|------------------|----------------------|
| T7 | 0.2 |
| 2742 | 0.2 |
| 10X | 10 |
| 10X | 10 |
| H ₂ O | 78.2 |
| Taq | 0.4 |
| DNA | 1 μ l of Diluted |
| | 100 μ l |

| HTPB411 | |
|------------------|--------------|
| T7 | 0.2 |
| 2741 | 0.2 |
| 10X | 10 |
| 10X | 10 |
| H ₂ O | 78.2 |
| Taq | 0.4 |
| DNA | 1 of diluted |
| | 100 μ l |

4 Tubes of each

HTPAN08 / HTPBY11

139

Run Program #66.

1/19/95

| | | |
|------|-----------|-------|
| 95°C | 5min | } 30x |
| 95°C | 30sec | |
| 55°C | 20sec | |
| 72°C | 1min | |
| 72°C | 7 1/2 min | |
| 4°C | Hold | |

Run 5ul of PCR on 1% TAE gel with 1kb ladder



HTPB411 Did Not
PCR - Need to
get more
clones.

HTPAN08 look
good.

HTPAN08 - precipitate with equal Vol
13% PEG / 16M NaCl
mix well

Spin 10min

Pour off supernatant

Wash pellet 1ml 70% EtOH


Spin 5min

Pour off supernatant

Allow pellet to dry at RT 15min.

Resuspend in 100ul TE

Run 2ul on 1% gel with 1kb ladder



Don + Know about
upper band 2.5kb - but
band looks good
Do TNT 1/20.

1/19/95

Clean up Borling pups to submit
for sequencing

Add TE to GS all samples to 400ul.

2X Phenol / SEVAG (1:1)

Extract with Equal Volume

2X SEVAG Extract with Equal Volume

Add 1/10 vol (40ul) 3M NaOAcate pH 5.5

2 vol (800ul) 100% ethanol

Mix well

Sit on ice 10min

Spin 15min

Remove supernatant

Wash pellet 100ul 70% ethanol

Spin 5min

Remove supernatant

allow pellet to dry at RT 5min

Resuspend pellet in 100ul TE

Run Dueton gel with 1kb ladder



(A) 1-4 HTPAN08 516bp + PAZ

(B) 5-8 HTPAN08 1856bp + PAZ

(C) 9-12 HTPAN08 516bp + PD10

(D) 13-16 HTPAN08 1856bp + PD10

(E) 17-19 HTPAN08 516bp + PQE10

Submit for sequencing
FP16 & RPO6.

IRIS Names:

(A) FAS51BP1PAZRP06/FP

(B) FAS185B1PAZRP06/FP

(C) FAS51BP1PD10RPO6/FP

(D) FAS185B1PD10RPO6/FP

(E) FAS51BP1PQE10RPO6/FP

FP16

HTPANO8 (HTPB11)

141

1/19/95
1/24/95

circulate 5ml (BTBmgs) (Ken)
with HTPB411 + PQE60


HTPANO8 185 + PQE constructs

pg 135

Destained gel

HTPB11 + PQE60

HTPANO8 185 + PQE60




Rem more of sample
on gel - 10%

looks like #4 induced

stain - Destain

1/20/95



Can't really
tell
- try inducing more

PQ6 / RAB
dp / RAB
OG / RAB
rd / RAB
ld / RAB
ile

1/20/95

DO TNT with PCR &
T₇ fragments (pg 139)

| | |
|----------------------------|--------------|
| Rabbit Reticulocyte Lysate | 12.5 |
| T ₇ Buffer | 1 |
| T ₇ Polymerase | 0.5 |
| Amino Acid Omet | 0.5 |
| ³⁵ S methionine | 2 |
| RNA _{in} | 0.5 |
| PCR Product | 4 |
| H ₂ O | 4 |
| | <u>25 ul</u> |

Incubate 30°C 2 hours

Heat 90°C 5 min

Quick Spin

Add 20 ul of 2X Association Buffer
Add 5 ul of TNTRun 15 ul on 15% Acrylamide gel
with C₄ labeled Rainbow marker

Run gel 150V 1 1/2 hrs

a until dry front near to bottom

Cut off stacking gel

Fix gel in Fixative

10% H₂AC, 30% MeOH
for 20 min at 37°C

Run gel fix

Amplify in 30 ml Amplify

for 20 min at 37°C

Dry gel 1 1/2 hrs at 80°C

Place on film

leave at RT over weekend

1/22/95

Incubate 150 ul T₇ Amp / T₇ Amp Km
with cultures 30°C 0.5 hr



(pg 149)

IL-6 PQECs / PD10

pg 70

143

1/19/95



IL-6 PD10

looks like correct clones

| IL-6 PQE | | | | | IL-6 PD10 | | | | |
|----------|--------|--------|--------|---|-----------|-------|--------|-------|--|
| 1-A2 | 6-B0 | 11-D12 | 14-G5 | | 1-A2 | 6-C4 | 11-E5 | 16-G1 | |
| 2-A3 | 7-C1 | 12-E5 | 17-G8 | | 2-A7 | 7-C5 | 12-E11 | 17-H2 | |
| 3-A7 | 8-C3 | 13-E8 | 18-H5 | | 3-B1 | 8-C8 | 13-F3 | 18-H3 | |
| 4-A10 | 9-C9 | 14-E11 | 19-H8 | | 4-B4 | 9-C9 | 14-F9 | 19-H5 | |
| 5-A11 | 10-D10 | 15-F1 | 20-F10 | 7 | 5-B11 | 10-D7 | 15-G8 | 20-H7 | |

1/23/95

Inoculate 200ul LB + Amp/Kan
with Cultures 1-20 of each
Incubate o/n at 37°C w/ aeration to DO
Mini induction

1/24/95

To 200ul of fresh LB + Amp/Kan
add 50ul of fresh Cultures

HTPB4.11 + HTPAN 308

149

PA42

Develop 1/23/95

1/23/95

Qiagen Maxi Preps of:

HUEAA88

HTPAN08504 51bp + PAZ #1

+ PAZ #3

+ PAZ #4

185bp + PAZ #1

+ PAZ #3

+ PAZ #4

+ PAZ #5

+ PQE 60 4-1

+ PQE 60 4-2

HTPB411515 + PAZ 1


Make Glycerol stocks of all

800ul Bacterial Stock + 800ul 30% Glycerol

pg42

HTPB411 + HTPAN08

Develop TNT.



Qiagen Maxi Preps of:

HMEAA88

HTPAN08S04 51bp + PAZ #1

+ PAZ #3

+ PAZ #4

185bp + PAZ #1

+ PAZ #3

+ PAZ #4

+ PAZ #5

+ PQE60 4-1

+ PQE60 4-2

HTPB411S15 + PAZ 1


Make Glycerol stocks of all

800ul Bacterial Stock + 800ul 30% Glycerol

HTPB411 + HTPAN08

py42

1/23/95



HTPAN08504
51bp

HTPAN08504
188bp

TNT 4/20/95
PCR of PA2 clones
HTPAN08504 51bp
HTPAN08504 188bp
T2 + 3 Xba primers

15% Gel
1/31/95
ON - 80C

HTPB411S15 + PA2


+PQE60 4-2
+PQE60

0.1 0.1 0.1 0.1 0.1 0.1

PAJ42

HTPB411 + HTPANJ08

1/23/95




HTPB411S15 + PAZ + PÄE60 4-2

Make Glycerol stocks of all
800ul Bacterial Stock + 800ul 30% Glycerol

1/23/95

Spin Culture 4.5 K 20 min
 Pour off Supernatant
 Resuspend pellet in 10 ml P1 + RNase
 let sit RT 5 min
 Add 10 ml P2 - mix by swirling
~~the~~ looked for cell lyses
 Add 10 ml P3 -
 Incubate on ice 30 min
 Spin 8 K 30 min
 Prepare Column
 apply 10 ml QBT
 allow to flow through
 Apply Supernatant to Column through
~~the~~ Kim wipe
 Allow to flow through
 Wash Column 30 ml Wash Buffer
 Elute DNA 15 ml TE
 Add isopropanol 0.7 vol (10.5 ml)
 Mix well
 Spin 8 K 30 min
 Pour off Supernatant
 Wash pellet 10 ml 70% Ethanol (-20°C)
 Spin 8 K 15 min
 Pour off Supernatant
 Allow pellet to dry at RT for 10 min
 Resuspend pellet in 400 µl TE pH 7.6
 Run 1 µl on gel with 1 kb ladder
 PA2 & PA2eo

looks good



HTPAN08 / HTPB44

151

1/23/9

Read OD₂₆₀/280
Dilute 1:200 in H₂O

| Sample ID | abs 260.0 nm | abs 280.0 nm | bkg abs 320.0 nm | 260.0 nm 280.0 nm | 280.0 nm 260.0 nm | | 40x2 |
|------------------------|--------------|--------------|------------------|-------------------|-------------------|------------|--------|
| 1 HMEABE | 0.0466 | 0.0293 | 0.0014 | 1.6230 | 0.6162 | 0.4 ug/ml | 188 ug |
| 2 HTPAN08 PAZ 10.1271 | 0.0721 | 0.0021 | -0.0021 | 1.7411 | 0.5743 | 1.2 ug/ml | 520 ug |
| 3 HTPAN08 PAZ 30.1188 | 0.0663 | -0.0002 | -0.0002 | 1.7859 | 0.5599 | 1.2 ug/ml | 480 ug |
| 4 HTPAN08 PAZ 40.0983 | 0.0547 | -0.0023 | -0.0023 | 1.7855 | 0.5664 | 1.0 ug/ml | 400 ug |
| 5 HTPAN08 PAZ 80.0919 | 0.0533 | -0.0003 | -0.0003 | 1.7201 | 0.5814 | 0.4 ug/ml | 368 ug |
| 6 " 30.0033 | 0.0001 | -0.0053 | -0.0053 | 1.5858 | 0.6306 | X | |
| 7 " 40.1070 | 0.0630 | -0.0008 | -0.0008 | 1.6905 | 0.5915 | 1.1 ug/ml | 440 ug |
| 8 " 50.1304 | 0.0771 | 0.0035 | 0.0035 | 1.7244 | 0.5799 | 1.3 ug/ml | 520 ug |
| 9 HTPAN08 PAZ 40.0967 | 0.0605 | 0.0045 | 0.0045 | 1.6471 | 0.6071 | 0.9 ug/ml | 358 ug |
| 10 HTPAN08 PAZ 10.0481 | 0.0287 | 0.0002 | 0.0002 | 1.6798 | 0.5953 | 0.48 ug/ml | 192 ug |
| 11 HTPAN08 PAZ 10.0502 | 0.0281 | -0.0017 | -0.0017 | 1.7434 | 0.5736 | 0.50 ug/ml | 200 ug |

Dilute DNA to 200-250 ng/ml
Submit for sequencing
for HTPAN08 - FR16 & R206
for HMEABE - R201, R202, R204, R205 & R206
for HTPB44

Inoculate 300 ml LB + Amp + Kan
with 5 ml HTPAN08 504 185 bp + P6560
O/N culture (41-1)

Inoculate at 37°C w/aeration
until OD₆₀₀ ~ 0.4-0.6

Add IPTG to 2 mM (10 ml of 100 mM IPTG)
Inoculate at 37°C w/aeration
5 hours

Spin culture 8K 20 min
Resuspend in 30 ml 6M GdnHCl pH 8
Store 4°C till tomorrow

BBB/Chloride

pg 2
Biodata
1/20

HTPANOS 7 HTPBY 11

1/24/95

Inoculate 200 μ l LBT Amp Kanw/ \oplus clones of $\frac{7}{11}$

HTPANOS 51bp + PD10 - 3

HTPANOS 185bp + PD10 - 17

HTPANOS 51bp + PQE60 - 17

Incubate w/ aeration 1 hour - 37°C

Add 100mM IPTG to 2mM \rightarrow 25ml

Incubate 37°C w/ aeration 5 hrs

HTPANOS 51bp + PQE60 - None (w/ 10)

Spin Cultures

Add 20ml H₂O to Resuspend pellet

Add 20ml 2X Dissociation Buffer

Heat 100°C 5 min


Run on 10% Acrylamide gel w/

Low mol. Weight Marker

150 V 1 1/2 hrs

Stain 50 min

DE stain 30 min



Not enough gels to run remaining
 samples, but does not look
 like anything induced

looks like from sequences,
 there is a problem with
 3'bp end -
 remove primer

Ag 5
 9
 book
 #27

**This Page is Inserted by IFW Indexing and Scanning
Operations and is not part of the Official Record**

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

- ☐ BLACK BORDERS
- ☐ IMAGE CUT OFF AT TOP, BOTTOM OR SIDES
- ☐ FADED TEXT OR DRAWING
- ☐ BLURRED OR ILLEGIBLE TEXT OR DRAWING
- ☐ SKEWED/SLANTED IMAGES
- ☐ COLOR OR BLACK AND WHITE PHOTOGRAPHS
- ☐ GRAY SCALE DOCUMENTS
- ☐ LINES OR MARKS ON ORIGINAL DOCUMENT
- ☒ REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY
- ☐ OTHER: _____

IMAGES ARE BEST AVAILABLE COPY.

As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.