





Steven M. Ruben
Appl. No. 10/662,429

Department _____
Subject _____
Name <u>AND SIM #6</u>
Address <u>5/5/94 - 5/5/94</u>
 43-648
Computation Notebook
Dennison Stationery Products Co., Framingham, MA 01701

75 Sheets x 9 1/4" Quad.
0 73333 43648 8

Ruben EXHIBIT #89

BEST AVAILABLE COPY

Department _____
Subject _____
Name ANN NIM #6
Address 515/94 - 5000
 43-648
Computation Notebook
Denslson Stationery Products Co., Framingham, MA 01701

0 73333 43648 8
7 Sheets
x 9"
Quad.

Ruben EXHIBIT 2089
Ruben v. Wiley et al.
Interference No. 105,077
RX 2089

pg 136 Bld # 127
AMK # 5

10/13/94

inoculate 100 ml TB + Amp
Single colony of HSC120
and HTPANOS
Incubate 37°C w/ aeration
overnight

10/14/94

TIP - 500 - Maxi Prep.
Spin cells 5K 15 min
Pour off Supernatant
Resuspend pellet in 10 ml PI + RNase
Add 10 ml P2
Let sit at RT 5 min
Add cold P3 - 10 ml
Incubate on ice 20 min
Spin 9K 30 min 4°C.
~~Transfer supernatant to~~

Equilibrate Tip - 500 with 10 ml QBT
buffer.
Transfer Supernatant to Column.
Strain through Kim wipe
Wash tip 2x 30 ml QF Buffer
Elute with 15 ml QF buffer
Add 10.5 ml isopropanol & mix
Sit on ice 10 min
Spin 9K 30 min
Pour off Supernatant
Wash pellet 15 ml 70% Etanol
Let pellet Air dry.
Dissolve pellet in 400 µl TE

54

54 Maigi Pip HTRAND & HSUSHZO

6/14/94

Read OD 260/280

Dilute 1:200

Sample ID	abs	abs	bkg abs	260.0 nm	280.0 nm	
	260.0 nm	280.0 nm	320.0 nm	280.0 nm	260.0 nm	
1 HSEB H2O	0.1716	0.1051	0.0058	1.6688	0.5992	1.7 ug/ml
2 HTANOB	0.1617	0.0983	0.0060	1.6862	0.5931	1.6 ug/ml

Run test on gel with pBSK.



from Glycerol Stocks -80°C

9/16/20

Pg 136

~~Amount to be paid is \$100.00~~
~~for the purchase of 100 shares of~~
~~the common stock of the~~
~~company at \$1.00 per share~~
~~plus \$5.00 per share for~~
~~brokerage fees and commissions~~
~~totaling \$105.00~~

7/20/94

incubate 3 hrs at 37°C

Heat 70°C 20 min.

Store 4°C

Use 5 μ l to transform 20 μ l SAR -Plate 5 μ l

7/21/94

Plates look very dense

Do PCR on Rescue.

H50
H505H20P04

0.15

0.15

3.2

3.2

23.1

0.2

2

32

M13R

10x dNTP

10x PCR

H₂O

Tag

Rescue

HE2

HE20142P03

0.15

0.15

3.2

3.2

23.1

0.2

2

32

PCR Program #69

95°C - 5 min

95°C 20 sec

55°C 20 sec

72°C 1 min

72°C 7.5 min

4°C

30X

⊕ Controls
plasmid DNARun 10 μ l on gel with ⊕ Controls
and 1 Kb ladderA. P. L.
7/21/94

7/26/94

PCR - Protein Exp

25 ng/lx HE9MF73S07 & S05 - XFEZ
 25 ng/lx HTPAN08S04 - TNF α - Fac.

(pg 54)

①

HE9MF73S07	1
2245 3' Bgl II Stop	1
2242 5' Bam HI Start	1
10x dNTP	5
10x PCR	5
H ₂ O	36.8
Tag	0.2
	50

②

HE9MF73S05	1
2245 3' Bgl II Stop	1
2242 5' Bam HI Start	1
10x dNTP	5
10x PCR	5
H ₂ O	36.8
Tag	0.2
	50

③

HTPAN08S04	1
2244 5' Nco Start 251	1
2241 3' Bam HI Stop	1
10x dNTP	5
10x PCR	5
H ₂ O	36.8
Tag	0.2
	50

④

HTPAN08S04	1
2243 5' Nco Start 185	1
2241 3' Bam HI Stop	1
10x dNTP	5
10x PCR	5
H ₂ O	36.8
Tag	0.2
	50

⑤

HTPAN08S04	1
2239 5' Bam HI Start 188	1
2238 3' Hind III Stop	1
10x dNTP	5
10x PCR	5
H ₂ O	36.8
Tag	0.2
	50

⑥

HTPAN08S04	1
2240 5' Bam HI Start 251	1
2238 3' Hind III Stop	1
10x dNTP	5
10x PCR	5
H ₂ O	36.8
Tag	0.2
	50

11 tubes ⑤ DNA 1 tube ⑥ control.

7/26/94

- ⊖ Control - No DNA.
 ⊕ Control DNA ⊕ Amplified M13F4R

PCR Program # 0158

95°C 5 min

95°C 30 sec


55°C 30 sec

72°C 1 min

72°C 7.5 min

30X

Run 5ul on gel.




Add Equal Volume 13% PEG / 1.6M NaCl
 -20°C O/N

7/27/94

Spin 15 min
 1X 70% wash
 Dry pellet
 Resuspended in 50ul H₂O
 Run 1ul on gel

7/27/94



2-6 Lark 1000

Digest w/ appropriate
Enzymes

7/28/94

Digest w/ appropriate Enzymes

#(2) A

DNA	10
10X #3	30
BamHI	1
H ₂ O	16
	<hr/> 30

#(2) B

DNA	10
10X #3	30
EcoRI	1
H ₂ O	16
	<hr/> 30

(3) A

DNA	10
10X 4	3
Nco	1
H ₂ O	16
	<hr/> 30

(3) B

DNA	10
10X 4	3
Bam	1
H ₂ O	16
	<hr/> 30

(4) A

DNA	10
10X #4	3
Nco	1
H ₂ O	16
	<hr/> 30

(4) B

DNA	10
10X #4	3
Bam	1
H ₂ O	16
	<hr/> 30

7/28/94

(5) A
 DNA 10
 10X#2 3
 Bam 1
 H₂O 16
 30

(5) B
 DNA 10
 10X#2 3
 H_{III} 1
 H₂O 16
 30

(6) A
 DNA 10
 10X#2 3
 Bam 1
 H₂O 16
 30

(6) B
 DNA 10
 10X#2 3
 H_{III} 1
 H₂O 16
 30

Incubate 37°C 3 hrs.
 Add Alternate Enzyme. Incubate 3 hrs.
 Repeat

Ligations:

#2 - Bam HI / Bgl II - pQE 60 & 70. Bam / Bgl II
 #3 - Nco / Bam - pQE60 > Nco / Bam
 #4 - Nco / Bam - pQE60 > Nco / Bam
 #5 - Bam / H_{III} - pDIO > Bam / H_{III}
 H_{II} Bam / H_{III} - pDIO > Bam / H_{III}



7/29/94

Run on 0.8% LMP gel with 1 Kb ladder.

Cut out of gel -

Gene Clean fragments

7/29/94

Add 500 μ l NaI

Heat 55°C 5min

mix well

Add 5 μ l of Glass milk

Let sit RT 2min w/ occasional mixing

Spin 5sec

Remove Supernatant

Resuspend in 500 μ l Wash Buffer

3x Spin 5sec

Remove Supernatant

Spin 5sec

Remove all supernatant as possible

Add 30 μ l of TE

Heat 55°C 5min

Spin 5sec


Transfer Supernatant to fresh tube

Store -20°C to do ligations

page 145

Protein Express.

7/29/94



Add 500 μ l NaI
Heat 55°C 5 min
mix well
Add 5 μ l of Glass milk
Let sit RT 2 min w/ occasional mixing
Spin 5 sec
Remove Supernatant
3x [Resuspend in 500 μ l Wash Buffer
Spin 5 sec
Remove Supernatant
Spin 5 sec
Remove all supernatant as possible
Add 30 μ l of TE
Heat 55°C 5 min
Spin 5 sec
Transfer Supernatant to fresh tube
Store -20°C to do ligations

pg 145

PROTEIN Expression

145

pg 145

7/29/94

Set up ligations

DNA	2
Vector	1
10x Buffer	2
H ₂ O	14
1x DNA ligase	1
	20µl

Store 4°C
Over Weekend

DNA fragment			Vector	
1	HE94173505	Bam/Bgl	PQE60	Bgl II/Bam
2	3' Bgl Stop		↓	Bam/Bgl II
3	5' Bam start		PQE70	Bgl II/Bam
4			↓	Bam/Bgl II
5		Bgl Bam	PQE60	Bgl II/Bam
6			↓	Bam/Bgl II
7			PQE70	Bgl II/Bam
8			↓	Bam/Bgl II
9	HTPANC504 5' Nco 251	Nco/Bam	PQE60	Nco/Bam
10	3' Bam Stop	Bam/Nco	↓	Nco/Bam
11	HTPANC504 5' Nco 185	Nco/Bam		Nco/Bam
12	3' Bam Stop	Bam/Nco	↓	Nco/Bam
13			PQE60	Bgl II/Bam
14			↓	Bam/Bgl II
15			PQE70	Bgl II/Bam
16			↓	Bam/Bgl II
17			PQE60	Bam/Nco/Bam
18				
19	1µl of DNA fragment.			

8/1/94

Transform M15 cells

100ul of Ligations
 1ul of 10 ng/ul pBSK
 1ul of 10 ng/ul pD10

Use M15 Chemically Competent Cells

DNA or ligations + 100ul Cells

Incubate on ice 1 hr.

Heat 42°C 45 sec

Put on ice

Add 400ul LB

Incubate 37°C 1 hr.

Plate 100ul into LB + Amp + Kan plates

Incubate 37°C O/N.

8/2/94

pick 12 colonies of 1-8 into 96 well dish LB + Amp + Kan

pick 24 colonies of 9-12 into 96 well dish LB + Amp + Kan

(100x)		(50x)		(50x)	
PCR #2A+B (1-8)		#3 (9+10)		#4 (11+12)	
#2245	0.4	#2244	0.4	#2243	0.4
#2242	0.4	#2241	0.4	#2241	0.4
10x dNTP	3.0		3.0		3.0
10x PCR	3.0		3.0		3.0
Taq	0.2		0.2		0.2
H ₂ O	23		23		23
	30ul		30		30

Protein Expression

147

8/2/94


PCR Program #69.

95°C 5min
95°C 20sec } 30x
55°C 20sec }
72°C 1min }
72°C 7 1/2 min.


Run 12 ul on 1% TAE gel with 1 Kb ladder

#2 - NO (+) - Whole gel Blank -

Re set up ligation again - Make more
PCR product (omit fragment)



#3 -
HTPANSOY
5' Nco 251
3' Bam Stop



#4
HTPANSOY
5' Nco 185
3' Bam

Protein Expression

8/2/94

PCR Program #69.

95°C	5min
95°C	20sec
55°C	20sec
72°C	1min
72°C	7 1/2 min

 } 30x

Run 12 ul on 1% TAE gel with 1 kb ladder

#2 - No ⊕ - Whole gel Blank -
Re set up ligation again - Make more
PCR product (omit fragment)



#3 -
HTPANSO4
5' Nco 251
3' Bam Stop



#4
HTPANCBSO4
5' Nco 185
3' Bam

8/2/94

Inoculate 3ml LB + Amp^r Kan
 into
 incubate 37°C O/N 1st 40 of each.

8/3/94

Add 4ml LB + Amp^r Kan
 Add 100 mM IPTG to 2 mM (140ul)
 incubate 37°C 5 hrs
 Spin 5 min - 1 ml
 Remove supernatant
 Resuspend 100ul H₂O
 Add 100ul 2X Sample Buffer



Heat 100°C for 5 min
 Chill on ice
 Spin 5 min

Run 10ul on 10% Acrylamide gel
 with Rainbow Marker.
 Forget to have a control of commercial
 so will re-run 1 gel with controls
 and selected samples.



Run 150 V 1 hr
 Stain 30 min
 DE stain O/N at RT

8/4/94

HTPANOSY 5' Nco 251 Start 35' Bam Stop.



HTPANOSY 5' Nco 181 Start 3' Bam Stop



**This Page is Inserted by IFW Indexing and Scanning
Operations and is not part of the Official Record**

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

- ☐ **BLACK BORDERS**
- ☐ **IMAGE CUT OFF AT TOP, BOTTOM OR SIDES**
- ☐ **FADED TEXT OR DRAWING**
- ☐ **BLURRED OR ILLEGIBLE TEXT OR DRAWING**
- ☐ **SKEWED/SLANTED IMAGES**
- ☐ **COLOR OR BLACK AND WHITE PHOTOGRAPHS**
- ☐ **GRAY SCALE DOCUMENTS**
- ☐ **LINES OR MARKS ON ORIGINAL DOCUMENT**
- ☒ **REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY**
- ☐ **OTHER:** _____

IMAGES ARE BEST AVAILABLE COPY.

As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.