





Steven M. Ruben
Appl. No. 10/662,429

Department	_____
Subject	_____
Name	<u>ANN KIM #5</u>
Address	<u>2/4/94 - 5/19/94</u>
 43-648	
Computation Notebook	
<small>Dennison Stationery Products Co., Framingham, MA 01701</small>	
	75 Sheets 11 1/2" x 9 1/2" 4x4 Quad.
0 73333-43648 8	

BEST AVAILABLE COPY

Ruben EXHIBIT #88

Department _____
Subject _____
Name ANN KIM #5
Address 2/4/94 - 5/19/94
 43-648
Computation Notebook
Dennison Stationery Products Co., Framingham, MA 01701

75 Sheets
11" x 9"
4x4 Quad.
0 73333 43648 B

Ruben EXHIBIT 2088
Ruben v. Wiley et al.
Interference No. 105,077
RX 2088

HIPANOB Screening 10

7

Thursday 2/3/94

~~WEEK~~ HIPANOB Screening

HIPANOB

HTP Lib from Scott Weissner

1×10^9 pfu/ml

12 plates - 7.2 ml LE 392 $D_{500} = 1.0$
+ 48 μ l 1:100 dilution

incubate 37°C 5 hrs
cool 4°C 1 hr.
do lifts

FRIDAY 2/4/94

Denature 0.5N NaOH , 1.5M NaCl

Neutralize 0.5M Tris pH 8 , 1.5M NaCl

let air dry

Monday 2/7/94

Primmbridge in $2 \times \text{PIPES buffer}$

42°C incubator

Fragment - HASK38, HESAY29, HTPANOS - RI/No 13

FRIDAY 2/4/94


Fragment Prep

Digest : HASK38
HESAY29
HTPANOS

With EcoRI /XhoI

DNA	10 μ l
10x #2	5 μ l
10x BSA	5 μ l
RI	2
Xho	2
H ₂ O	24 μ l
	40 μ l

Incubate 37°C, 2 hrs
Run 1.5 μ l




Run on low melt gel + cut out

Monday 2/7/94

Rem on 2nd LMPol
with 1 Kb

Rem 80V 2hrs - 2 1/2 hrs



Cost and bands:

HALSK38	RI/xho	~600 bp.
HEBA 29	RI/xho	1.1 Kb
	RI/xho	1.3 Kb
HTPA 03	RI/xho	1.1 Kb

Use TRANS from Random Prime
Rymo

HIV-1 Screening

(pg 8)

25

Tuesday 2/8/94

Pour off Probe

Wash - 0.2XSSC ; 0.1% SDS

Wash 3 x 50°C
0.2XSSC ; 0.1% SDS
30min

Put on film
- 80°C ON.

Wednesday 2/9/94

Develop film.

Total of 37 (+) clones

- 1 - 4
- 2 - 2
- 3 - 4
- 4 - 3
- 5 - 4
- 6 - 3
- 7 - 4
- 8 - 5
- 9 - 0
- 10 - 0
- 11 - 3
- 12 - 5

Pick plug into
600 500 ul SM Buffer
store 4°C

Do 2° Screening on all (+) clones

plates 1ul → 300ul LE392
plate onto 100mm plate

HTIPAN08 Screening
2/15/94 Tuesday

lul of phage \rightarrow 300ul LE392 Δ low = 1.0.

incubate 37°C 15 min
Plate onto NZ4 Plate (100mm)
in 5ml LB + 0.75% Agarose
melted 55°C

incubate 37°C 95 hrs
Place plate at 4°C.

2/16/94 Wednesday

Do plate lifts
Denature: 2 min
Neutralize: 5 min

Phybridge @ 42°C in 2x P/PES

Random Prime fragment

Low Melt. qd fragment 2ul
5x Buffers 10ul
34ul

Heat 100°C 5 min
Quench. Chill. Quench again

Add:

5x dATP Buffer 10ul
 α -³²P dATP 5ul
Klenow 1ul
60ul

incubate 37°C 15 min

1+ IPANOS Screening

31

Wednesday 2/16/94

-1.0

nm)
35

Put through G 25 column

Count 1 ul
2.23 x 10⁶ counts/ul

Use 20ul → 40ul Hyb buffer
Boil probe 100°C 5 min in 50ul
Salmon Sperm DNA - Quick Chill

Membrane 42°C ON

Thursday 2/17/94

Wash 0.2 x SSC; 0.1% SDS

Wash 3 x 55°C
0.2 x SSC; 0.1% SDS
30 min / wash

Put on film -80°C

Develop after 3 hrs

FRIDAY 2/18/94

⊕ Clones

~~1/1~~
~~2/2~~
3/4

~~4/3~~
~~5/4~~
6/5

~~7/4~~
~~8/5~~
9/6

10/8

37 ⊕ clones

1000 clones

On Mem

Rescue

	1-3	2-1	3-4	5-2	7-1	8-2	11-1
12-3	1-2	2-2	4-1	5-3	7-2	8-3	11-2
12-4	1-1	3-1	4-2	5-4	7-3	8-4	11-3
12-5	1-4	3-2	4-3	6-2	7-4	6-1	12-1
		3-3	5-1	6-3	8-1	8-5	12-2

H12A006 2nd Screening

2/22/94 Tuesday

D₅ min. ReconsPhage 100 μ lEXA-~~st~~ 4 μ lX-1 Blue OD₆₀₀ = 1.0 - 200 μ l

Incubate 37°C 15 min

Add 3 ml LB

Incubate 37°C 2 hrs with aeration

Heat 70°C 20 min
Spin 3K 10 minTransfer to fresh tubes
Store at 4°C

2/23/94 Wednesday

5 μ l \Rightarrow 200 μ l SOLR OD₆₀₀ = 1.0

Incubate 37°C 15 min

Plate 100 μ l onto LB Amp + IPTG/x

Incubate 37°C O/N

DO PCR on Phage mid - (see pg 59)

2/24/94 Thursday

Plates over grown

Re-plate -

Add 10 mM HgSO₄ to 1 ml

Incubate 37°C 15 min

Plate 75 μ l on LB Amp + IPTG/x

Incubate 37°C O/N

(pg 59)

H1rANOS

2° Screening

57

(pg 32)

FRIDAY

2/25/94

Pick clones of each.
into LBT Amp
Incubate 4 hrs
Put at 4°C

- 96 Well Dish
37°C w aeration

Monday

2/26/94

Incubate 37°C 2 hrs
Do PCR

M13F	0.1
M13R	0.1
Culture	1
HIPAN 1575	0.1
10x dNTP	2.1
10x PCR	2.1
Tag	0.15
H ₂ O	15.35
	21

Use HIPANOS
Plasmid DNA
as (+) Control.

PCR Program 69.

Tuesday

3/27/94

Run 10 µl on gel with
1 kb ladder

HIPANOB Screening

pg 32

Wednesday

59
2/23/14

PCR Phagemid


Primer = HTPANOB GGS

Primer	0.1
M13F	0.1
10x PCR	2.1
10x dN	2.1
Phagemid	1
H ₂ O	15.45
Taq	0.15
	21.

⊕ Control
HTPANOB plasmid
DNA.

PCR Program #69

Run 10ul + 1Kb Ladder.



see pg 32

lecta
ml
np

Min

3
1
5
4
2
8
1
2
1
2
53
1
25
56

~~3/1/94~~ Monday

Pg

(56)


3/1/94 PABed Tuesday

inoculate into 3rd TB + Amp for
 Boring mini prep
 incubate 37°C w/aeration

3/2/94 Wednesday

Do boring mini prep

Run Gels on 1% gel w. in Kb ladder
 and pBSK



3/3/94 Thursday

PCR Plasmid DNA

Dilute 1/2000 in H₂O use 1ul in PCR

MRF	0.1
HPAN 157-S	0.1
10x dNTP	2.1
10x PCR	2.1
Taq B	0.15
H ₂ O	15.45
DNA	1
	21ul

① Control
+ 10x dNTP
+ plasmid

Pg 67

HT...N Screening
Pg 60

Thursday

67
3/3/94

PCR Program #69.

Run on



• have
Sequenced
Reverse

Make internal
primer

2, 3, 4, 7, 9, 13, 14, 15,
22, 23, 24, 25, 28, 33

2	S01	13	S06	24	S11
3	S02	14	S07	25	S12
4	S03	15	S08	28	S13
7	S04	22	S09	33	S14
9	S05	23	S10		

HTPAN 03, 04, 13 & 14 look good.

Sequence, with gene specific Primers
3' → 5'

Forward Sequence did not match

Re screen

68

HTPANOS SCREENING

4/19/94 Tuesday

Plate 20 plates
~40,000 plaques / plate 1×10^7 pfs/ml.0.8 ul phage into 6 ml LE392 $OD_{600} = 1$ Incubate 37°C 15 min

Put into melted LB + Agarose (0.75%)

Plate onto 150 mm NZY Plates

Incubate 37°C 5-6 hrsPut at 4°C O/N.

4/20/94 Wednesday

Denature 2 min

Neutralize 5 min

Ph Prehybridize 2 hr -42°C
at 2x PIPESRandom Prime
HTPANOS xho/RI

DNA	1 ul
5x Primers	10 ul
H ₂ O	23 ul
	<u>34</u>

Heat 100°C 5 min

Quick Chill Spin

Chill on ice

Add	10 ul	5x Buffer	dATP
	5 ul	$\gamma\text{-P}$ dATP	
	1 ul	Klenow	

MPANOS

Screening

69

Wednesday 9/20/94

Incubate 37°C 10 min

Put through G-25 column

Count cpl.

$\sim 1.3 \times 10^4$ 30 $\mu\text{l.}$

Boil 5 min

Chill on ice

Add to 2X PIPES H_2O

Incubate 42°C O/N

Thursday 9/21/94

Wash 0.2X SSC, 0.1% SDS.

Wash 3X at 65°C
0.2X SSC, 0.1% SDS.

Put on film -80°C
Signal not very strong - so
leave until Monday

Monday 9/25/94

Develop film

HTP - 24
HTN - 13

70

TRAN 08

4/25/94 Monday

Pick plugs into 300ul of
3m.

Incubate at Room 16°C ON.

4/26/94 Tuesday

PCR Room Phase

M13F	0.2
157-S	0.2
10X dN	2.5
10X PCR	2.5
H ₂ O	14.4
Taq	0.2
Phase	5
	25

Program # 58



Do 2° screens of ⊕ clones

PA
PA

HTPAN 08

129

(pg 70)

Wednesday

4/27/94

Use small NZY plates -

* No plates so had to make own
NZY Plates *

Thursday

4/28/94

LE392.0000 = 0.5 - 100 μ l
+ 1.500 Dilution of phage.

Incubate 37°C 6 hrs
Cool at 4°C ON

FRIDAY

4/29/94

Do plate lifts

place orientation marks.

Monday

4/5/94

Denature 2min 0.5N NaOH / 1.5M NaCl

Neutralize 5min 0.5M Tris pH 8.0 / 1.5M NaCl

Wet on 2x SSC 5min

Allow filters to dry

133

H7PANOB

5/3/94 Tuesday

b

Probe Blots

80°C 1 hour

make probe

H7PANOB xhs /RI

DNA	1ul
5x PRIMERS	10ul
H ₂ O	23ul
	<hr/> 34ul

Heat 100°C 5 min

Quick Spin

Place on ice

Add

5x dATP Buffer	10ul
³² P dATP	5ul
Isolator	1ul
	<hr/> 5ul

Incubate 37°C 10 min

put through Gels 6-25 column

Count 1ul

1.04 x 10⁶ counts /ul ~ 50ul~~used 200~~

Use 30ml Hyb Buffer

Prehybridize Blots	1 hr in
in 8 PIPES	Buffer at 42°C

Pour off Prehyb add the Hyb buffer -

Tuesday

5/3/94

Heat

Probe + 50ul Salmon Sperm DNA
at 100°C 5min

Quick Chill

Add 70ul to Hyb buffer w/ lifts

Incubate 42°C O/N w/ lifts

Slight Agitation

Wednesday 5/4/94

Wash Blots -

Rinse 1x 0.2xSSC / 0.1% SDS

Wash 3x in

0.2xSSC + 0.1% SDS
at 65°C

to

No Film Cassettes

so store at RT in 0.2xSSC + 0.1% SDS

Thursday

5/5/94

Put on film

-80°C for 1 hrs

Develop film

Picked 4 clones

H. PANOS

5/6/94

Thursday PEIDAY

Have full length

found that 504, 513, 514

fit into a contig &

When Blasted against HGS/TIGR
Database returned a clone

which when blasted against

HGS/TIGR EST. Produced a match
to HTPAN.

Sequence to get full length

5/8/94

Make Oligos to sequence

More 3' - forward Forward end

5/10/94

Monday

Submit for sequencing

5/11/94

Tuesday

Transform Df/5x low pressure
& Do maxi Plasm.

5/12-5/13 off

- Carrie took plasmid out of
37°C to 4°C

(pg 136)

HTPAN 08

(10/10/132)

5/16/99 Monday

Pick & Clous of lawn

SO₄, S₁₃ & S₁₄ into LB Amp
in 1/2 well each

Incubate 37°C ON

5/17/99 Tuesday

Run PCR using M13 F & Rev.

		25X
M13 F & M13 R	by vol	0.1
10x dNTP		2
10x PCR		2
H ₂ O		13.7
Topo		0.2
Antisense		2
		25
		50
		50
		342.5
		50
		—

Run PCR Program #69

95°C 5min

95°C 20sec

55°C 20sec


72°C 1min

72°C 7.5min

4°C

20X

Run 10µl on gel with 1kb ladder



Tuesday 5/17/94

Inoculate 400 ml
TB + Amp with
S04 - 3
S13 - 5

Incubate 37°C w/ aeration

Wednesday 5/18/94


Main Prep

Spin Culture 5K 10 min
Pour off
Resuspend in 15 ml Rinsuspension
Buffer
Add 15 ml Cell Lysis Buffer
Add 15 ml Neutralization Buffer
Spin 8K 10 min
Transfer to fresh Tube
Add Isopropanol to 0.6 volumes
Mix
Spin 9K 30 min
Pour off
Resuspend pellet in 1 ml TE
Transfer to 2 - 1.5 ml eppendorf
- 2x Phenol / SEVAG
- 2x SEVAG
Ethanol ppt 1/10 vol 3M NaAc.
2x vol Ethanol
Inc 15 min
Spin 15 min
Pour off supernatant
2x 70% Ethanol wash
1x 100% Ethanol wash
Dry pellet & Resuspend
Read OD 260/280
Run 2 ul on gel with markers

138

HTPAN08

5/18/94 Wednesday

5/19/94 Thursday
Read O.D.₂₆₀ / 280

1:100 Dilution

	abs 260.0 nm	abs 280.0 nm	big abs 320.0 nm	260.0 nm 280.0 nm	280.0 nm 260.0 nm
HTPAN08	0.0011	0.0007	0.0005	4.0362	0.2478
HTPAN08S04	0.0679	0.0344	-0.0035	1.8790	0.5322
HTPAN08S13	0.1344	0.0674	-0.0049	1.9259	0.5192
	0.6327	0.3200	-0.0027	1.9689	0.5079

0.34 ug/lul.
0.6 ug/lul.
3.16 ug/lul.

Submit for Sequencing

HTPAN08

Fc

Rc

R157S

RP01b

R318S

RP03b

R166S

RP05b

R338S

RP06b

F488A

RP07b

HTPAN08S04

Fc

Rc

RP01

RP03b

RP05

RP06

RP07

HTPAN08S04 pol RP10

RP10

HTPAN08S13

Fc

Rc

RP01

RP03b

RP05

RP06

RP07

HTPAN08S13 pol RP09

RP09

RP10

Pg 53 Book #168
AMC #6

Sample ID	abs 260.0 nm	abs 280.0 nm	bkg abs 320.0 nm	260.0 nm		280.0 nm			
				260.0 nm	280.0 nm	260.0 nm	280.0 nm		
1mer 1HE3629	0.0558	0.0297	-0.0001	1.8768	0.5328	0.46ug/ml	73.4 pmol/l		
1mer 2HPR480	0.0508	0.0336	-0.0005	1.5008	0.6683	0.42ug/ml	67 pmol/l		

2/2/94

2/7/94

Sample ID	abs 260.0 nm	abs 280.0 nm	bkg abs 320.0 nm	260.0 nm		280.0 nm			
				260.0 nm	280.0 nm	260.0 nm	280.0 nm		
1HTAN08F 3185	0.0473	0.0364	0.0044	1.2674	0.7890	0.39ug/ml	59.1 pmol/l		
2HTAN08R 4050	0.0350	0.0176	0.0054	1.7582	0.5687	0.29ug/ml	46.1 pmol/l		
3HTAN08R 665	0.0499	0.0270	0.0009	1.8805	0.5318	0.41ug/ml	65.7 pmol/l		

2/11/94

Sample ID	abs 260.0 nm	abs 280.0 nm	bkg abs 320.0 nm	260.0 nm		280.0 nm			
				260.0 nm	280.0 nm	260.0 nm	280.0 nm		
1mer 1HMAA 51-2055	0.0430	0.0324	0.0065	1.4044	0.7121	0.22ug/ml	37 pmol/l		
1mer 2HMTN 17-2595	0.0384	0.0237	-0.0000	1.6219	0.6166	0.19ug/ml	30 pmol/l		

2/15/94

Sample ID	abs 260.0 nm	abs 280.0 nm	bkg abs 320.0 nm	260.0 nm		280.0 nm			
				260.0 nm	280.0 nm	260.0 nm	280.0 nm		
1 5512	0.0300	0.0194	0.0023	1.6231	0.6161				
2 5513	0.0324	0.0164	0.0011	2.0444	0.4892				

3/4/94

Sample ID	abs 260.0 nm	abs 280.0 nm	bkg abs 320.0 nm	260.0 nm		280.0 nm			
				260.0 nm	280.0 nm	260.0 nm	280.0 nm		
1HBE669 95	0.0551	0.0346	0.0015	1.6215	0.6187	0.3ug/ml	50 pmol/l		
2HBE668 97	0.0586	0.0288	0.0021	2.0539	0.4889	0.3ug/ml	51.6 pmol/l		

150

3/8/94

Sample ID	abs 260.0 nm	abs 280.0 nm	bkg abs 320.0 nm	260.0 nm	280.0 nm
				280.0 nm	260.0 nm
1153					
25677 HASK POI	0.0691	0.0450	-0.0006	1.5287	0.6541
35676 HTPANOS	0.0935	0.0524	0.0017	1.8117	0.552089 pmol/ul
45675 HEBAY POI	0.1129	0.0588	0.0023	1.9576	0.5108 102 pmol/ul
55678 HCAN POI	0.0750	0.0528	0.0022	1.4427	0.6931 79 pmol/ul
	0.0827	0.0505	0.0027	1.6746	0.5972 87 pmol/ul

3/10/94

				260.0 nm	280.0 nm	
Sample ID	abs 260.0 nm	abs 280.0 nm	bkg abs 320.0 nm	280.0 nm	260.0 nm	
5697	1HE9MF3 142	0.0857	21mer 0.0477	0.0010	1.8121	0.5519 81 pmol/ul
5699	2HE7SE29 231	0.0604	20mer 0.0385	-0.0009	1.5565	0.6425 60 pmol/ul
5688	3HTPANOS 485	0.0536	21mer 0.0283	-0.0017	1.8443	0.5422 51 pmol/ul
5685	4HPLB52 222	0.0551	20mer 0.0375	-0.0010	1.4573	0.6882 55 pmol/ul
5686	5HPRAL36 551	0.0268	18mer 0.0179	-0.0000	1.4972	0.6679 30 pmol/ul

3/10/94

Sample ID	abs	abs	bkg abs	260.0 nm	280.0 nm
	260.0 nm	280.0 nm	320.0 nm	280.0 nm	260.0 nm
5756 1HTPAUCB 124A	0.0467	0.0270	0.0018	1.7815	0.5613 61 pmol / ul

3/17/94

Sample ID	abs	abs	bkg abs	260.0 nm	280.0 nm	
	260.0 nm	280.0 nm	320.0 nm	280.0 nm	260.0 nm	
5780	1740186	0.0746	0.0417	-0.0003	1.7828	0.5609 19mer 70 pmol/ul
5781	2741087	0.0770	0.0533	-0.0005	1.3833	0.7228 19mer 81 pmol/ul
5782	3741030	0.0762	0.0508	-0.0021	1.4809	0.6752 19mer 80 pmol/ul
5783	4741034	0.0847	0.0564	-0.0024	1.4806	0.6754 21mer 80 pmol/ul
5		0.0900	0.0600	0.0000	0.0000	
5764	6741036	0.0463	0.0310	-0.0028	1.4579	0.6859 20mer 46 pmol/ul
6785	7741037	0.0732	0.0437	-0.0028	1.6350	0.6116 19mer 77 pmol/ul
5786	8741039	0.0658	0.0527	-0.0019	1.2388	0.8073 19mer 69 pmol/ul

3/18/94

Sample ID	abs	abs	bkg abs	260.0 nm	280.0 nm
	260.0 nm	280.0 nm	320.0 nm	280.0 nm	260.0 nm
5787 1	-0.0007	-0.0016	0.0005	0.6323	1.8786
5788 2741040 R 2010	0.0538	0.0362	-0.0022	1.4574	0.6882 59 pmol/ul
5799 3741042 R 2010	0.0137	0.0056	-0.0021	2.0592	0.4856 15.2
5800 4741047 R 2010	0.0111	0.0045	-0.0031	1.8718	0.5343 11.1 pmol/ul
5809 5HETAS6 R 2010	0.0061	0.0020	-0.0034	1.7501	0.5714 64 pmol/ul
5810 6HETAS6 R 2010	0.0377	0.0129	-0.0036	2.5057	0.3991 39 pmol/ul
5811 7HAFBOS R 2010	0.0244	0.0129	-0.0043	1.6674	0.5998 24.4 pmol/ul

**This Page is Inserted by IFW Indexing and Scanning
Operations and is not part of the Official Record**

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

- ☐ **BLACK BORDERS**
- ☐ **IMAGE CUT OFF AT TOP, BOTTOM OR SIDES**
- ☐ **FADED TEXT OR DRAWING**
- ☐ **BLURRED OR ILLEGIBLE TEXT OR DRAWING**
- ☐ **SKEWED/SLANTED IMAGES**
- ☐ **COLOR OR BLACK AND WHITE PHOTOGRAPHS**
- ☐ **GRAY SCALE DOCUMENTS**
- ☐ **LINES OR MARKS ON ORIGINAL DOCUMENT**
- ☒ **REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY**
- ☐ **OTHER:** _____

IMAGES ARE BEST AVAILABLE COPY.

As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.